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
## Intended use

Prepared sample tubes are used with needles for venous blood collection. The prepared sample tubes are used for collecting, transporting and processing blood samples in clinical laboratories. The products are intended for use in a professional environment by qualified medical personnel.

## Product description

Prepared sample tubes consist of a plastic tube with a flat or round bottom in various sizes and a colour-coded plastic screw cap with or without an integrated membrane or a push cap, and various additives (preparations). The selection of the additive depends on the analytical test method. It is specified by the manufacturer of the test reagents and/or the analysis equipment on which the test is conducted.

Colour codes of the prepared sample tubes, screw caps and push caps:

Additive	Letter code	Cap colour based on ISO 6710*	Cap colour based on BS 4851**
<b>Serum tube</b>			
Prepared sample tubes with serum	CAT	Red	White
Prepared sample tubes with serum gel	CAT	Brown	Brown
<b>Lithium heparin tube</b>			
Prepared sample tubes with lithium heparin	LH	Green	Orange
<b>EDTA tube</b>			
Prepared sample tubes with K3 EDTA	K3E	Violet	Red
Prepared sample tubes with K3 EDTA, 3 ml, membrane cap	K3E	Violet	Red
<b>Tube with glycolysis inhibitors</b>			
Prepared sample tubes with fluoride/heparin, 2 ml	FH	Grey	Yellow
<b>Tube with citrate</b>			
Prepared sample tube with citrate	9NC	Blue	Green

\*DIN EN ISO 6710: Single-use containers for human venous blood specimen collection

\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, revoked

## SARSTEDT Prepared sample tubes Serum CAT / Serum Gel CAT

The prepared sample tube Serum CAT contains a synthetic granulate coated with a clotting activator (silicate) and is used for serum collection\*\*. The synthetic granulate settles between the blood clot and serum during centrifugation. The prepared sample tube Serum Gel CAT contains a polymer-based gel, in addition to the synthetic granulate coated with a clotting activator (silicate), and is used for serum collection\*\*. The gel settles between the blood clot and serum during centrifugation and forms a stable separation layer during transport, storage and analysis. Serum is used as sample material in routine biochemical and immunological testing and in serology. After collecting blood, let the blood coagulate in the prepared sample tube Serum CAT / Serum Gel CAT for 30 minutes in an upright position to ensure a well-formed horizontal separation layer after centrifugation. The recommended time is based on an intact coagulation process. Blood from patients with coagulation disorders due to illness or under anticoagulant therapy requires more time for complete coagulation.

## SARSTEDT Prepared sample tubes Lithium Heparin LH

The prepared sample tube Lithium Heparin LH contains a synthetic granulate coated with the anticoagulant lithium heparin and is used for collecting whole blood and plasma\*\*. The synthetic granulate settles between the blood cells and plasma during centrifugation. Anticoagulation with heparin occurs through activation of antithrombin.

**NOTE: No lithium determinations may be carried out from the prepared sample tube Lithium Heparin LH.**

## SARSTEDT Prepared sample tubes EDTA K3E

The prepared sample tube EDTA K3E contains the anticoagulant K3 EDTA and is used for collecting whole blood\*\*. K3 EDTA is in a spray dose form. EDTA whole blood is used as sample material for haematological analyses. Blood smears should be made within four hours of blood collection. Anticoagulation occurs through calcium ion complex formation by EDTA. Blood from the prepared sample tube EDTA K3E can also be used for routine immunohaematology analyses and tests for infectious diseases. The user is responsible for validating the suitability of the sample material for these analyses with the appropriate test reagents / analysis equipment incl. storage conditions.

### SARSTEDT Prepared sample tubes Fluoride Heparin FH

The prepared sample tube Fluoride Heparin FH contains the anticoagulant heparin as well as fluoride for glycolysis inhibition and is used for collecting plasma\*\*. Fluoride heparin is in a spray dose form. Glycolysis inhibition through fluoride alone starts approx. 2 hours after filling the blood collection tube and shows its full effect after approx. 4 hours. On average, therefore, up to 6 % of glucose breaks down after a few hours and up to 10-15 % after 24 hours in whole blood.

**NOTE: Fluoride can cause increased haemolysis. Further information on substances that can cause interference can be found in the instructions for use of the relevant test reagent manufacturer.**

### SARSTEDT Prepared sample tubes Citrate 9NC

The prepared sample tube Citrate 9NC contains the anticoagulant trisodium citrate and is used for collecting whole blood and plasma\*\*. The mixture ratio of citrate to blood is 1:9 – 1 part by volume citrate and 9 parts by volume blood. Correct filling is absolutely necessary for the analysis. Citrate plasma is used as sample material in routine haemostaseology tests. Anticoagulation occurs through calcium ion complex formation by citrate.

**\*\*Suitability of the sample material depends on the analyte and test reagent / analysis equipment. Corresponding recommendations of the manufacturer must be observed.**

### Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures of your facility. In case of exposure to biological specimens or of a needlestick injury, consult a doctor, as HIV, HCV, or HBV or other infectious diseases can be transmitted. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects (needles) for blood collection in suitable sharps containers.
4. Under- or overfilling the prepared sample tubes leads to an incorrect ratio of blood to preparation/additive and can give rise to false analysis results.
5. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.

### Storage

Store the products at room temperature.

### Transport

The products with screw caps are primary receptacles according to the ADR (Packing instruction P650) and IATA regulations.

### Limitations

1. The storage period and temperature of a filled sample tube depend on the stability of the analytes for investigation. The assessment is made by the laboratory, or information is obtained from the manufacturer's instructions for use for the test reagents / analysis device.
2. Even when plasma or serum is separated by centrifuging the prepared sample tubes and/or a barrier is present, not all cells are necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of the analytes. The stability of the analyte should be evaluated with respect to the storage container and the conditions of the relevant laboratory.
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use for the test reagents / analysis device.

### Collecting and handling specimens

**READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW BLOOD.**

#### Equipment required for specimen collection:

1. All required prepared sample tubes, labelled by additive.
2. Gloves, gown, goggles or other suitable protective gear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) Luer needle for collecting venous blood. Possibly also non-prepared collection receptacles with a multi-adapter.
5. Disinfectant to clean the collection site (follow the guidelines of the facility for sterile specimen collection for preparing the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

## Handling:

COLLECT THE SPECIMEN IN LINE WITH THE ESTABLISHED GUIDELINES AND PROCEDURES IN YOUR FACILITY.

1. Open the prepared sample tube and set the cap aside.
2. Place your specimen material in the prepared sample tube. Ensure that the prepared sample tube is not filled with more or less than the nominal volume.
3. Close the prepared sample tube with the corresponding cap. Invert all samples with an additive gently several times.

## Notes on coagulation of prepared sample tubes with serum / serum gel:

After collecting blood, let the prepared sample tubes with serum/serum gel coagulate for 30 minutes in an upright position to ensure a clear separation layer after centrifugation. The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to illness or under anticoagulant therapy requires more time for complete coagulation.

## Centrifugation

### WARNING!

**SARSTEDT prepared sample tubes are designed for 3,500 x g. Sample tubes with a diameter of 8.5 mm (1.2 ml) are an exception, as they are currently validated up to 2,500 x g. Only use suitable tube holders or inserts. Centrifuging prepared sample tubes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the prepared sample tubes, which can release potentially infectious material.**

Centrifuge inserts should be selected according to the size of the prepared sample tubes used. The relative centrifugal force has the following ratio to the selected rpm:

$$RCF = 11.2 \times r \times (\text{rpm}/1,000)^2$$

RCF: "Relative centrifugal force", also called "g-force",

"rpm": "Revolutions per minute",

"r": "Radius of rotation from the centre of the centrifuge to the bottom of the prepared sample tube" (cm).

Prepared sample tubes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.

Prepared sample tubes with gel are designed solely for use in centrifuges with swing-out rotor. Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended.

Prepared sample tubes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user themselves.

Ensure that the prepared sample tubes fit into the centrifuge inserts.

Prepared sample tubes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

**CAUTION! Do not remove broken prepared sample tubes by hand.**

*Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.*

Prepared sample tubes	Relative centrifugal force (g)			
	2000 x g	2500 x g	3000 x g*	3500 x g*
Serum	10 min	10 min	6 min	4 min
Serum Gel	15 min	10 min	4 min	4 min
Lithium Heparin	10 min	10 min	7 min	7 min
K3 EDTA*	n.v.	n.v.	7 min	6 min
Citrate	9 min	8 min	7 min	6 min
Fluoride/Heparin	9 min	8 min	7 min	6 min

n.v. = not validated

\* Applies to all prepared sample tubes except those with 8.5 mm diameter (1.2 ml).

## Information on separation material

The flow properties of the separation material (gel) are temperature-dependent. If prepared sample tubes are cooled before or during centrifugation, a continuous layer of the separation gel is no longer ensured. Prepared sample tubes must not be centrifuged again / re-centrifuged.

## Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled prepared sample tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

## Product-specific standards and guidelines as amended/updated

CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.


CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

\*CLSI (Clinical and Laboratory Standards Institute)

## Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE mark
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Keep away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.


All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, изтерлен.

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CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.  
 CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.  
 CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.  
 CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.  
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".  
 \*CLSI (Clinical and Laboratory Standards Institute)

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
\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, atšauktas

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- CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
- CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.
- CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.
- CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
- WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
- \*CLSI (Clinical and Laboratory Standards Institute)

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
\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, geri çekildi

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