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SARSTEDT Micro-Probengefäße

DE	Gebrauchshinweis – SARSTEDT Micro-Probengefäße	2 - 6
EN	Instructions for use – SARSTEDT Micro sample tubes	7 - 11
BG	Инструкции за употреба – SARSTEDT микропробетки за проби	12 - 16
CS	Návod k obsluze – Mikrozkušavky na vzorky SARSTEDT	17 - 21
DA	Brugsanvisning – SARSTEDT mikroprøverør	22 - 26
EL	Οδηγίες χρήσεως – Μικροφιαλίδιο δειγμάτων SARSTEDT	27 - 31
ES	Instrucciones de uso – Microtubos SARSTEDT	32 - 36
ET	Kasutusjuhend – Mikrokatsutid SARSTEDT	37 - 41
FR	Mode d'emploi – Microtubes SARSTEDT	42 - 46
HR	Uputa za upotrebu – SARSTEDT Mikropruvete za uzorke	47 - 51
HU	Használati utasítás – SARSTEDT mikro-mintaedélyek	52 - 56
IT	Istruzioni d'uso – Micro-provette SARSTEDT	57 - 61
KO	사용 설명서 – SARSTEDT 마이크로 검체 용기	62 - 66
LT	Naudojimo instrukcijos – SARSTEDT mikromėgintuvėliai	67 - 71
LV	Lietošanas norādes – SARSTEDT paraugu mikrotrauki	72 - 76
NL	Gebruiksaanwijzing – SARSTEDT micro-monsterbuisjes	77 - 81
NO	Bruksanvisning – SARSTEDT mikro-prøverør	82 - 86
PL	Instrukcja obsługi – Mikropróbówki SARSTEDT	87 - 91
PT	Instruções de utilização – Microtubos de amostras SARSTEDT	92 - 96
RO	Instrucţiuni de utilizare – Microtuburi pentru recoltare SARSTEDT	97 - 101
RU	Инструкция по применению – Микропробирки SARSTEDT	102 - 106
SK	Návod na Použitie – SARSTEDT mikroskúšavky na vzorky	107 - 111
SL	Navodila za uporabo – Mikropruvete za vzorce SARSTEDT	112 - 116
SV	Bruksanvisning – SARSTEDT preparerade mikrorör	117 - 121
TH	คำแนะนำในการใช้งาน – หลอดเก็บตัวอย่างขนาดเล็กของ SARSTEDT	122 - 126
TR	Kullanım için talimatlar – SARSTEDT mikro test tüpleri	127 - 131
ZH	使用说明 – SARSTEDT 微量样本管	132 - 136

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Intended use

Micro sample tubes are used with Micro needles as a system for venous blood collection. The Micro sample tubes are used for collecting, manual in-house transport and processing of blood samples in the clinical laboratory. The products are intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

Micro sample tubes consist of a flat-bottomed plastic receptacle and a colour-coded plastic screw cap or a soft push cap with loop and various additives (preparations). The selection of the additive depends on the analytical test method. It is specified by the manufacturer of the test reagents and/or the analysis equipment on which the test is conducted.

Colour codes of the Micro sample tube screw caps or soft push caps:

Additive	Letter code	Cap colour based on ISO 6710	Cap colour based on BS 4851*
Serum tube			
Micro sample tube with serum	CAT	red	white
Micro sample tube with serum gel	CAT	brown	brown
Lithium heparin tube			
Micro sample tube with lithium heparin	LH	green	orange
EDTA tube			
Micro sample tube with K3 EDTA	K3E	violet	red
Micro sample tube with K3 EDTA incl. membrane closure	K3E	violet	red
Tube with glycolysis inhibitors			
Micro sample tube with fluoride heparin/glucose	FH	grey	yellow
Citrate tube			
Micro sample tube with citrate 9NC	9NC	blue	green

*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry.

SARSTEDT Micro sample tube Serum/Serum gel

Micro sample tubes with serum contain a clotting activator and are used for collecting serum. After centrifugation, the gel in the Micro sample tube Serum/Gel forms a layer separating the coagulum and the serum. Micro sample tubes with serum/serum gel are intended for clinical chemistry tests from serum.

SARSTEDT Micro sample tube Lithium heparin gel

Micro sample tubes with lithium heparin are coated with the anticoagulant lithium heparin. They are used for plasma recovery and are intended for clinical-chemical examinations.

SARSTEDT Micro sample tube K3 EDTA

Micro sample tubes with K3 EDTA are coated with the anticoagulant K3 EDTA and are intended for haematological examinations.

SARSTEDT Micro sample tube Fluoride/heparin

Micro sample tubes with fluoride/heparin are coated with the anticoagulant heparin and contain fluoride for glycolysis inhibition. They are used for glucose determination.

NOTE: Fluoride (used for measuring glucose) can cause an increase in haemolysis. Further information on substances that can cause interference can be found in the respective instructions for use of the assay manufacturer.

SARSTEDT Micro sample tube Citrate 9NC

Micro sample tubes with citrate 9NC contain the anticoagulant citrate and are intended for haemostasis investigations.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures of your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects (needles) for blood collection in suitable sharps containers.
4. If blood is collected from an intravenous (IV) line, it must be ensured that the line is properly flushed in accordance with the institution (=cleaned of the IV solution) before beginning blood collection. Properly flushing the line prevents faulty laboratory results.
5. Under- or overfilling the Micro sample tubes leads to an incorrect ratio of blood to the preparation/additive and can result in incorrect analysis results.
6. Do not use Micro sample tubes after the expiry date. The expiry date of the Micro sample tubes is the last day of the month and year indicated.

Storage

Store the products at room temperature.

Transport

The products with screw cap are primary receptacles according to the ADR (packing instruction P650) and IATA regulations.

Limitations

1. If blood samples are to be stored in Micro sample tubes, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. Even if plasma or serum is separated by centrifuging the Micro sample tubes and/or a barrier is present, all cells are not necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of the analytes. The stability of the analyte should be evaluated with respect to the storage container and the conditions of the respective laboratory
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW BLOOD.

Materials needed for collecting a specimen/blood:

1. All required Micro sample tubes, labelled by additive.
2. Gloves, gown, goggles or other suitable protective wear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. Micro needle for venous blood collection.
5. Disinfectant to clean the collection site (follow guidelines of the facility for sterile specimen collection for preparing the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

Recommended order of collection:


1. Tubes for coagulation tests of serum with/without gel
2. Tubes with heparin
3. Tubes with EDTA
4. Tubes with glycolysis inhibitors
5. Containers with other additives

NOTE: Follow the regulations of your facility for the order of collection.


Venous blood collection in drip technique using a Micro needle

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING BLOOD COLLECTION.

- 1 Unscrew the screw cap of the Micro sample tube and set it aside.
- 2 For blood collection with a Micro needle, collect the venous blood in the Micro sample tube using the drip technique. Stop blood collection when the Micro sample tube is filled to the fill mark.
- 3 After blood collection, first remove the Micro needle
- 4 ...and discard it.
- 5 Close the Micro sample tube with the attached screw cap.

- 1 Open the soft cap of the Micro sample tube.
- 2 For blood collection with a Micro needle, collect the venous blood in the Micro sample tube using the drip technique. Stop blood collection when the Micro sample tube is filled to the fill mark.
- 3 After blood collection, first remove the Micro needle
- 4 ...and discard it.
- 5 Close the Micro sample tube with a soft cap.

Notes on coagulation/ on Micro sample tubes Serum:

After collecting blood, let the blood coagulate in the **Micro sample tubes Serum/Serum Gel** for 30 minutes in an upright position to ensure a clear layer of separation after centrifugation.

The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or coagulation-inhibiting treatment requires more time for complete coagulation.

Centrifugation

WARNING!

Centrifuging Micro sample tubes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Micro sample tubes, which can release potentially hazardous substances.

Centrifuge inserts are to be selected according to the size of the Micro sample tubes used. The relative centrifugal force has the following relationship to the RPM (revolutions per minute) set:

$$RCF \text{ (g force)} = 11.2 \times r \times (\text{RPM}/1000)^2,$$

"RCF": "Relative centrifugal force",

"RPM": "Revolutions per minute" (min^{-1}),

"r" [in cm]: "Centrifuge radius from the centre of the centrifuge to the bottom of the Micro sample tube.

Micro sample tubes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.

Micro sample tubes with gel are designed solely for use in centrifuges with swing-out rotor. Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended. In addition, centrifugation should be conducted no later than 2 hours after collection, as this could otherwise yield false results.

Micro sample tubes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user.

It must be ensured that the Micro sample tubes fit well in the centrifuge inserts. Micro sample tubes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

CAUTION! Do not remove broken Micro sample tubes by hand.

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

Image	Preparation	Times and g values (=RCF*)				
		Standard recommendation		Alternative recommendation		
	Serum	5 min	10,000 x g	or	10 min	
	Lithium heparin		2,000 x g			
	Glucose fluoride		10,000 x g			
	Serum gel ¹		1,500 x g			
	Citrate	10 min				

¹ For tubes prepared with gel, we recommend using swing-out rotors.

*RCF: Relative centrifugal force

Centrifugation at 20 °C

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled Micro sample tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

DIN EN ISO 6710: Single-use containers for human venous blood specimen collection

CLSI* C38 "Control of Preanalytical Variation in Trace Element Determinations", Approved Guideline.

CLSI* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline .

CLSI* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

Article number

Batch number

Use by

CE marking

In vitro diagnostic device

Follow the instructions for use

If reused: Risk of contamination

Store away from sunlight

Store in a dry place

Manufacturer

Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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