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Intended use

Prepared tubes for HbA_{1c} testing using Tina-quant® are used together with an End-to-End capillary for capillary blood collection for manual or automated measuring of the long-term marker for blood glucose (Tina-quant® HbA_{1c}). Prepared tubes for HbA_{1c} testing using Tina-quant® are used for the collection, transport and processing of capillary blood specimens in a clinical laboratory. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

Prepared tubes for HbA_{1c} testing using Tina-quant® consist of a plastic tube with a flat bottom and a black plastic membrane screw cap and are used together with the enclosed 10 µl End-to-End capillary. The prepared sample tubes contain 1,000 µl of the haemolysis solution Tina-quant® HbA_{1c}.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils according to the guidelines and procedures in your facility. In case of direct contact with biological specimens or of a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. The products are intended for single use. Dispose of all sharp objects (needles, adapters) for blood collection in suitable containers for hazardous biological waste.
4. Under or overfilling the prepared sample tubes leads to an incorrect ratio of blood to preparation/additive and can give rise to false analysis results.
5. The blood collected with the prepared tubes and then processed is not intended to be re-injected into the human body.
6. Do not use the products after the expiry date. The expiry date of the product is the last day of the month and year indicated.

Storage

Store the products at room temperature.

Limitations

1. If blood samples are stored in the prepared sample tubes for HbA_{1c} testing using Tina-quant®, the stability of the HbA_{1c} should be evaluated by the relevant laboratory or taken from the instructions for use of the analysis device manufacturer.
2. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.


Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW BLOOD.

Preparation for drawing capillary blood and material needed:

1. The required prepared tubes for HbA_{1c} testing using Tina-quant® and the corresponding End-to-End capillaries.
2. Gloves, gown, goggles or other suitable protective gear to protect from bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet or incision lancet.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples).
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

General instructions:

- 1 Prepared tube for HbA_{1c} testing using Tina-quant® consisting of a prepared sample tube and an End-to-End capillary.
- 2 Clean the puncture site with disinfectant. Puncture the skin in a way that ensures good blood flow. Discard the first drop. Then hold the End-to-End capillary horizontally with one end in the middle of the drop of blood and completely fill the End-to-End capillary without air bubbles. Use a swab to wipe any blood off the outside of the capillary.
- 3 Open the prepared sample tube and place the End-to-End capillary filled with blood in it.
- 4 Place the End-to-End capillary filled with blood in the prepared sample tube. Close the prepared sample tube and swirl it gently. Swirl until all the blood from the capillary has gone into solution.

Use the prepared tubes for HbA_{1c} testing using Tina-quant® in accordance with the instructions for use of the test manufacturer and the Standard Operating Procedure (SOP) of your facility. The prepared tubes for HbA_{1c} testing using Tina-quant® can now be used for automated analysis. For use in devices, follow the respective device manufacturer's instructions.

Note:

For the stability of HbA_{1c}, follow the instructions for use of the test reagents / analysis device manufacturer.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled capillary blood collection systems must be disposed of in suitable containers for hazardous biological waste.
4. Dispose of potentially contaminated material according to the regulations and guidelines of the facility.

Product-specific standards and guidelines as amended/updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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CLSI* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection" , Approved Standard.

CLSI* GP42 " Collection of Capillary Blood Specimens" , Approved Standard.

CLSI* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva:World Health Organization, 2002" .

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