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SARSTEDT Schlauch-Segment-Öffner

DE	Gebrauchshinweis – SARSTEDT Schlauch-Segment-Öffner	2 - 4
EN	Instructions for Use – SARSTEDT tube segment opener	5 - 7
BG	Инструкции за употреба – SARSTEDT Отварачка за сегменти на маркуч	8 - 10
CS	Návod k obsluze – Pomůcka na propíchnutí segmentu hadičky SARSTEDT	11 - 13
DA	Brugsanvisning – SARSTEDT slangesegment-åbner	14 - 16
EL	Οδηγίες χρήσεως – Συσκευή ανοίγματος συνοδών σωληνίσκων ασκών αίματος SARSTEDT	17 - 19
ES	Instrucciones de uso – Abridor de segmentos de tubo SARSTEDT	20 - 22
ET	Kasutusjuhend – Voolikusegmeni avaja SARSTEDT	23 - 25
FR	Mode d'emploi – Perce-tubulure SARSTEDT	26 - 28
HR	Uputa za upotrebu – SARSTEDT otvarač segmenta crijeva	29 - 31
HU	Használati utasítás – SARSTEDT tömlő szegmens nyitó	32 - 34
IT	Istruzioni d'uso – Pungitubo per segmenti SARSTEDT	35 - 37
KO	사용 설명서 – SARSTEDT 튜브 세그먼트 오픈너	38 - 40
LT	Naudojimo instrukcijos – SARSTEDT žarnelės segmento atidarytuvas	41 - 43
LV	Lietošanas norādes – SARSTEDT caurulītes segmenta atvērējs	44 - 46
NL	Gebruiksaanwijzing – SARSTEDT-slangsegmentopener	47 - 49
NO	Bruksanvisning – SARSTEDT slangesegmentåpner	50 - 52
PL	Instrukcja obsługi – Otwieracz do drenów SARSTEDT	53 - 55
PT	Instruções de utilização – Abridor de segmentos de tubo flexível SARSTEDT	56 - 58
RO	Instrucțiuni de utilizare – Deschizător de segment cu tub SARSTEDT	59 - 61
RU	Инструкция по применению – Открыватель сегментных трубок SARSTEDT	62 - 64
SK	Návod na Použitie – SARSTEDT otvárač hadičkového segmentu	65 - 67
SL	Navodila za uporabo – Segmentni odpiráč cevi SARSTEDT	68 - 70
SV	Bruksanvisning – SARSTEDT segmentöppnare	71 - 73
TH	คำแนะนำในการใช้งาน – ฝ่ายเปิดส่วนสายท่อของ SARSTEDT	74 - 76
TR	Kullanım için talimatlar – SARSTEDT tüp segment açıcı	77 - 79
ZH	使用说明 – SARSTEDT 血液片段段样管打开器	80 - 82

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Intended use


The tube segment opener is used for the hygienically clean transfer of blood (erythrocyte concentrate) from the segments of a blood bag into a secondary container (e.g. a Coombs tube) or directly into the test card of a bedside test.

The products are intended for use in a professional transfusion medicine environment by qualified medical and laboratory personnel.

Product description

The tube segment opener is made of plastic and contains a disposable medical needle made of stainless steel. Depending on the design, they are available for tube outer diameters of 11,5 mm, 12 mm or tube inner diameters of 8.2 mm and can be pre-assembled in a tube with a round or pointed bottom. The tube segment opener is a disposable product and has been specially developed for the transfer of erythrocyte concentrates into a secondary tube or a bedside test card.

Product overview:

Type	Description
A	Tube segment opener
B	Tube segment opener for bedside test
C	Tube segment opener with round bottom tube pre-assembled
D	Tube segment opener with pointed bottom tube pre-assembled

Safety information and warnings

1. Check the tube segment opener before use.
2. Do not use a damaged tube segment opener.
3. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
4. Handle all biological samples and sharp blood collecting utensils (needles) in accordance with the guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
5. Dispose of used and contaminated tube segment openers and secondary tubes in suitable disposal containers for biological hazardous substances.
6. Injection and subcutaneous administration are not allowed.
7. Reusing the product can lead to infections, injuries or death.
8. Do not use the tube segment opener after the expiry date. The expiry date of the tube segment opener is the last day of the month and year indicated.

Storage

Store the tube segment opener at room temperature.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE YOU START WITH THE APPLICATION.

Material required for the transfer:


1. Use gloves, gown, goggles or other suitable protective wear to protect against bloodborne pathogens or potentially infectious materials.
2. If necessary, Secondary tube or bedside test
3. Blood bag
4. Sharps disposal container for the safe disposal of used material.

Blood transfer with tube segment opener

General instructions

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING USE.


1. Observe the medical guidelines for transfusions, as well as your facility's procedures and guidelines.
2. Select a suitable secondary tube for the prescribed request.
3. Select a suitable tube for the application (depending on the version with \varnothing 12 or 11.5 mm) and press the tube segment opener onto it (fig. 1). Depending on the design, the tube segment opener is already pre-assembled with a tube. The tube used is also described below as a secondary tube.
4. Select blood bag.
5. Remove the tube segment opener from the carton.
6. Remove a barcode from the blood bag and apply it to the secondary tube.
7. Tear off the outer tube segment of the blood bag at the perforation line.
8. Hold the tube segment opener in one hand. Insert the tube segment at a slight angle (fig. 2a).
9. Push the tube segment forward (fig. 2b) without squeezing the tube (fig. 2a).
10. Now press the end of the tube once to release a small drop of blood (fig. 3).
11. Avoid dripping by decreasing pressure before detaching the tube segment from the tube segment opener (fig. 4).
12. Pull the tube segment opener out of the secondary tube and dispose of it in a sharps disposal container (fig. 5).

Tube segment opener for bedside test cards

General instructions

1. Observe the medical guidelines for transfusions, as well as your facility's procedures and guidelines.
2. The bedside test must be performed by physician.
3. Prepare bedside test.
4. Prepare a blood bag.
5. Remove the tube segment opener from the carton.
6. Take a barcode from the blood bag and affix it to the bedside test card.
7. Tear off the outer tube segment of the blood bag at the perforation line.
8. Hold the tube segment opener in one hand. Insert the tube segment at a slight angle (fig. 1a).
9. Push the tube segment forward (fig. 1b) without squeezing the tube.
10. Pierce the test card membrane with the needle end of the tube segment opener (fig. 2).
11. To dispense a drop of blood into the test card chamber, gently squeeze the tube segment opener.
12. Avoid unintended dripping by decreasing pressure before detaching the tube segment from the tube segment opener (fig. 3).
13. Pull the tube segment opener out of the secondary tube and dispose of it in a sharps disposal container (fig. 4).

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated tube segment openers must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended / updated

DIN 13097-4; Medical needles - Part 4: Ground types, requirements and testing

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: risk of contamination
	Keep away from sunlight
	Store in a dry place
	Do not use if package is damaged
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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