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SARSTEDT Thrombo Plus®

DE	Gebrauchshinweis – SARSTEDT Thrombo Plus®	2 - 4
EN	Instructions for use – SARSTEDT Thrombo Plus®	5 - 7
BG	Инструкции за употреба – SARSTEDT Thrombo Plus®	8 - 10
CS	Návod k obsluze – SARSTEDT Thrombo Plus®	11 - 13
DA	Brugsanvisning – SARSTEDT Thrombo Plus®	14 - 16
EL	Οδηγίες χρήσεως – SARSTEDT Thrombo Plus®	17 - 19
ES	Instrucciones de uso – SARSTEDT Thrombo Plus®	20 - 22
ET	Kasutusjuhend – SARSTEDT Thrombo Plus®	23 - 25
FR	Mode d'emploi – SARSTEDT Thrombo Plus®	26 - 28
HR	Uputa za upotrebu – SARSTEDT Thrombo Plus®	29 - 31
HU	Használati utasítás – SARSTEDT Thrombo Plus®	32 - 34
IT	Istruzioni d'uso – Thrombo Plus® SARSTEDT	35 - 37
KO	사용 설명서 – SARSTEDT Thrombo Plus®	38 - 40
LT	Naudojimo instrukcijos – SARSTEDT Thrombo Plus®	41 - 43
LV	Lietošanas norādes – SARSTEDT Thrombo Plus®	44 - 46
NL	Gebruiksaanwijzing – SARSTEDT Thrombo Plus®	47 - 49
NO	Bruksanvisning – SARSTEDT Thrombo Plus®	50 - 52
PL	Instrukcja obsługi – SARSTEDT Thrombo Plus®	53 - 55
PT	Instruções de utilização – Thrombo Plus® SARSTEDT	56 - 58
RO	Instrucțiuni de utilizare – SARSTEDT Thrombo Plus®	59 - 61
RU	Инструкция по применению – SARSTEDT Thrombo Plus®	62 - 64
SK	Návod na Použitie – SARSTEDT Thrombo Plus®	65 - 67
SL	Navodila za uporabo – Thrombo Plus® SARSTEDT	68 - 70
SV	Bruksanvisning – SARSTEDT Thrombo Plus®	71 - 73
TH	คำแนะนำในการใช้งาน – หลอดเก็บตัวอย่าง Thrombo Plus® ของ SARSTEDT	74 - 76
TR	Kullanım için talimatlar – SARSTEDT Thrombo Plus®	77 - 79
ZH	使用说明 – 莎斯特 (SARSTEDT) Thrombo Plus®	80 - 82

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Intended use

Thrombo Plus® is a prepared sample tube that is used, following the addition of EDTA blood (20 µl, 50 µl or 100 µl), for manually counting platelets according to common haematological practice. Thrombo Plus® is used for collecting capillary or venous EDTA blood and processing EDTA blood samples for *in vitro* diagnostic tests. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

Thrombo Plus® consists of a plastic tube with a round bottom, a red push cap with an integrated applicator and a 2.5 mg Hg salt/ml solution, pH 2.5. The nominal volume is 2 ml.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from blood and possible exposure to potentially infectious specimen material and transmissible pathogens.
2. Handle all biological specimens and sharp blood collection utensils according to the guidelines and procedures in your facility. In case of direct contact with biological specimens or of a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects for capillary blood collection in appropriate disposal containers.
4. The product is intended for single use. Dispose of all sharp objects (needles, lancets) for blood collection in suitable containers for hazardous biological waste.
5. Under- or overfilling the prepared sample tube leads to an incorrect ratio of blood to preparation/additive and can give rise to false analysis results.
6. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.

Storage

Store the product at room temperature.

Limitations

If blood samples are stored in the prepared sample tube, the stability of the analytes should be evaluated by the relevant laboratory or taken from the instructions for use of the analysis device manufacturer.


Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW BLOOD.

Equipment required for specimen collection:

1. All required Thrombo Plus® sample tubes and 20 µl, 50 µl or 100 µl of blood from an EDTA-coated End-to-End capillary or EDTA blood from a venous blood collection tube, e.g. S-Monovette® EDTA.
2. Gloves, gown, goggles or other suitable protective gear to protect from bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. Capillary holder, art. no. 95.1048, recommended for handling End-to-End capillaries.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

NOTE: Follow the regulations of your facility for the order of collection.

- 1
 - a Pipette 20 µl, 50 µl or 100 µl of EDTA blood from a venous EDTA blood collection tube or
 - b hold a 20 µl, 50 µl or 100 µl EDTA-coated End-to-End capillary horizontally or tilted slightly to collect the drops of blood. Stop collecting blood when the End-to-End capillary is completely filled with blood with no air bubbles.
- 2
 - a Unscrew the Thrombo Plus® sample tube.
 - a Put 20 µl, 50 µl or 100 µl of EDTA blood pipetted from a venous EDTA blood collection tube into the Thrombo Plus® sample tube or
 - b place the filled End-to-End capillary in the Thrombo Plus® sample tube.

The sample is haemolysed immediately after the blood is added.

Important: Do not put the EDTA blood into the Thrombo Plus® sample tube until at least 2-3 minutes after collecting the sample.
- 3
 - a Close the Thrombo Plus® sample tube and invert it carefully so that any End-to-End capillary used is emptied.
- 4
 - a The haemocytometer (e.g. from Thoma, Neubauer or others) is filled using the applicator integrated in the stopper.

Then the platelets are counted manually according to common haematological practice and following the instructions for use of the haemocytometer.

Note: For counting with normal values, a quantity of 20 µl of blood (factor 5,000) is sufficient. If low platelet counts are anticipated, larger quantities of blood, for example 50 µl (factor 2,000) or 100 µl (factor 1,000) can be used.

Calculation

Counting result x factor = platelet count / µl blood

Solution	Blood	Factor
2 ml	20 µl	5.000
2 ml	50 µl	2.000
2 ml	100 µl	1.000

Normal value: 150,000 - 400,000/µl blood

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled capillary blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

*CLSI (Clinical and Laboratory Standards Institute)

Hazard statements

H302 Harmful if swallowed.
H412 Harmful to aquatic life with long lasting effects.

Precautionary statements

P264 Wash hands thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
P273 Avoid release to the environment.
P301+P312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P330 Rinse mouth.
P501 Dispose of contents/container according to the regulations.

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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