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Intended use


Micro sample tubes for blood sugar testing from haemolysate are used with End-to-End capillaries and a haemolysing solution. The End-to-End capillary is used for capillary blood collection. Micro sample tubes are used for holding an End-to-End capillary filled with capillary blood and for the haemolysis and processing of blood samples in clinical laboratories. The products are intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

Micro sample tubes consist of a plastic vessel with a conical bottom, with or without a skirt depending on the version, and a colour-coded soft cap with loop, with and without a writable lid, various haemolysing solutions and a separate End-to-End capillary prepared with K2 EDTA.

Colour codes of the Micro sample tube soft caps:

Cap colour / pre-dosed with	Dosage/capillary	Vessel type
Yellow	Haemolysing solution (SARSTEDT) 1,000 µl/20 µl	1
Colourless	Haemolysing solution (DiaSys) 1,000 µl/20 µl	1
Colourless	Haemolysing solution (Roche Gen.2) 500 µl/20 µl	2
Yellow	Haemolysing solution (e.g. ESAT) 1,000 µl/20 µl	2

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures of your facility. In case of direct contact with biological specimens or of a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects (needles) for blood collection in suitable sharps containers.
4. Under- or overfilling the Micro sample tubes leads to an incorrect ratio of blood to preparation/additive and can result in false analysis results.
5. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.

Storage

Store the products at room temperature.

Limitations

1. If blood samples are to be stored in Micro sample tubes, the stability of the glucose should be evaluated by the relevant laboratory or taken from the instructions for use of the analysis device manufacturer.
2. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.

Collecting and handling specimens





READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW BLOOD.

Materials needed for collecting a specimen/blood:

1. All required Micro sample tubes and the corresponding K2 EDTA-coated End-to-End capillaries.
2. Gloves, gown, goggles or other suitable protective gear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet or incision lancet.
5. Disinfectant to clean the collection site (follow guidelines of the facility for sterile specimen collection to prepare the collection site).
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION AND WORK WITH BLOOD SAMPLES.

- 1 Clean the puncture site with disinfectant. Puncture the skin in a way that ensures good blood flow. Discard the first drop. Then hold the End-to-End capillary horizontally with one end in the middle of the drop of blood and completely fill the End-to-End capillary without air bubbles. Use a swab to wipe any blood off the outside of the capillary with a swab.
- 2 Open the Micro sample tube.
- 3 Put the End-to-End capillary filled with blood into the Micro sample tube.
- 4 Close the Micro sample tube and swirl it gently.

Note:

Use Micro sample tubes for testing blood sugar in accordance with the instructions for use of the test manufacturer and the Standard Operating Procedure (SOP) of your facility. The Micro sample tubes for testing blood sugar can be used for automated analysis. For use in devices, follow the respective device manufacturer's instructions.

For the stability of the glucose, follow the instructions for use of the test reagents / analysis device manufacturer.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled Micro sample tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

CLSI* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
 CLSI* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
 *CLSI (Clinical and Laboratory Standards Institute)

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CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
CLSI* H21 "Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays", Approved Guideline.
WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
*CLSI (Clinical and Laboratory Standards Institute)

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CLSI* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

CLSI* H21 " Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays" , Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002" .

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