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Intended use

The Micro sample tubes, brilliant cresyl blue are used for collecting, transporting and processing EDTA blood samples in the clinical laboratory. The product is intended for use in a professional environment by laboratory personnel for reticulocyte counting.

Product description

The Micro sample tube, brilliant cresyl blue consists of a plastic tube with a conical bottom, blue screw cap and a 1% brilliant cresyl blue solution. The nominal volume is 100 µl.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures of your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects (needles) for blood collection in suitable sharps containers.
4. If blood is collected from an intravenous (IV) line, it must be ensured that the line is properly flushed in accordance with the facility (=cleaned of the IV solution) before beginning blood collection. Properly flushing the line prevents faulty laboratory results.
5. Under- or overfilling the Micro sample tubes leads to an incorrect ratio of blood to the preparation/additive and can result in false analysis results.
6. Do not use the product after the expiry date. The expiry date is the last day of the month and year as indicated.

Storage

Store the product at room temperature.

Transport

The product is a primary receptacle according to the ADR (Packing instruction P650) and IATA regulations.

Limitations

1. If blood samples are to be stored in Micro sample tubes, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.

Specimen collection and handling

READ ALL OF THIS DOCUMENT BEFORE YOU START HANDLING THE PRODUCT.


Preparation for handling the Micro sample tube, brilliant cresyl blue, and material needed:

1. Select the Micro sample tube, brilliant cresyl blue incl. 100 µl from an EDTA-coated End-to-End capillary or 100 µl of EDTA blood from a reaction tube for the required sample.
2. Gloves, gown, goggles or other suitable protective gear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
5. Dry, clean swabs.
6. Plasters.
7. Sharps disposal container for the safe disposal of used material.

NOTE: Follow the regulations of your facility for the order of collection.

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION OR WORK WITH BLOOD SAMPLES.

- 1 a) Pipette 100 µl EDTA blood from an EDTA tube.
b) Hold a 100 µl EDTA-coated End-to-End capillary horizontally or tilted slightly to collect the drops of blood. Stop collecting blood when the End-to-End capillary is completely filled with blood with no air bubbles.
- 2 Screw open the Micro sample tube, brilliant cresyl blue.
a) Transfer Pipette 100 µl from an EDTA tube into the Micro sample tube, brilliant cresyl blue.
b) Transfer filled End-to-End capillary into the Micro sample tube, brilliant cresyl blue.
- 3 Close the Micro sample tube, brilliant cresyl blue and invert.
- 4 Wait for a reaction period of 15-30 minutes.
- 5 Invert the Micro sample tube, brilliant cresyl blue again.
- 6 Carry out the smear.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled Micro sample tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline .
 CLSI* H20 "Reference Leukocyte (WBC) Differential Count (Proportional) and Evaluation of Instrumental Methods", Approved Standard.
 *CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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