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EN Instructions for Use
Pressure filtration system
Filtropur L

Order No.	Description	Pore size	Sales unit
83.3944	Filtropur L 0.2 S with integral prefilter inlet: hose nipple (stepped hose nipple with 6 - 12 mm outer diameter) outlet: hose nipple (stepped hose nipple with 6 - 12 mm outer diameter)	0.2 µm	50/case 10/inner packaging individually wrapped sterile
83.3944.001	Filtropur L 0.2 LS with integral prefilter inlet: Luer-Lock outlet: hose nipple (stepped hose nipple with 6 - 12 mm outer diameter)	0.2 µm	50/case 10/inner packaging individually wrapped sterile
83.1850	Membrane pump for Filtropur L 0.2		1/case

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Pressure filtration system Filtropur L

Filtropur L 0.2 S and L 0.2 LS are sterile, individually wrapped filtration units for single laboratory use only and must not be applied to filter liquid for in-vivo use.

A polyethersulfone membrane with non-detectable non-specific protein binding characteristics and a glass fibre prefilter (which is 100 % free of binding agents) are sealed in a polycarbonate housing.

Eight symmetrically aligned ventilation holes enable **automatic ventilation** of the Filtropur L unit. A hydrophobic PTFE membrane protects the openings against the penetration of water.

Automatic ventilation guarantees **permanent use of the whole filtration area**.

Filtropur L pressure filtration units are available in the following versions:

- Connector: inlet hose nipple or luer lock
- With integral glass fibre prefilter
- 20 cm² filtration area and automatic ventilation
- Connector: outlet protected by sterile filling bell
- Polycarbonate housing with polyethersulfone membrane filter
- The pressure filtration unit is marked with "IN" = inlet and "OUT" = outlet. The sterile bell is positioned at the outlet.

Technical specifications

Housing material	polycarbonate
Filter material	polyethersulfone
	glass fibre prefilter
	(100 % free of binding agents)
Filter diameter	62 mm
Filter area	20 cm ²
Bubble Point	>= 3.2 bar
Hold-up volume before Bubble Point	approx. 1.5 ml
Housing burst pressure	> 5 bar
Flow direction	unidirectional
Max. recommended inlet pressure	3 bar
Range of filtration	100 ml - 10l
Liquid supply	peristaltic or pressure pump (diaphragm pump)

Instructions for Use

- The liquid to be filtered can be supplied via a pressure tank or peristaltic pump. Specifically for this range of application we offer a self-regulating diaphragm pump with electronic speed control and a set of connecting hoses from silicone.
- Connect the hose from your pressure source to the inlet side of the Filtropur L (IN) without taking the sterile side out of the packaging.
- Remove the Filtropur L unit from the packaging and mount the unit above the filtrate receptacle using a stand.
- Start filtration by slowly increasing the applied pressure. The entrapped air escapes through the PTFE membrane by automatic venting. As soon as the polyethersulfone membrane is saturated you can increase the inlet pressure to its maximum (set to 2.5 bar!).
- After filtration, the integrity of the filter can be tested as follows: Connect the filter to a pressure source and adjust the inlet side pressure to a constant 3.0 bar. Hold the Filtropur L unit under water. Air bubbles will appear from the automatic vents but none should at the outlet. This testing method is quite sensitive and detects the smallest defects with suitable accuracy.

Notes:

Avoid unnecessary air intake during the filtration by ensuring, for example, that the inlet end of the intake tubing is kept below the surface of the liquid being filtered, otherwise a pocket of air may build up above the membrane, resulting in partial or complete "blocking" of the membrane. Should this occur, however, pump all of your solution out of the system and then continue filtration as described in point 4.

Chemical compatibility

Filtropur L is not designed for organic solvents, alcohols, acetone, surface-active substances (emulsifying agents) and strong acids and bases.