


Salivette® VD

DE: Für den Arzt und das Labor · EN: For the doctor and the lab · BG: За лекаря и лабораторията · CS: Pro lékaře a laboratoře · DA: Til lægen og til laboratoriet · EL: Για το γιατρό και το εργαστήριο · ES: Para el médico y el laboratorio · ET: Arstile ja laboratooriumile · FR: Pour le médecin et le laboratoire · HR: Za liječnika i laboratorij · HU: Az orvosok és laboratóriumok részére · IT: Per il medico e il laboratorio · JP: 医師、および研究室用 · KO: 의사 및 검사실 사용 용도

DE: Gebrauchshinweis · EN: Instructions for Use · BG: Инструкции за употреба · CS: Návod k obsluze · DA: Brugsanvisning · EL: Οδηγίες χρήσεως · ES: Instrucciones de uso · ET: Kasutusjuhend · FR: Mode d'emploi · HR: Uputa za upotrebu · HU: Használati utasítás · IT: Istruzioni d'uso · JP: 取扱説明 · KO: 사용 설명서

5 DE: 2 min ohne Kaubewegung EN: 2 min without chewing BG: 2 мин. без дъвчене CS: 2 min bez žvýkání DA: 2 min. uden at tygge EL: 2 λεπτά χωρίς μάσημα ES: 2 min sin movimientos ET: 2 minutit ilma närimiseta FR: 2 min sans mastiquer HR: 2 min bez žvakanja HU: 2 percig rágomozgás nélkül IT: 2 min senza masticare JP: 咀嚼運動なしで2分間 KO: 씹지않고 2분

10 DE: Rückstände EN: Residues BG: Остатъци CS: Zbytky DA: Rester EL: Υπολείμματα ES: Residuos ET: Proov FR: Échantillon HR: Ostaci HU: Maradványok IT: Residui JP: 残留物 KO: 잔물

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EN: Collection and storage of saliva samples for screening for viral material, e.g. nucleic acids of SARS-CoV-2. The saliva sample should be taken at least 60 minutes after a meal (liquid/solid food intake) or after medications have been taken orally in order to avoid the saliva being contaminated by potential interferents. The swab should be kept in the mouth (e.g. in the cheek) for 2 minutes without chewing. The amount of saliva collected should be at least 500 µl. If an extremely small amount of saliva is produced, leave the swab in the mouth for longer. Saliva sample: The storage duration and temperature will depend on the shelf-life for the parameters to be tested. Therefore, and since the onset of bacterial growth in saliva can be expected after just a few hours at room temperature, it is recommended that the Salivette® be sealed tightly after sampling and either centrifuged and analysed within 4 hours or placed immediately in a refrigerator before being processed further. According to packaging instruction P650, shipping container (article No. 78.898) with absorbent liner plus screw cap (article No. 65.679) and a shipping box/bag are recommended for the shipment of Salivette® samples. The cotton roll (swab) is a natural product. Discoloration and darker inclusions or areas are natural and unavoidable.

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