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Intended use

The S-Sedivette® needles and adapters are used together as a system for venous blood collection. They are used for collecting, transporting and processing venous blood for manual *in vitro* diagnostic determination of the erythrocyte sedimentation rate (ESR*). The measurement results are obtained according to the Westergren method.

The product is intended for use in a professional environment by qualified medical and laboratory personnel.

*ESR = Erythrocyte sedimentation rate

Product description

The S-Sedivette® blood collection system consists of a plastic container, a plunger, a plunger rod and a colour-coded plastic screw cap with an integrated membrane developed especially for connecting it to an S-Monovette® (safety) needle or (Safety-) Multify® needle, and a citrate buffer preparation (0.105 mol/l pH 5.5 trisodium citrate/citric acid solution), equivalent to a 3.2% trisodium citrate/citric acid buffer solution in a mixed ratio of 1:4 (1 part citrate + 4 parts blood) and a glass sphere as a mixing aid.

The volume and permissible tolerances of the additive and the ratio of blood to additive comply with the requirements and recommendations of the international DIN EN ISO 6710 standard "Single-use containers for venous blood specimen collection" and the Clinical and Laboratory Standards Institute's Approved Standards (CLSI). The S-Sedivette® is internally sterile.

S-Sedivette® screw cap colour codes:

Additive	Letter code	Cap colour DIN EN ISO 6710	Cap colour based on BS 4851*
S-Sedivette®	4NC	black	violet

*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, revoked.

SARSTEDT S-Sedivette® Blood Collection System Citrate 4NC

The S-Sedivette® contains the anticoagulant trisodium citrate and is used for collecting whole blood. Trisodium citrate is present as a 0.105 mol/l pH 5.5 trisodium citrate/citric acid buffer solution (3.13% trisodium citrate solution; often rounded up to 3.2%) and amounts to 20% of the nominal volume of an S-Sedivette®. The mixture ratio of citrate to blood is 1:4 – 1 part by volume citrate and 4 parts by volume blood. Correct filling is absolutely necessary for the analysis.

This citrate whole blood is used as specimen material for routine tests, determining the erythrocyte sedimentation rate.

Anticoagulation occurs through calcium ion complex formation by citrate.

Please observe the NOTES under sample collection and handling for the S-Sedivette® and butterfly needle combination.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from blood and possible exposure to potentially infectious specimen material and transmissible pathogens.
2. Treat all biological specimens and sharp blood collection utensils (needles, adapters) according to the guidelines and procedures in your institution. Seek a doctor in the event of direct contact with biological samples or a needle-stick injury, as HIV, HCV, HBV and other infectious diseases can be transmitted in this way. Use safety needles/Safety-Multify® needles with an integrated needle guard. Always follow the safety guidelines and procedures of your facility.
3. The product is intended for single use. Dispose of all sharp objects (needles, adapters) for blood collection in appropriate hazardous biological waste containers.
4. If blood is collected through an intravenous (IV) line, it must be ensured that the line is properly flushed in accordance with facility guidelines (= cleaned of the IV solution) before beginning blood collection. Properly flushing the line prevents faulty analysis results.
5. Underfilling the S-Sedivette® will result in an incorrect ratio of blood to preparation/additive and can produce incorrect analysis results.
6. The blood collected and processed with the S-Sedivette® is not intended to be re-injected into the human body.
7. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.
8. Liquid preservatives and anticoagulants are clear. Do not use if the liquids are cloudy or precipitation is present.

storage

Store the product at room temperature until use.

Limitations

1. It should be stored at room temperature. The ESR should be determined within the first 4 hours after collecting the blood. The specimen can be stored in a refrigerator (4 °C) for a longer period (max. 24 hours). The specimen must then be brought to room temperature before being used.
2. Homogenise the blood sample by swirling it carefully (inverting it 5 times) before measuring the erythrocyte sedimentation rate. The air bubble/glass bead must always travel to the other end. Do not shake quickly to homogenise the sample.
3. Measuring requires an ambient temperature of 18–25 °C and protection against vibrations, draughts and direct sunlight must be ensured before measuring.
4. The measuring principle complies with CLSI H2-A5.
5. The measurement accuracy is +/- 1 mm of the measurement path.
6. The mixing ratio of 1:4 has a direct effect on the analysis result and must be adhered to.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING VENOUS BLOOD COLLECTION.

Equipment required for specimen collection

1. All necessary S-Monovette® blood collection systems (including the S-Sedivette®).
2. (Safety) needles or (Safety-)Multify® needles
3. Gloves, gown, goggles or other suitable protective gear to protect from bloodborne pathogens or potentially infectious materials.
4. Labels for patient identification.
5. Disinfectant to clean the collection site (follow the guidelines of the facility for sterile specimen collection for preparing the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean, disposable swabs.
7. A tourniquet.
8. Plaster or bandage.
9. Sharps disposal container for the safe disposal of used material.

Recommended order of collection:

1. Blood culture vials.
2. Tubes with citrate or clotting activator (with/without gel, Serum CAT/Serum Gel CAT).
3. Tubes with clotting activator or citrate.
4. Tubes with heparin with/without gel.
5. Tubes with EDTA with/without gel.
6. Tubes with glycolysis inhibitors.
7. Tubes with other additives.

NOTE: When starting collection with a citrate tube or another tube with a liquid additive in combination with a butterfly needle (e.g. (Safety-)Multify® needle), collection with an empty tube (e.g. S-Monovette® neutral Z) first is recommended to avoid underfilling of the tube due to filling of the tubing of the butterfly needle.

This can ensure the correct filling of the tube and thus a correct mixing ratio (liquid to blood).

NOTE: Follow the regulations of your facility for the order of collection.

Preventing reflux

S-Sedivettes contain chemical additives. To prevent reflux, the following should be noted:

1. Use only the materials recommended in the instructions for use.
2. Position the patient so that there is easy access to the vein, and the arm or other puncture site on the patient is directed downwards.
3. Hold the S-Sedivette® with the screw cap pointing upwards and, if using a (Safety-) Multify® needle, also below the level of the arm/puncture site.
4. It must be ensured that the contents of the S-Sedivette® do not touch the screw cap or the end of the needle during the venipuncture.
5. Loosen the tourniquet while the blood is flowing into the S-Sedivette®.

Blood collection

Collection techniques

A distinction is made between two collection techniques: the aspiration technique and the vacuum technique.

The aspiration technique enables gentle blood collection with a continuous, slow blood flow, by controlled pulling on the plunger rod. In this way, the blood flow can be adapted immediately to all vein conditions and circumstances as well as reducing haemolysis.

The vacuum technique using a pre-evacuated tube allows blood to be collected according to the negative pressure principle with a continuous rapid blood flow. This technique is therefore ideal for good vein conditions and easy collection conditions.

Handling of blood collection: see also the handling video for the S-Monovette® in the aspiration technique or vacuum technique with the safety needle or Safety-Multifly® needle: [www.sarstedt.com/en/download/videos/Diagnostics/Venous blood](http://www.sarstedt.com/en/download/videos/Diagnostics/Venous%20blood)

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING BLOOD COLLECTION.


1. Select a suitable S-Sedivette® for the required sample.
2. Select a site for the venipuncture. Please follow the guidelines of the institution for preparing the venipuncture site.

Blood collection with the aspiration technique

Important! Use the S-Sedivette® only for blood collection, not for injection.

Use only with these needles (S-Monovette® (safety) needles, or (Safety-)Multifly® needles/adapters for the S-Monovette®).

Use with a safety needle is shown here:

1. Immediately prior to the venipuncture, put the S-Sedivette® into the integrated holder of the safety needle **A** and secure it by **twisting clockwise slightly B**.
Apply a venous tourniquet if necessary (max. 1 minute). Prepare the puncture site with a suitable disinfectant. Do not touch the puncture site after disinfection. Position the patient so that there is easy access to the vein, and the arm or other puncture site on the patient is directed downwards.

2. Pull the needle's protective cap to remove it. Puncture the vein, slowly retract the plunger rod and loosen the tourniquet while the blood is flowing into the S-Sedivette® **C**. Wait until the flow of blood stops to allow correct filling.

3. Move the S-Sedivette® by **twisting it anticlockwise slightly D** and then pull (disconnect) it from the safety needle **E**. The safety needle remains in the vein.

For multiple specimens – as described above – the subsequent S-Monovettes are connected to the safety needle and additional specimens are taken.

Completing blood collection: First disconnect the S-Sedivette® **D** + **E**, then withdraw the safety needle from the vein. Using one hand, either activate the needle guard on a stable, level surface – until the needle snaps into the needle guard with a **perceptible and audible "click"** – or by pressing the index finger against the lower end of the needle guard, see also instructions for use for the S-Monovette® (safety) needle or (Safety-)Multifly® needle.

Immediately mix the blood sample thoroughly. Invert the S-Sedivette® slowly five times. The air bubble/glass bead must always travel to the other end. Do not shake quickly.

For transport and centrifugation, retract the plunger rod until the plunger snaps audibly into place **F** (plunger clicks in the bottom of the S-Sedivette®) and break off the plunger rod **G**.


1. Press a clean, dry swab against the puncture site until bleeding stops.
2. As soon as coagulation has set in, you may apply a dressing if desired.
3. Dispose of the safety needle in a suitable container for hazardous biological waste.

Blood collection with the vacuum technique

Important! Use the S-Sedivette® only for blood collection, not for injection.

Use only with these needles (S-Monovette® (safety) needles, or (Safety-)Multifly® needles/adapters for the S-Monovette®).

Use with a safety needle is shown here:

1 To begin blood collection gently, SARSTEDT recommends collecting the first S-Monovette® using the aspiration technique, see VENIPUNCTURE TECHNIQUE Aspiration technique **A – E**. Then continue with the vacuum technique.

Before collecting blood, use the plunger rod to pull the plunger to the bottom of the S-Sedivette® and allow it to snap into place audibly **A**. Then the plunger rod is broken off **B**.

2 Slide the prepared S-Sedivette® into the integrated holder of the S-Monovette® safety needle and secure it by **twisting it clockwise slightly C + D**.

3 Wait until the blood stops flowing to allow correct filling.

4 Move the S-Sedivette® by **twisting it anticlockwise slightly E** and then pull (disconnect) it from the safety needle **F**. The safety needle remains in the vein.

For multiple specimens – as described above – the subsequent S-Monovettes are connected to the safety needle and additional specimens are taken.

Completing blood collection:

First disconnect the S-Sedivette® **E + F**, then withdraw the safety needle from the vein. Using one hand, activate the needle guard either on a stable, level surface – until the needle snaps into the needle guard with a **perceptible and audible "click"**– or by pressing the index finger against the lower end of the needle guard.

Immediately mix the blood sample thoroughly. Invert the S-Sedivette® slowly five times. The air bubble/glass bead must always travel to the other end.


Do not shake quickly.

1. Press a clean, dry swab against the puncture site until bleeding stops.
2. As soon as coagulation has set in, you may apply a dressing if desired.
3. Dispose of the safety needle in a suitable container for hazardous biological waste.

Note for barcode labelling:

Place the label on the bottom of the S-Sedivette®.

Labels which are attached too high or askew will disrupt the measurement as they will cover the blood column in the measurement area.

Performing the ESR measurement

General information

Homogenise the blood sample by swirling it carefully before measuring the erythrocyte sedimentation rate. We recommend using the Sarmix® rotation mixer.

Automated determination:

The Sediplus® is available for automated determination. Observe the Instructions for Use of the Sediplus® S 2000 NX (Art. no.: 90.1092). Online instructions for use: Sediplus® S 2000 NX; GB 557 at: www.sarstedt.com/en/download/instructions-for-use

Manual determination:

For manually determining the ESR, use the SARSTEDT ESR rack with graduations Art. no.: 90.1090. Observe the Instructions for Use of the SARSTEDT ESR rack. Online instructions for use: SARSTEDT Manual erythrocyte sedimentation – ESR rack & sedimentation pipettes; GB 562 at www.sarstedt.com/en/download/instructions-for-use

1. Invert the sample carefully.
2. The measurement system (Sediplus® or the SARSTEDT sedimentation rack) must stand on a horizontal base.
To perform the measurement, place the filled S-Sedivette® vertically in the measurement system (Sediplus® or the SARSTEDT sedimentation rack). Observe the corresponding Instructions for Use.
3. ESR determination should be completed within 4 hours of blood collection.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. Dispose of potentially contaminated material according to the regulations and guidelines of the facility.

Product-specific standards and guidelines as amended/updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

CLSI* H2-A5 2011: Procedures for the Erythrocyte Sedimentation Rate Test, 5th Edition Item Order Code: H02ASE.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

Standard for containers for venous blood collection:

DIN EN ISO 6710: Single-use containers for human venous blood specimen collection

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ICSH recommendations for measurement of erythrocyte sedimentation rate. J. Clin. Pathol. 1993 46:198-203

Comparative Study SARSTEDT Blood Sedimentation Systems (see at: www.sarstedt.com/produkte/diagnostik/blutsenkung/blutentnahmesysteme/ Literatur)

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

Order no.

Batch number

Use by

CE mark

In vitro diagnostic devices

Follow the instructions for use

Risk of contamination if reused

Keep away from sunlight

Store in a dry place

Manufacturer

Country of manufacture

Sterile fluid path

Single sterile barrier system

Single sterile barrier system with protective packaging outside

Do not resterilise

Technical modifications reserved

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
 CLSI* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.
 CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
 CLSI* H2-A5 2011: Procedures for the Erythrocyte Sedimentation Rate Test, 5th Edition Item Order Code: H02ASE.
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

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- CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
- CLSI* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.
- CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
- CLSI* H2-A5 2011: Procedures for the Erythrocyte Sedimentation Rate Test, 5th Edition Item Order Code: H02ASE.
- WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

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CLSI* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection" , Approved Standard.

CLSI* GP41 " Collection of Diagnostic Venous Blood Specimens" , Approved Standard.

CLSI* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

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WHO/DIL/LAB/99.1 Rev02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002" .

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