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SARSTEDT Blutgaskapillaren

DE	Gebrauchshinweis – SARSTEDT Blutgaskapillaren	1 - 4
EN	Instructions for use – SARSTEDT blood gas capillaries	5 - 7
BG	Инструкции за употреба – капилари за кръвни газове SARSTEDT	8 - 10
CS	Návod k obsluze – kapiláry pro stanovení krevních plynů SARSTEDT	11 - 13
DA	Brugsanvisning – SARSTEDT Kapillærør til blodgasanalyse	14 - 16
EL	Οδηγίες χρήσεως – Τριχοειδή αερίων αίματος SARSTEDT	17 - 19
ES	Instrucciones de uso – Capilares para gasometría SARSTEDT	20 - 22
ET	Kasutusjuhend – veregaasikapillaarid SARSTEDT	23 - 25
FR	Mode d'emploi – Tubes capillaires pour gaz du sang SARSTEDT	26 - 28
HR	Uputa za upotrebu – SARSTEDT kapilare za analizu plinova u krvi	29 - 31
HU	Használati utasítás – SARSTEDT vérgáz kapillárisok	32 - 34
IT	Istruzioni d'uso – Capillari per emogasanalisi SARSTEDT	35 - 37
KO	사용 설명서 – SARSTEDT 혈액 가스 모세관	38 - 40
LT	Naudojimo instrukcijos – SARSTEDT kraujo dujų kapiliarai	41 - 43
LV	Lietošanas norādes – SARSTEDT asins gāzu kapilārcaurulītes	44 - 46
NL	Gebruiksaanwijzing – SARSTEDT-bloedgascapillairen	47 - 49
NO	Bruksanvisning – SARSTEDT blodgaskapillarer	50 - 52
PL	Instrukcja obsługi – kapilary do gazometrii SARSTEDT	53 - 55
PT	Instruções de utilização – Capilares de gases sanguíneos SARSTEDT	56 - 58
RO	Instrucțiuni de utilizare – capilare pentru gaze sanguine SARSTEDT	59 - 61
RU	Инструкция по применению – Капилляры для газов крови SARSTEDT	62 - 64
SK	Návod na Použitie – SARSTEDT kapiláry na analýzu krvných plynov	65 - 67
SL	Navodila za uporabo – kapilare za plinsko analizo krvi SARSTEDT	68 - 70
SV	Bruksanvisning – SARSTEDT blodgaskapillärer	71 - 73
TH	คำแนะนำในการใช้งาน – หลอดเก็บตัวอย่างเลือดจากเส้นเลือดฝอยของ SARSTEDT	74 - 76
TR	Kullanım için talimatlar – SARSTEDT Kan gazı için kapilerler	77 - 79
ZH	使用说明 – SARSTEDT 血气毛细管	80 - 82

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Intended use

The blood gas capillaries are coated with calcium balanced heparin and are intended for taking a capillary blood sample for point of care testing (POCT) and an immediate determination of blood gases, pH or electrolytes in a blood gas analysis device. The blood gas capillaries are used with the corresponding quick-release caps and mixing wires. Using the mixing magnet, the blood sample is optimally mixed with the anticoagulant following the drawing of capillary blood. The products are intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The blood gas capillaries are made of unbreakable PET plastic and are prepared with calcium balanced heparin. They are available in volumes ranging from 100 to 175 µl.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures in your institution. In case of direct contact with biological specimens, consult a doctor due to the risk of transmission of HIV, HCV, HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects (needles) for blood collection in appropriate hazardous biological waste containers.
4. Underfilling the blood gas capillary leads to an incorrect ratio of blood to preparation/additive and can result in faulty analysis results.
5. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.
6. Please note the safety data sheet at <https://www.sarstedt.com/en/download/safety-data-sheets/>.

Storage

Store the products at room temperature.

Collecting and handling specimens


READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW CAPILLARY BLOOD.

Preparation for drawing capillary blood and material needed:

1. Blood gas capillary.
2. Gloves, gown, goggles or other suitable protective gear to protect against pathogens or potentially infectious materials transmitted by blood.
3. Labels for specimen identification.
4. (Safety) lancet.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples).
6. Dry, clean swabs.
7. Plaster.
8. Sharps disposal container for the safe disposal of used material.

Capillary technique with the blood gas capillary

- 1 Prepare the blood gas capillary by placing an appropriate quick-release cap on one end.
- 2 Now put a mixing wire into the blood gas capillary.
- 3 Clean the puncture site with disinfectant. Puncture the skin in a way that ensures good blood flow. Discard the first drop. Take the cap off the capillary. Then hold the capillary horizontally with one end in the middle of the drop of blood and completely fill the capillary without air bubbles.
- 4 Close both ends of the capillary tightly with the caps.
- 5 Move the mixing wire back and forth 10-15 times over the entire length of the capillary with the aid of the mixing magnet to mix the blood with the anticoagulant.
- 6 Mix the sample again immediately before the analysis. Then position the mixing wire at the end of the capillary.
- 7 Remove both quick-release caps.
- 8 Let the blood sample be suctioned off by the device. The blood gas analysis should be completed within 15 minutes after drawing capillary blood.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and burned.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

US standards / ISO standards (as amended/updated)

Product-specific standards:

CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI C46 Blood Gas and pH Analysis and Related Measurements; Approved Standard

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use before
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Keep away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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