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## Intended use


The manual ESR\* racks and the SARSTEDT sedimentation pipettes are used for the manual, *in vitro* diagnostic determination of the erythrocyte sedimentation rate (ESR) from anticoagulated citrate whole blood (mixing ratio 1:4; 1 part citrate + 4 parts blood). The measurement results are obtained according to the Westergren method. Depending on the blood collection system used, venous or capillary blood is used. Capillary blood is used with the Microvette® CB 200 ESR and venous blood with the S-Sedivette®, the S-Monovette® ESR or the ESR screw-cap tube. The products are intended for use in a professional environment by qualified medical and laboratory personnel.

\*ESR = Erythrocyte sedimentation rate

## Product description

ESR racks (also called sedimentation racks) are made of plastic and, depending on the design, are matched to the blood collection tubes S-Sedivette®, S-Monovette® ESR, the screw-cap tube ESR or Microvette® CB 200 ESR. They are available with and without a graduated scale on the rear wall and have 10 measuring stations each. The SARSTEDT sedimentation pipettes made of plastic, with and without graduation, are used as accessories. The S-Monovette® ESR or the screw-cap tube ESR in combination with the SARSTEDT sedimentation pipette with graduation are suited for use with the SARSTEDT sedimentation rack without rear wall. Depending on the version, the SARSTEDT sedimentation racks with scaled rear wall fit the S-Sedivette® as well as the combination: S-Monovette® ESR or to the screw-cap tube ESR with the sedimentation pipette with 0 mark (no graduation) or the Microvette® CB 200 ESR with its sedimentation capillary.

## Product overview

Type	Art. no.	Name	For use with
<b>A</b>	90.1090	SARSTEDT ESR rack with graduated scale	Suitable for the S-Sedivette®
<b>B</b>	90.1060	SARSTEDT ESR rack with graduated rear wall	Suitable for the S-Monovette® ESR or the screw-cap tube ESR in combination with sedimentation pipette with 0 mark (variant E)
<b>C</b>	90.1060.062	SARSTEDT ESR rack without rear wall	Suitable for the S-Monovette® ESR or the screw-cap tube ESR in combination with sedimentation pipette with printed graduated scale (variant F)
<b>D</b>	90.1091	SARSTEDT ESR rack for the Microvette® CB 200 ESR	Suitable for the Microvette® CB 200 ESR with sedimentation capillary
<b>E</b>	86.1996	SARSTEDT sedimentation pipettes with 0 mark	Suitable for SARSTEDT ESR rack with graduated rear wall (variant B)
<b>F</b>	86.1996.062	SARSTEDT sedimentation pipettes with printed graduated scale	Suitable for SARSTEDT ESR racks without rear wall (variant C)

## Safety information and warnings

1. Follow the directions for handling.
2. General precautions: Use gloves and general personal protective equipment to protect yourself from possible exposure to potentially infectious specimen material and transmitted pathogens.
3. Handle all biological specimens and collection utensils according to the guidelines and procedures in your facility. In the event of contact with biological samples, consult a doctor, as this can transmit infectious diseases. Always follow the safety guidelines and procedures of your facility.
4. The sedimentation pipettes are intended for single use. Dispose of all products in the hazardous biological waste containers.
5. The sedimentation pipettes must not be used after the expiry date. The expiry date of the product is the last day of the month and year indicated.

## Storage

Store the products at room temperature.

## Physical-chemical principle

Due to the higher density of erythrocytes compared to plasma, they slowly form sediment in anticoagulated blood. However, the sedimentation of erythrocytes is influenced by various factors that can accelerate or slow down the rate of sedimentation. In addition to the difference in density, the most important influencing variables include erythrocyte aggregation and its electrical charge. The negative charge of the erythrocytes causes them to repel one another and the sedimentation takes place only slowly. A changed composition of the plasma proteins, especially the increase in the concentration of acute phase proteins in inflammatory processes, reduces the zeta potential (repulsive effect of the negative charge) and thus promotes erythrocyte aggregation, which leads to a faster sedimentation behaviour of the blood cells and consequently to higher measured values for ESR\*. In contrast, the increase in concentration of negatively charged plasma proteins leads to slower sedimentation of the erythrocytes and thus to lower ESR measured values.

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Reference values should be set at the local level in accordance with the recommendations for reference values. In view of the progressive increase in ESR with age, separate values should be set for each decade of adult life for men and women. Several other variables influence the ESR and can therefore affect the reference values, e.g. haemoglobin levels, medications, menstrual cycle, pregnancy and smoking. The table "Reference values for the Westergren ESR method" lists reference values for the ESR, which can be used as a guide for establishing local reference values.

## Reference values for the Westergren ESR method\*\*

Age (years)	Mean ESR			
	Male	Female	upper limit of normal	
			Male	Female
18–30	3.1	5.1	< 7.1	< 10.7
31–40	3.4	5.6	< 7.8	< 11.0
41–50	4.6	6.2	< 10.6	< 13.2
51–60	5.6	9.4	< 12.2	< 18.6
60–70	5.6	9.4	< 12.7	< 20.2
> 70	5.6	10.1	< 30	< 35

A normal ESR does not exclude non-inflammatory organ diseases, organ dysfunction and malignant tumours.\*\*\*

\*\*CLSI H2-A5 2011: Procedures for Erythrocyte Sedimentation Rate Test.

\*\*\*Bruhn et al.; Lehrbuch der Labormedizin; (Textbook of Laboratory Medicine) Schattauer, 1999; ISBN 3-7945-1856-X

## Limitations

- The blood samples should be stored at room temperature. The ESR should be determined within the first 4 hours after collecting the blood.  
The specimen can be stored in a refrigerator (4 °C) for a longer period (max. 24 hours). The specimen must then be brought to room temperature before being used.
- The blood sample must be thoroughly homogenised by careful swirling before measuring the erythrocyte sedimentation rate.
- The measurement requires an ambient temperature of 18–25 °C and must be protected from vibrations, draughts and direct sunlight or in accordance with country and laboratory-specific specifications.
- The measuring principle complies with CLSI H2-A5.
- The measuring accuracy is +/- 1 mm of the measuring path.
- The SARSTEDT ESR racks can be cleaned with ethanol/isopropanol-based disinfectants.
- The mixing ratio of 1:4 has a direct effect on the analysis result and must be adhered to.
- The rack may be used only if the scale cover is not deformed and the scale is clearly legible.

## Collecting and handling specimens

READ ALL OF THIS DOCUMENT BEFORE STARTING THE ESR MEASUREMENT.

### Equipment required for the ESR determination

1. Sample tube filled with blood sample (S-Sedivette®, S-Monovette® ESR, screw-cap tube ESR or Microvette® CB 200 ESR)\*.
2. SARSTEDT ESR rack matching the sample tube and, if necessary, a suitable SARSTEDT sedimentation pipette.
3. Gloves, gowns, eye protection or other suitable protective clothing to protect against pathogens or potentially infectious materials transmitted by sample material.
4. Alarm clock or stopwatch.

\*For blood collection and sample transport, the instructions for use of the respective sample tube (S-Sedivette®, S-Monovette® ESR, screw-cap tube ESR or Microvette® CB 200 ESR) must be observed.

### How to determine the ESR

General information


The different variants of the SARSTEDT ESR racks and SARSTEDT sedimentation pipettes are described below:

- A: Use of SARSTEDT ESR rack, Art. no.: 90.1090
- B: Use of SARSTEDT ESR rack, Art. no.: 90.1060 with SARSTEDT sedimentation pipettes, Art. no.: 86.1996
- C: Use of ESR rack, Art. no.: 90.1060.062 with SARSTEDT sedimentation pipettes, Art. no.: 86.1996.062
- D: Use of SARSTEDT ESR rack, Art. no.: 90.1091

#### A: How to use the SARSTEDT ESR rack Art. no.: 90.1090

The SARSTEDT ESR rack (Art. no.: 90.1090) is suitable only for use with the S-Sedivette® (Art. no.: 06.1690.xxx).



1. Carefully mix the sample by inverting it (Fig. 1).
2. The SARSTEDT ESR rack must stand on a horizontal base. To perform the measurement, place the filled S-Sedivette® vertically in the SARSTEDT ESR rack. The zero point of each sample must be set with the knurled screw (arrow) (Fig. 2).
3. Start the measurement by setting and starting an alarm clock or a stopwatch for your measurement times. Optionally set to 30 minutes (1st measured value) and 1 hour (2nd measured value) or 1 hour (1st measured value) and 2 hours (2nd measured value).
4. After the time +/- 1 minute has elapsed, the distance in millimetres between the bottom of the plasma layer and the upper end of the column with the sedimented erythrocytes is read off from the graduated scale on the rear wall. Make sure that no leukocytes (buffy coat) are contained in the erythrocyte column. Record the numeric value.

**B: How to use the SARSTEDT ESR rack Art. no.: 90.1060 with SARSTEDT sedimentation pipettes Art. no.: 86.1996**

The SARSTEDT ESR rack (Art. no.: 90.1060) is suitable only for use with the S-Monovette® ESR (Art. no.: 05.1079.xxx) or the screw tube ESR (Art. no.: 47.410.xxx) and the SARSTEDT sedimentation pipette (Art. no.: 86.1996).

- Carefully mix the samples by inverting them (Fig. 1).
- Open the S-Monovette® ESR/the screw-cap tube ESR by turning the screw cap anticlockwise (Fig. 2).
- Insert the sedimentation pipette into the S-Monovette® ESR/screw tube ESR. Carefully push down the sedimentation pipette. The sedimentation pipette will then fill with blood. Screw in the plunger sleeve of the sedimentation pipette with the S-Monovette® ESR/the screw-cap tube ESR by turning it clockwise (Fig. 3).
- The sedimentation pipette is properly screwed together when the filling mark on the sedimentation pipette has exactly reached the 0 mark (Fig. 4).
- For the measurement, place the S-Monovette® ESR/the screw-cap tube ESR with the now screwed in sedimentation pipette vertically on the SARSTEDT ESR rack with graduated rear wall. The SARSTEDT ESR rack must stand on a horizontal base. The zero point of each sample must be checked and, if necessary, adjusted by rotating the plunger sleeve of the sedimentation pipette (Fig. 5).
- Start the measurement by setting and starting an alarm clock or a stopwatch for your measurement times. Optionally set to 30 minutes (1st measured value) and 1 hour (2nd measured value) or 1 hour (1st measured value) and 2 hours (2nd measured value).
- After the time +/- 1 minute has elapsed, the distance in millimetres between the bottom of the plasma layer and the upper end of the column with the sedimented erythrocytes is read off from the graduated scale on the rear wall. Make sure that no leukocytes (buffy coat) are contained in the erythrocyte column. Record the numeric value.


**C: How to use the SARSTEDT ESR rack Art. no.: 90.1060.062 with SARSTEDT sedimentation pipettes Art. no.: 86.1996.062**

The SARSTEDT ESR rack (Art. no.: 90.1060.062) is suitable only for use with the S-Monovette® ESR (Art. no.: 05.1079.xxx) or the screw tube ESR (Art. no.: 47.410.xxx) and the SARSTEDT sedimentation pipette (Art. no.: 86.1996.062).


- Carefully mix the samples inverting them (Fig. 1).
- Open the S-Monovette® ESR/the screw-cap tube ESR by turning the screw cap anticlockwise (Fig. 2).
- Insert the sedimentation pipette into the S-Monovette® ESR/screw tube ESR. Carefully push down the sedimentation pipette. The sedimentation pipette will then fill with blood. Screw in the plunger sleeve of the sedimentation pipette with the S-Monovette® ESR/the screw-cap tube ESR by turning it clockwise (Fig. 3).
- The sedimentation pipette is properly screwed together when the filling mark on the sedimentation pipette has exactly reached the 0 mark (Fig. 4).
- For the measurement, place the S-Monovette® ESR/the screw-cap tube ESR with the now screwed in sedimentation pipette vertically on the SARSTEDT ESR rack without rear wall. The SARSTEDT ESR rack must stand on a horizontal base. The zero point of each sample must be checked and if needed adjusted by rotating the plunger sleeve of the sedimentation pipette (Fig. 5).
- Start the measurement by setting and starting an alarm clock or a stopwatch for your measurement times. Optionally set to 30 minutes (1st measured value) and 1 hour (2nd measured value) or 1 hour (1st measured value) and 2 hours (2nd measured value).
- After the time +/- 1 minute has elapsed, the distance in millimetres between the bottom of the plasma disc and the upper end of the column with the sedimented erythrocytes is read off from the graduated scale printed on the sedimentation pipette. Make sure that no leukocytes (buffy coat) are contained in the erythrocyte column. Record the numeric value.

**D: How to use the SARSTEDT ESR rack Art. no.: 90.1091**

The ESR rack (Art. no.: 90.1091) is suitable only for use with the Microvette® CB 200 ESR (Art. no.: 18.1325.xxx).

- Carefully mix the sample by inverting it (Fig. 1).
- Open the Microvette® CB 200 ESR by pulling the stopper out upwards. Insert the sedimentation capillary into the Microvette® CB 200 ESR. Carefully push down the sedimentation capillary. The sedimentation capillary will then fill with blood (Fig. 2).
- The SARSTEDT ESR rack must stand on a horizontal base. To perform the measurement, place the Microvette® CB 200 ESR with sedimentation capillary vertically in the SARSTEDT ESR rack. The zero point of each sample must be exactly set with the knurled screw (arrow) (Fig. 3).
- Start the measurement by setting and starting an alarm clock or a stopwatch for your measurement times. Optionally set to 30 minutes (1st measured value) and 1 hour (2nd measured value) or 1 hour (1st measured value) and 2 hours (2nd measured value).
- After the time +/- 1 minute has elapsed, the distance in millimetres between the bottom of the plasma layer and the upper end of the column with the sedimented erythrocytes is read off from the graduated scale on the rear wall. Make sure that no leukocytes (buffy coat) are contained in the erythrocyte column. Record the numeric value.

**Disposal**

- The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
- Disposable gloves prevent the risk of infection.
- Contaminated or filled SARSTEDT sedimentation pipettes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
- Dispose of potentially contaminated material according to the regulations and guidelines of the facility.

## Product-specific standards and guidelines as amended / updated

CLSI\*\* H2-A5 2011: Procedures for Erythrocyte Sedimentation Rate Test.

Thomas, L.: Labor und Diagnose. Digital edition of Labor & Diagnose 2020; Chapter 19.3.2 Blutkörperchensenkungs-Reaktion (BSR) [Erythrocyte Sedimentation Reaction (ESR)] ([www.labor-und-diagnose-2020.de/k19.html#\\_idTextAnchor8638](http://www.labor-und-diagnose-2020.de/k19.html#_idTextAnchor8638))

ICSH recommendations for measurement of erythrocyte sedimentation rate. J. Clin. Pathol. 1993 46:198-203

\*\*CLSI (Clinical and Laboratory Standards Institute)

## Key for symbols and labels:

Article number

Batch number

Use by

CE mark

*In vitro* diagnostic device

Follow the instructions for use

If reused: Risk of contamination

Keep away from sunlight

Store in a dry place

Manufacturer

Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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