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## SARSTEDT CSF Röhren

<b>DE</b>	Gebrauchshinweis – SARSTEDT CSF Röhren	2 - 4
<b>EN</b>	Instructions for use – SARSTEDT CSF tubes	5 - 7
<b>BG</b>	Инструкции за употреба – SARSTEDT ЦСТ епруветки	8 - 10
<b>CS</b>	Návod k obsluze – Zkumavky CSF SARSTEDT	11 - 13
<b>DA</b>	Brugsanvisning – SARSTEDT CSF-rør	14 - 16
<b>EL</b>	Οδηγίες χρήσεως – Σωληνάρια SARSTEDT ENY	17 - 19
<b>ES</b>	Instrucciones de uso – Tubos CSF SARSTEDT	20 - 22
<b>ET</b>	Kasutusjuhend – SARSTEDTI CSF-katsutid	23 - 25
<b>FR</b>	Mode d'emploi – Tubes CSF SARSTEDT	26 - 28
<b>HR</b>	Uputa za upotrebu – SARSTEDT CSF epruvete	29 - 31
<b>HU</b>	Használati utasítás – SARSTEDT CSF csövek	32 - 34
<b>IT</b>	Istruzioni d'uso – Provette LCR SARSTEDT	35 - 37
<b>KO</b>	사용 설명서 – SARSTEDT CSF 튜브	38 - 40
<b>LT</b>	Naudojimo instrukcijos – SARSTEDT CSF mėgintuvėliai	41 - 43
<b>LV</b>	Lietošanas norādes – SARSTEDT CSF stobriņi	44 - 46
<b>NL</b>	Gebruiksaanwijzing – SARSTEDT CSF-buisjes	47 - 49
<b>NO</b>	Bruksanvisning – SARSTEDT CSF rør	50 - 52
<b>PL</b>	Instrukcja obsługi – próbówki na PMR SARSTEDT	53 - 55
<b>PT</b>	Instruções de utilização – Tubos de LCR SARSTEDT	56 - 58
<b>RO</b>	Instrucțiuni de utilizare – Tuburi pentru LCR SARSTEDT	59 - 61
<b>RU</b>	Инструкция по применению – SARSTEDT Пробирки для спинномозговой жидкости (СМЖ)	62 - 64
<b>SK</b>	Návod na Použitie – skúmavky SARSTEDT CSF	65 - 67
<b>SL</b>	Navodila za uporabo – epruvete CSF SARSTEDT	68 - 70
<b>SV</b>	Bruksanvisning – SARSTEDT CSF-rör	71 - 73
<b>TH</b>	คำแนะนำในการใช้งาน – หลอดเก็บตัวอย่าง CSF ของ SARSTEDT	74 - 76
<b>TR</b>	Kullanım için talimatlar – SARSTEDT BOS Tüpleri	77 - 79
<b>ZH</b>	使用说明 – SARSTEDT CSF 脑脊液样品管	80 - 82

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
Lewczuk et al, Effect of sample collection tubes on cerebrospinal fluid concentrations of tau proteins and amyloid beta peptides, ClinChem 52 (2):332-4 2006.

Perret-Liaudet et al, Risk of Alzheimer's disease biological misdiagnosis linked to cerebrospinal collection tubes, J Alzheimers Dis 31(1):13-20 2012

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
## Intended use

The CSF tubes are specimen receptacles for CSF (cerebrospinal fluid) that have been designed for the collection, transport, processing and storage of CSF for *in vitro* diagnostics. The specimen material in the tubes can be sent for automatic or manual specimen processing. The specifications of the device manufacturer must be observed. The products are intended for use in a professional environment by healthcare professionals and qualified medical and laboratory personnel.

## Product description

CSF tubes consist of a clear receptacle with a graduated scale and a round bottom, and are offered in different variants. The versions vary in capacity and dimensions, and are available with different imprints depending on the variant. The products are sealed with a screw cap and individually sterile packed.

## Product overview

## Tubes and caps

Type	Description
A	CSF set
B	CSF tube 5 ml
C	CSF tube 10 ml

## Safety information and warnings

**Do not use the product if the packaging is damaged.**

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from possible exposure to potentially infectious sample material and transmissible pathogens, as well as hazardous substances.
2. Handle all biological specimens and sharp collection utensils according to the established guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Check each individual package for damage before use and do not use if packaging is damaged.
4. Dispose of all sharp/pointed objects (needles) for sample collection in suitable containers for hazardous biological waste.
5. The products are intended for single use. Dispose of all products in the hazardous biological waste containers.
6. Do not use the product after the expiry date. The expiry date is the last day of the month and year as indicated.

## Storage

Store the products at room temperature.

## Transport

The CSF tubes (5 ml & 10 ml) together with an ADR/IATA (packaging instruction P650) compliant secondary container, incl. suction insert and type-tested outer packaging, comply with the ADR P650/IATA regulations.

## Limitations

1. The storage time and temperature of a filled CSF tube depend on the expiry date of the parameters to be examined. Clarify the time within which the analysis should be carried out with the responsible laboratory or refer to the analysis manufacturer's instructions for use.
2. Check the CSF specimen quality. Note the information provided by the test manufacturer regarding suitability for haemolysis and precipitates for clinical tests.

## Collecting and handling specimens

**READ THIS DOCUMENT IN FULL BEFORE STARTING SAMPLE COLLECTION.**

### Equipment needed to collect a specimen:

1. Appropriate protective clothing and general equipment for CSF collection in line with the established hygiene guidelines at your facility.
2. All required specimen receptacles as per the request.
3. Labels for specimen identification.
4. Sharps disposal container for the safe disposal of used equipment.

**NOTE: Check whether specific pre-analytical handling processes are necessary according to the test manufacturer and follow these steps.**

**COLLECT THE SPECIMEN IN LINE WITH THE ESTABLISHED GUIDELINES AND PROCEDURES IN YOUR FACILITY.**

## Handling of the CSF tube

1. To deposit the sterile CSF tube, open the blister pack over a sterile cloth. The transparent film must be facing upwards.
2. Open the CSF tube by turning the screw cap anticlockwise.
3. Place the screw cap on a sterile cloth with the opening facing upwards.
4. Allow the cerebrospinal fluid to drain directly into the CSF tube.
5. Put the cap on and carefully close the CSF tube.
6. Place the CSF tube in an upright position according to the intended storage conditions for the specimen material.
7. When using the CSF set, the sequence can be ticked on the printed boxes in the case of multiple collections.
8. Affix the specimen identification.
9. Store the tube closed according to the storage conditions specified for your sample material.

## Centrifugation

**IMPORTANT!**  
**SARSTEDT CSF tubes (5 ml & 10 ml) are designed for a maximum of 4,000 x g (RCF - relative centrifugal force – g-force). Use only suitable tube holders or inserts. Centrifugation of a cracked CSF tube, or centrifugation using excessive centrifugal force can cause breakage of the CSF tube, releasing potentially hazardous substances.**

Centrifuge inserts should be selected according to the size of the CSF tube used. The relative centrifugal force has the following ratio to the selected rpm:

$$\text{RCF (g-force)} = 11.2 \times r \times (\text{rpm}/1,000)^2$$

"RCF": "Relative centrifugal force",

"rpm": "Revolutions per minute",

"r" [in cm]: Centrifuge radius from the centre of the centrifuge to the bottom of the CSF tube.

CSF tubes should be centrifuged in accordance with the centrifugation conditions set out below. If other conditions are used, they must be validated by the respective user.

Ensure that the CSF tube fits snugly in the centrifuge adapter. Tubes that protrude beyond the adapter can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

**CAUTION! Do not remove a broken CSF tube by hand.**

*Information on disinfecting the centrifuge can be found in the instructions for use of the centrifuge.*

## Centrifugation recommendation

The recommended RCF (g-force) for the centrifugation of CSF specimens with precipitates is 2,000 x g for 10 minutes.

## Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled CSF tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. Dispose of potentially contaminated material according to the regulations and guidelines of your facility.

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
Lewczuk et al., Effect of sample collection tubes on cerebrospinal fluid concentrations of tau proteins and amyloid beta peptides, ClinChem 52 (2):332-4 2006.

Perret-Liaudet et al., Risk of Alzheimer's disease biological misdiagnosis linked to cerebrospinal collection tubes, J Alzheimers Dis 31(1):13-20 2012

Leitlinien für Diagnostik und Therapie in der Neurologie, Lumbalpunktion und Liquordiagnostik, Entwicklungsstufe [Guidelines for Diagnostics and Therapy in Neurology, Lumbar Puncture and CSF Diagnostics, Development Stage]: S1, AWMF (Association of Scientific Medical Societies in Germany) registration number: 030/141

## Key for symbols and labels:

	Article number		Do not use if package is damaged
	Batch number		Sterilisation by irradiation
	Use by		Single sterile barrier system with external protective packaging
	CE mark		Do not resterilise
	<i>In vitro</i> diagnostic device		
	Follow the instructions for use		
	If reused: Risk of contamination		
	Keep away from sunlight		
	Store in a dry place		
	Manufacturer		
	Country of manufacture		

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority. SARSTEDT AG & Co. KG

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