

The S-Monovette® RNA Exact stabiliser solution maintains sample integrity over multiple freeze/thaw cycles without loss of quality

Schulz D.², Reiter J.¹, Richter S.¹, Kämper M.¹

¹ SARSTEDT AG & Co. KG, Research & Development, Nümbrecht, Germany

² SARSTEDT AG & Co. KG, Marketing & Product Management, Nümbrecht, Germany



Introduction

RNA molecules are considered highly unstable, especially since the sample can easily become contaminated with ubiquitous RNases. Even in living organisms, RNA can easily be destroyed by oxidative reactions, spontaneous hydrolysis and enzymes. It is therefore important to stabilise samples that are taken for subsequent RNA analysis against degradation.

Due to the instability of RNA molecules, it is common practice to freeze both stabilised and non-stabilised samples at -80°C . It is also customary to thaw the samples in order to take an aliquot for sample isolation and then re-freeze them. Even when samples have been gathered together for a test series, it is customary to freeze the samples for a short period and then thaw them for joint processing. However, this standard procedure runs up against the fact that repeated freeze/thaw cycles are detrimental to sample integrity and ought to be avoided.

We have therefore investigated the performance of the S-Monovette® RNA Exact stabiliser solution in repeated freeze/thaw cycles.

Materials & Methods

3 samples were taken from each of 8 test subjects using the S-Monovette® RNA Exact in accordance with the instructions for use. The reference specimen (T0) was frozen at -80°C after having been stored for at most 1 hour. The other 2 samples were frozen, and afterwards thawed during the day and re-frozen at -80°C at night over a period of 3 days. These blood samples have therefore been frozen and thawed out a total of 4 times prior to isolation.

The RNA was isolated using the NucleoSpin® RNA Blood Midi Kit (REF 740210.20) manufactured by MACHEREY-NAGEL GmbH & Co. KG, in accordance with the manufacturer's instructions. The concentration of the isolated RNA was measured and the quality control performed photometrically using the Eppendorf SE BioPhotometer®. The sample integrity (RIN) was determined using an Agilent 2100 Bioanalyzer.

Result

The purity of the samples subjected to multiple freeze/thaw cycles, was shown using OD 260/280 ratios, which is a measure of contamination with proteins. OD 260/280 ratios were comparable to the controls, at an average of 2.1. The mean of the photometrically determined concentration (yield) for the samples subjected to multiple freeze/thaw cycles was 72 ng/ μl . The RNA concentration of the control samples was 18.6% higher on average. Nevertheless, the sample integrity (RIN) of the samples subjected to multiple freeze/thaw cycles is unchanged compared to the control samples, at an average of 8.8.

	RIN	STDEV
Thawed 4x	8.8	0.2
Control	8.8	0.29

The sample integrity was also checked using quantitative RT PCR and shown using so-called $\Delta\Delta\text{CT}$ evaluation. Here, having taken into account isolation fluctuations above the average of the CT value of three housekeeping genes (Δ), the change in the samples subjected to multiple freeze/thaw cycles compared to the control samples is shown ($\Delta\Delta$). A value here that is as close to 0 as possible indicates that no changes have occurred in the initial concentration of this mRNA.




Figure:

Depiction of the conserved expression levels of genes after freeze-thaw cycles. RNA from S-Monovette® RNA Exact stabilized blood was either isolated directly after blood collection or after four freeze-thaw-cycles. cDNA was synthesized and expression levels of selected genes were compared to the levels obtained from the directly isolated sample with the ddCt-method.

Discussion

The experiment shows that repeatedly thawing and re-freezing samples stabilised with RNA Exact can result in a slight loss in yield. Thawing and re-freezing 4 times resulted in an average loss of 18.6% in the 8 tested donors.

The RNA integrity (RIN) and the ratios of the various mRNAs to one another (tested on 5 different genes) in samples stabilised with RNA Exact were not affected by the 4 freeze/thaw cycles performed.

The S-Monovette® RNA Exact stabiliser solution is able to ensure sample integrity over multiple freeze/thaw cycles without loss of quality. As it was already previously established that long-term storage of isolated RNA, even at -80°C, leads to a loss in quality, S-Monovette® RNA Exact now offers an alternative solution to sample storage, making it possible to process aliquots from the sample and re-freeze the rest, for example. This goes some way to making day-to-day work much easier.

For additional product or technical information, please e-mail us at marketing@sarstedt.com or visit www.sarstedt.com.

SARSTEDT AG & Co. KG

Sarstedtstraße 1
D-51588 Nümbrecht
Germany

www.sarstedt.com
info@sarstedt.com

Worldwide Support Offices

Australia

Tel: +61 8 8349 6555
Fax: +61 8 8349 6882
info.au@sarstedt.com

Austria

Tel: +43 2236 616 82
Fax: +43 2236 620 93
info.at@sarstedt.com

Belgium

Tel: +32 3 541 7692
Fax: +32 3 541 8103
info.be@sarstedt.com

Brazil

Tel: +55 11 4152 2233
Fax: +55 11 4152 3198
info.br@sarstedt.com

Canada

Tel: +1 514 328 6614
Tol free: 1 888 727 7833
Fax: +1 514 328 9391
info.ca@sarstedt.com

China

Tel: +86 21 5062 0181
Fax: +86 21 5058 0700
info.cn@sarstedt.com

Croatia

Tel: +385 95 36 77 030
Fax: +385 14 96 10 75
info.hr@sarstedt.com

Czech Republic

Tel: +420 281 021 491
Fax: +420 281 021 495
info.cz@sarstedt.com

Finland

Tel: +358 9 374 1044
Fax: +358 9 374 1176
info.fi@sarstedt.com

France

Tel: +33 3 84 31 95 95
Fax: +33 3 84 31 95 99
info.fr@sarstedt.com

Germany

Telefon +49 22 93 305-0
Telefax +49 22 93 305-3450
Service (0800) 0 83 305-0
info@sarstedt.com

Greece

Tel: +30 210 6038 274
Fax: +30 210 6038 276
info.gr@sarstedt.com

Hungary

Tel: +36 1 383 1216
Fax: +36 1 383 1213
info.hu@sarstedt.com

Ireland

Tel: +353 53 91 44922
Fax: +353 53 91 44998
info.ie@sarstedt.com

Italy

Tel: +39 045 8510 114
Fax: +39 045 8510 118
info.it@sarstedt.com

Japan

Tel: +81 3 5215 5400
Fax: +81 3 5215 6400
info.jp@sarstedt.com

Latvia

Tel: +371 6 731 0386
Fax: +371 6 704 0723
info.lv@sarstedt.com

Mexico

Tel: +52 55 8501 1577
Fax: +52 55 8501 1578
info.mx@sarstedt.com

Netherlands

Tel: +31 76 501 7550
Fax: +31 76 501 7626
info.nl@sarstedt.com

Norway

Tel: +47 64 856 820
Fax: +47 64 856 821
info.no@sarstedt.com

Poland

Tel: +48 22 722 0543
Fax: +48 22 722 0795
info.pl@sarstedt.com

Portugal

Tel: +351 21 915 6010
Fax: +351 21 915 6019
info.pt@sarstedt.com

Russia

Tel: +7 495 937 5228
Fax: +7 495 937 5228
info.ru@sarstedt.com

Slovakia

Tel: +421 2 682 45 933
Fax: +421 2 682 45 934
info.sk@sarstedt.com

Spain

Tel: +34 93 846 4103
Fax: +34 93 846 3978
info.es@sarstedt.com

Sweden

Tel: +46 42 19 84 50
Fax: +46 42 19 84 59
info.se@sarstedt.com

Switzerland

Tel: +41 81 750 1880
Fax: +41 81 750 1899
info.ch@sarstedt.com

Turkey

Tel: +90 216 290 18 65
Fax: +90 216 290 18 64
info.tr@sarstedt.com

United Arab Emirates

Tel.: +971 4 3888 080
Fax: +971 4 3888 282
info.ae@sarstedt.com

United Kingdom

Tel: +44 116 2359 023
Fax: +44 116 2366 099
info@sarstedt.co.uk

USA

Tel: +1 800 257 5101
Tel: +1 828 465 4000
Fax: +1 828 465 4003
customerservice@sarstedt.us