

# S-Monovette® ThromboExact

Specially designed for thrombocyte determination when pseudothrombocytopenia is suspected

**NEW!**



## For thrombocyte determination when pseudothrombocytopenia is suspected

The term pseudothrombocytopenia describes an erroneously low thrombocyte count. In contrast to a genuine reduction of blood platelets, pseudothrombocytopenia does not constitute a clinical picture but, instead, a mere preanalytical phenomenon that happens to occur when working with automated blood cell analysers<sup>1,2,3</sup>.

Generally, pseudothrombocytopenia is caused by thrombocyte aggregation, which makes precise counting of the thrombocytes impossible. Aggregation most frequently occurs when EDTA is used as an anticoagulant. However, patient-related aggregation can also occur with other additives used for anticoagulation, such as Heparin or Citrate.

Early detection of this phenomenon ensures that the diagnostic and therapeutic consequences of a thrombocytopenia misdiagnosis are avoided.

The following graph shows an example of multiple-intolerance reactions:




Fig. 1: Example of multiple-intolerance reactions (Rostock University Hospital)

According to special literature in this field, adding certain aminoglycoside antibiotics<sup>4,5</sup> or using a CTAD solution<sup>6,7</sup> may prevent in vitro aggregation of thrombocytes. Such inhibitors are, however, rarely used because they are typically expensive and alter the blood/additive dilution ratio.

The S-Monovette® ThromboExact provides a solution to the problems described above. The product is currently being tested at the universities of Munich and Rostock for the determination of an accurate thrombocyte count when EDTA- induced pseudothrombocytopenia and multiple-intolerance reactions occur.

Initial data compiled from these studies are illustrated as follows:




Fig. 2: Thrombocyte values from EDTA blood and S-Monovettes ThromboExact analysed immediately and 4 hours after sampling Munich University, Sysmex XE-2100)

| Preparation    | EDTA  | Citrate | ThromboExact |     |     |     |     |     |      |
|----------------|-------|---------|--------------|-----|-----|-----|-----|-----|------|
|                |       |         | immediately  | 1 h | 2 h | 4 h | 5 h | 7 h | 12 h |
| Elapsed Time   | 2-4 h | 2-4 h   |              |     |     |     |     |     |      |
| Patient No. 12 | 126   | 266     | 265          | 273 | 258 |     |     |     | 254  |
| Patient No. 13 | 57    | 172     | 279          | 279 | 288 | 284 |     |     | 261  |
| Patient No. 14 | 25    | 313     | 378          | 398 | 393 | 384 |     |     | 376  |
| Patient No. 15 | 37    |         | 142          | 124 |     |     |     |     |      |
| Patient No. 16 | 77    | 72      | 157          | 149 |     | 156 |     |     | 152  |
| Patient No. 17 | 78    | 213     | 222          | 281 |     | 281 | 278 | 267 |      |
| Patient No. 18 | 12    | 20      | 142          | 150 | 144 | 147 |     |     | 142  |
| Patient No. 19 | 116   | 356     | 362          |     |     |     |     |     |      |

Tab. 1: Thrombocyte count (in 1,000/µl) from blood treated with EDTA, Citrate, and S-Monovettes ThromboExact at defined times after sampling (Rostock University, Sysmex XE 2100).

The preliminary results obtained from these studies to date are as follows:

- The data compiled for the S-Monovette® ThromboExact confirm a case of genuine thrombocytopenia for patient No. 3.
- The period of time after blood collection during which an EDTA-induced pseudothrombocytopenia manifests varies among individual patients. This becomes apparent when comparing the thrombocyte count obtained for patients 1, 2, 6, and 7 with that for patients 8, 9, and 10.
- Patients 13, 16, and 18 show examples of multiple-intolerance reactions.
- The data obtained for patients 12, 13, 14, 16, and 18 reveal that the S-Monovette® ThromboExact provides a correct thrombocyte count 12 hours after blood collection.

### Conclusion:

*Thrombocytes can still be correctly determined 12 hours after blood collection.*

*Intolerance reactions to EDTA and/or Citrate can be ascertained.*

*Corrections due to the dilution ratio are no longer necessary.*

*The specific nature of the S-Monovette® ThromboExact ensures assignment of the sample to the correct laboratory.*

## References

- 1:** E. Gowland, H. E. M. Kay, J. C. Spillman und J. R. Williamson: Agglutination of platelets by a serum factor in the presence of EDTA *J. Clin. Path.* **1969**, 22, 460-464.
- 2:** T. Nilsson und B. Norberg: Thrombocytopenia and pseudothrombocytopenia: A clinical and laboratory problem *Scand. J. Haematol.* **1986**, 37, 341-346.
- 3:** F. Silvestri, L. Virgolini, C. Savignano, F. Zaja, M. Velisig und M. Baccarani: Incidence and Diagnosis of EDTA-Dependent Pseudothrombocytopenia in a Consecutive Outpatient Population Referred for Isolated Thrombocytopenia *Vox Sang* **1995**, 68, 35-39.
- 4:** S. Sakurai, I. Shiojima, T. Tanigawa und K. Nakahara: Aminoglycosides prevent and dissociate the aggregation of platelets in patients with EDTA-dependent pseudothrombocytopenia *Br. J. Haematol.* **1997**, 99, 817-823.
- 5:** A. J. P. F. Lombarts, J. J. Zijlstra, R. H. M. Peters, C. G. Thomasson und P. F. H. Franck: Accurate Platelet Counting in an Insidious Case of Pseudothrombocytopenia *Clin Chem Lab Med* **1999**, 37, 1063-1066.
- 6:** G. Contant, M. Goualt-Heilmann und J. L. Martinoli: Heparin inactivation during blood storage: its prevention by blood collection in citric acid, theophylline, adenosine, dipyridamole *Thromb. Res.* **1983**, 31, 365-374.
- 7:** O. Ohnuma, Y. Shirata und K. Miyazawa: Use of theophylline in the investigation of pseudothrombocytopenia induced by edetic acid *J. Clin. Pathol.* **1988**, 41, 915-917.

## Ordering Information

| Order No.   | Description               | Packaging                          |
|-------------|---------------------------|------------------------------------|
| 05.1168.001 | S-Monovette® ThromboExact | 50 pcs./inner box<br>500 pcs./case |

SARSTEDT AG & Co.  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone +49 2293 3050  
Fax +49 2293 305-122  
info@sarstedt.com  
www.sarstedt.com