




Gentle to animals

Good for animals,
designed for your practice.



Treating animals well



“We use the S-Monovette® for large dogs. This allows us to fill the sample tube easily and cleanly even with just one hand.”

As an innovator with decades of experience in human pre-analytics, SARSTEDT offers a number of solutions that have also been proven in veterinary medicine where the focus is on safety, user friendliness and animal welfare.

With products such as the S-Monovette®, we set safety standards for closed venous blood collection from pets such as dogs and cats. We also offer sample collection solutions that meet the specific requirements of your smallest patients such as small mammals, birds or reptiles.


Test our thoughtfully designed systems and see for yourself the many practical benefits for your day-to-day work.

Blood collection with SARSTEDT. Good for animals, designed for your practice.

Pre-analytics engineered by SARSTEDT

The pre-analytical phase is the most important for every examination. Even with sophisticated analytical equipment, the result is only as good as the quality of the sample.

With our expertise in pre-analytics, we bring you devices that support your daily work. The goal is to give you more time for the animals and owners that come to your practice.



“If I get only a few drops of blood, the Micro needle helps me get as much as possible from the collection – it prevents coagulation at the needle aperture.”



“Using the Multifly® needle is an enormous relief for me. When the puncture is made, I have enough time to fill the required sample tubes.”



“We call the Multivette® the ‘vacuum’. We just connect it to the Luer needle and can fill two or three tubes with blood while we are in the vein.”

- ✓ Safe
- ✓ Gentle
- ✓ Hygienic

The convenient solution for every animal – together for greater animal welfare.

Many different patients come to your veterinary practice, all with their individual needs. SARSTEDT gives you suitable solutions for pre-analytics in veterinary medicine.

- The S-Monovette® – safe, gentle and hygienic blood collection from larger pets
- The Multivette® – collect blood gently and easily from small pets using only venous pressure
- The Micro needle – even in especially sensitive animals, collect every drop reliably, e.g. with a Micro sample tube or a Microvette®
- Many other practical systems and single products for collecting blood, urine and fecal samples

More about Animal Health Care from SARSTEDT



animalhealthcare.sarstedt.com



Closed blood collection from larger pets




The S-Monovette® blood collection system allows you to perform **closed blood collection with two techniques**. The aspiration technique enables gentle collection adapted to the weaker blood flow of fragile veins. When blood flow is strong, the dual system of the S-Monovette® also allows the conventional vacuum technique for blood collection.

The sustainable blood collection system

The great advantage of blood collection with an S-Monovette® is that the blood can be taken directly into the correctly prepared sample tube. No extra disposable syringe is needed to transfer the blood sample to other tubes. This saves resources and makes blood collection safer and more hygienic than alternative collection options.

The more sustainable cardboard packaging also keeps disposal costs and waste volume low.

Make the trip to the veterinarian a positive experience for your patients and their owners



Test now
for free!



animalhealthcare.
sarstedt.com





"I went to the doc again today, but this time something was different: Everything was completely relaxed.

Before, I just got pricked in the leg and this red liquid spilled out everywhere. Today, everything went very quickly without any mess. That must be because of the doc's new equipment.

I hope it's like this every time now!"

Lucky, male boxer



The S-Monovette®

Safety needle guard for use with one hand

Integral adapter provides secure, tight connection

Visible check for successful vein puncture using the first drop of blood that enters the S-Monovette® or the tubing (with the Safety-Multifly® needle)

The S-Monovette® is available for blood volumes from 1.2–4.5 ml in a range of additives

Gentle collection using the aspiration technique



Note: the S-Monovette® can also be used with the Safety-Multifly® needle.

Relaxed blood collection from small pets




The Multivette® 600 is designed and prepared **for small blood volumes** of 600 µl. Blood is collected in a nearly closed condition via natural venous pressure. This makes collection especially gentle and easy.

The blood sample can be centrifuged directly in the Multivette® 600. This eliminates sample transfer and saves time.

The small inner diameter makes it particularly easy to pipette blood off after centrifugation. Additionally, the sample can be securely closed for sending to the laboratory.


The Multivette® 600



"Today my can openers put me in this box again and took me away. When they finally let me out, they first rubbed me behind my ears.*

Then something magical happened: They held this thing against my fur, the color changed from white to red, they rubbed me again and I got to go home. It wasn't bad at all."


Bella, Bengal cat
**In the cat world, can opener means "human"*



Specially designed inner tube is filled from venous pressure alone.

Range of preparations for 600 µl blood samples

Successful drip method in especially sensitive animals



"Hi guys! My mom was worried because I haven't been feeling well for a few days. My vet first drew some blood to find out why. She had a special needle for me that I could hardly feel. Actually, I didn't really notice anything.

Mom was also much more relaxed than usual. That was most certainly because it was over so quickly – she can't stand the sight of blood..."

Coco, female Chihuahua



The Micro needle with micro sample tubes or Microvette®

Especially in very small animals and **difficult vein conditions**, every drop of blood counts. The Micro needle ensures that every drop of blood flows into the sample tube and does not coagulate beforehand. The appropriate sample tube can be selected according to the expected sample volume.

Micro sample tubes are suitable for volumes of 1.3 ml. For smaller patients, the Microvette® in volumes 100 – 500 µl is the right choice.

Optimal flow properties

Minimal blood loss with no residual volume in the holder

Direct filling eliminates blood loss due to sample transfer

Flexibility resulting from a range of volumes and additives



Test now for free!



animalhealthcare.sarstedt.com



SARSTEDT

SARSTEDT, Inc.

1025 St. James Church Road
P.O. Box 468
Newton, NC 28658 - 0468

Tel: +1 828 465 4000
Toll free: +1 800 257 5101
Fax: +1 828 465 4003

Customer Service Germany
Phone +49 0800 0 83 305 0

customerservice.us@sarstedt.com
www.sarstedt.com



Further information and free samples
for testing:
animalhealthcare.sarstedt.com

If you have any questions:
We'd be happy to help!

Visit our website: www.sarstedt.com