# Instruction manual

flexiPERM® slide, micro12, disc, conA, conB

After removing the flexiPERM®, hydrophobicity remains in the zones surrounding each of the culture chambers. Subsequent individual staining steps of the different culture chambers on one slide is therefore easier to perform.

After proper cleaning and sterilization, the flexiPERM® may be used repeatedly.

## Cleaning

After cell cultivation, immediately rinse the flexiPERM® in flowing water accompanied by slight rubbing between thumb and forefinger. Rinse the flexiPERM® in distilled water for a few moments.

If heavily soiled, the flexiPERM® may also be cleaned in diluted solution of HCl, NaOH or KOH, and then stored in sterile distilled water.

Longterm exposure to concentrated acids and bases destroys flexiPERM®.

Store in a plastic bag or 70 % pharmaceutical ethanol.

## Sterilizing

a) Ethanol bath - Prior to use, place flexiPERM® in a 70 % pharmaceutical ethanol solution 20 min and then place it in sterile distilled water for 20 min to dissipate residual ethanol. Gently shake off excess distilled water and allow to air dry in an open cell culture dish under a laminar flow hood.

b) Autoclaving - Immerse flexiPERM® in a suitable, vented container (e.g. glass jar) of distilled water and autoclave it under normal conditions (121°C, 20 min). Following sterilization and cooling, store the flexiPERM® in the closed capped container. Prior to use, gently shake off excess distilled water and allow the flexiPERM® to air dry in an open petri dish under a laminar flow hood.


## Notes

flexiPERM® turns turbid (milky) by absorption of water or ethanol. When rendered turbid by water, the flexiPERM® may be used again after drying, see "Sterilizing b)". When rendered turbid by ethanol place in sterile distilled water for 20 min and dry it as described above.

flexiPERM® should not come into contact with disinfectants and detergents.

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
# Instruction manual

flexiPERM® slide, micro12, disc, conA, conB

flexiPERM® is a reusable cell culture insert made of silicone.

flexiPERM® slide (1) and flexiPERM® micro12 (2) can be used in combination with a microscope slide (26 mm x 76 mm) thus creating 8 or 12 compartments, respectively.

flexiPERM® disc (3) can be used to subdivide tissue culture dishes and lumox® dishes 50 in four smaller growth areas.

flexiPERM® conA (4) and conB (5) are single chamber versions. The flexiPERM® itself constitutes the walls of the chambers, while the slide or dish constitutes the bottom.

When used in combination with microscope slides flexiPERM® conA and conB facilitate examination by micromanipulation (e.g. plant and animal physiology).

In addition to the surface of slides or dishes, other substrates (e.g. glass, plastic or foil) may be used as growth surface.

flexiPERM® is hydrophobic, non-toxic to cells and is resistant to heat, cold and almost all laboratory chemicals.

flexiPERM® is stable against diluted acids and bases and also in a temperature range from -20°C to +125°C. Concerning solvent compatibility ask for technical data sheet.

## Instructions for use

For adherence the smooth, polished under side of the flexiPERM® must be gently pressed against the dry slide or dish. Only new, clean substrate should be used. It should be ensured that the polished side of the flexiPERM® is in good condition, as scratches could affect the adherence of the flexiPERM® to the substrate. Do not use mechanically treated slides (e.g. Superfrost). When coated slides are required, the coating must be accomplished after having adhered the flexiPERM® onto the slide.

For better adherence of the cells to glass slides, the glass surface may be slightly abraded by treating it with 0.1 N NaOH (0.4%), followed by neutralization with 0.1 N HCl (0.36%). Rinse the slides in sterile distilled water and dry them prior using for cell culture.

Isolated, leak-proof chambers are formed upon pressing flexiPERM® on the desired substrate. Adherence can be checked by viewing the flexiPERM®/substrate contact from below in reflected light.

flexiPERM® could be used for longterm experiments approximately two weeks without becoming detached.

Upon completion of cell cultivation, and prior to fixation, the flexiPERM® can be detached from the substrate by lifting flexiPERM® at the flap.