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## SARSTEDT Safety-Multifly® Kanülen/Multifly® Kanülen

|           |                                                                                                       |         |
|-----------|-------------------------------------------------------------------------------------------------------|---------|
| <b>DE</b> | Gebrauchshinweis – SARSTEDT Safety-Multifly® Kanülen/Multifly® Kanülen                                | 2 - 4   |
| <b>EN</b> | Instructions for use – SARSTEDT Safety-Multifly® Needles/Multifly® Needles                            | 5 - 7   |
| <b>BG</b> | Инструкции за употреба – SARSTEDT Safety-Multifly® Канюли/Multifly® Канюли                            | 8 - 10  |
| <b>CS</b> | Návod k obsluze – SARSTEDT Jehel Safety-Multifly®/Jehel Multifly®                                     | 11 - 13 |
| <b>DA</b> | Brugsanvisning – SARSTEDT Safety-Multifly® Kanyler/Multifly® Kanyler                                  | 14 - 16 |
| <b>EL</b> | Οδηγίες χρήσεως – SARSTEDT Βελόνες Safety-Multifly®/Βελόνες Multifly®                                 | 17 - 19 |
| <b>ES</b> | Instrucciones de uso – SARSTEDT Agujas de seguridad Multifly®/Agujas Multifly®                        | 20 - 22 |
| <b>ET</b> | Kasutusjuhend – SARSTEDT Safety-Multifly® Kanüülid/Multifly® Kanüülid                                 | 23 - 25 |
| <b>FI</b> | Käyttöohje – SARSTEDT Safety-Multifly®-Kanyylit/Multifly®-Kanyylit                                    | 26 - 28 |
| <b>FR</b> | Mode d'emploi – SARSTEDT Aiguilles Safety-Multifly®/Multifly®                                         | 29 - 31 |
| <b>HR</b> | Uputa za upotrebu – SARSTEDT Kanile Safety-Multifly®/Kanile Multifly®                                 | 32 - 34 |
| <b>HU</b> | Használati utasítás – SARSTEDT Safety-Multifly® Kanülök/Multifly® Kanülök                             | 35 - 37 |
| <b>IT</b> | Istruzioni d'uso – SARSTEDT Aghi Safety-Multifly®/Aghi Multifly®                                      | 38 - 40 |
| <b>KO</b> | 사용 설명서 – SARSTEDT Safety-Multifly® 바늘/Multifly® 바늘                                                    | 41 - 43 |
| <b>LT</b> | Naudojimo instrukcijos – SARSTEDT „Safety-Multifly®“ Kaniulės/„Multifly®“ Kaniulės                    | 44 - 46 |
| <b>LV</b> | Lietošanas norādes – SARSTEDT Safety-Multifly® Kanīles/Multifly® Kanīles                              | 47 - 49 |
| <b>NL</b> | Gebruiksaanwijzing – SARSTEDT Safety-Multifly®-Naalden/Multifly®-Naalden                              | 50 - 52 |
| <b>NO</b> | Bruksanvisning – SARSTEDT Safety-Multifly®-Kanyler/Multifly®-Kanyler                                  | 53 - 55 |
| <b>PL</b> | Instrukcja obsługi – SARSTEDT Igły Safety-Multifly®/Igły Multifly®                                    | 56 - 58 |
| <b>PT</b> | Instruções de utilização – SARSTEDT Agulhas Safety-Multifly®/Agulhas Multifly®                        | 59 - 61 |
| <b>RO</b> | Instrucțiuni de utilizare – Canule SARSTEDT Safety-Multifly®/Canule SARSTEDT Multifly®                | 62 - 64 |
| <b>RU</b> | Инструкция по применению – SARSTEDT Безопасные иглы-бабочки Safety-Multifly® / иглы-бабочки Multifly® | 65 - 67 |
| <b>SK</b> | Návod na Použitie – SARSTEDT Safety-Multifly® Ihly/Multifly® Ihly                                     | 68 - 70 |
| <b>SL</b> | Navodila za uporabo – SARSTEDT Kanila Safety-Multifly®/Kanila Multifly®                               | 71 - 73 |
| <b>SV</b> | Bruksanvisning – SARSTEDT Safety-Multifly®-Kanyler/Multifly®-Kanyler                                  | 74 - 76 |
| <b>TH</b> | คำแนะนำในการใช้งาน – SARSTEDT หัวเข็ม Safety-Multifly®/หัวเข็ม Multifly®                              | 77 - 79 |
| <b>TR</b> | Kullanım için talimatlar – SARSTEDT Safety-Multifly® İğneler/Multifly® İğneler                        | 80 - 82 |

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CLSI\* GP41 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.

\*CLSI (Clinical and Laboratory Standards Institute)

### Normen zu Nadeln/Kanülen:

DIN 13097-4: Hypodermic needles - Part 4: Point geometry, requirements and testing.

DIN EN ISO 6009: Hypodermic needles for single use - Colour coding for identification.

DIN EN 10088-1: Stainless steels - Part 1: List of stainless steels.

DIN EN ISO 9626: Stainless steel needle tubing for the manufacture of medical devices - Requirements and test methods.

### Normen zu Etiketten:

DIN EN 1041: Information supplied by the manufacturer of medical devices

DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

### Normen zu Sterilisation:

DIN EN ISO 11135: Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices

DIN EN 556-1: Sterilization of medical devices - Requirements for medical devices to be designated „STERILE“ - Part 1: Requirements for terminally sterilized medical devices

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DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

### Systembezogene Normen:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes

MDSAP: Medical Device Single Audit Program gemäß ISO 13485

DIN EN ISO 14971: Medical devices - Application of risk management to medical devices

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 SARSTEDT AG & Co. KG  
Sarstedtstr. 1  
D-51588 Nümbrecht  
www.sarstedt.com

## Intended use

The Safety-Multify® Needle / Multify® Needle is a sterile butterfly needle developed for single use. The Safety-Multify® Needle / Multify® Needle, together with a SARSTEDT S-Monovette®, is a system for venous blood collection. Additionally, the Safety-Multify® Needle / Multify® Needle can be used for a short infusion. The needle protector of the Safety-Multify® Needle is used immediately after blood collection of the last S-Monovette® / short infusion to cover the used needle and thus reduce potential accidental needle prick injuries. The products are to be used by trained personnel.

## Product description

The Safety-Multify® Needle / Multify® Needle is a sterile, pyrogen-free disposable product packed individually in a blister pack. Safety-Multify® Needles / Multify® Needles consist of a medical 3/4 inch (19 mm) long disposable stainless steel needle with a triple bevel, a flexible hose (80 mm or 200 mm) and an integrated needle holder (also called multi-adapter). There is a variant of the Safety-Multify® Needles / Multify® Needles without an integrated needle holder (multi-adapter) with a hose length of 240 mm. The products are available with the needle diameters 20, 21, 23, and 25 gauge (G). The needle protector of the Safety-Multify® Needles is used to cover the needle after use. The products contain only components that are DEHP and latex free.

## Safety instructions and warnings

Do not use the product if the blister packaging is damaged or the needle is bent or broken.

1. Follow the directions for handling.
2. Wear gloves and general personal protective equipment to protect yourself from blood and possible exposure to pathogens transmitted by blood.
3. Do not use force to pull the butterfly needle or prematurely activate the safety mechanism of the Safety-Multify® Needle (needle protector), as this would impair the functioning of the Safety-Multify® Needle.
4. Before using, check each individual packaging for damage and do not use if it is damaged.
5. Dispose of used, contaminated Safety-Multify® Needles / Multify® Needles immediately in suitable sharps disposal containers and dispose of equipment for blood collection (e.g. lancets, adapters) in sharps disposal containers as well.
6. After activating the safety mechanism of the Safety-Multify® Needle, do not open it again. Do not use force.
7. During use and disposal, always keep your hands behind the needle.
8. Do not bend the needle.
9. Injections and subcutaneous use are not allowed.
10. Treat all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures in your institution. In case of exposure to biological specimens or a needle prick injury, consult a doctor, as HBV, HCV, or HIV or other infectious diseases can be transmitted. Use needles with an integrated needle protector. The safety guidelines and procedures of your institution must always be followed.
11. Reusing the product can lead to infections, injuries or death.
12. When removing or disposing of the needle, check for possible leakage of blood or medication. It is therefore recommended to dispose of the Safety-Multify® Needle / Multify® Needle with the tip of the needle pointing up with the needle protector activated (Safety-Multify® Needle).
13. Do not cover the safety mechanism of the Safety-Multify® Needle with tape.
14. Before using for a short infusion, ensure that all the air is first removed (priming).

## Collecting specimens and handling

**READ THIS DOCUMENT COMPLETELY BEFORE STARTING A VEIN PUNCTURE.**

**Material required for blood collection that is not included:**

1. Use gloves, gown, goggles or other suitable protective gear to protect against pathogens or potentially infectious materials transmitted by blood.
2. Use disinfectant to clean the blood collection site. Follow your institution's guidelines for sterile sample collection to prepare the blood collection site and the directions for handling the S-Monovette®. Do not use any cleaning materials containing alcohol if the specimens are to be used for a blood alcohol test.
3. Clean, dry disposable swabs as used in your institution
4. A tourniquet.
5. A sharps disposal container.


## Preparation for blood collection

Ensure that the following materials are at hand before blood collection:

1. See above: Material required for blood collection that is not included.
2. All necessary S-Monovettes.
3. Labels for patient identification.

## Venous puncture with the Safety-Multify® Needle / Multify® Needle:

1. Apply tourniquet. Prepare the puncture site with a suitable disinfectant. Do not touch the puncture site after disinfection. Position the patient to allow easy access to the vein.
2. Open the blister packaging of the Safety-Multify® Needle / Multify® Needle aseptically using the peel-off technique with thumb and index finger. The transparent film must be facing upwards. After the disinfectant has dried, remove the protective sleeve from the needle and ensure that the bevelled side of the Safety-Multify® Needle/Multify® Needle is facing upwards. Perform venous puncture by pressing the coloured wings together upwards (1), then sliding the S-Monovette® into the multi-adapter of the Safety-Multify® Needle / Multify® Needle and securing it with a **¼ turn clockwise (connect)** (2). Draw blood with the S-Monovette® by carefully withdrawing the plunger up to the bottom of the S-Monovette® (3). Loosen the tourniquet while drawing blood into the first S-Monovette®.
3. After completing blood collection, first twist the S-Monovette® **¼ turn counterclockwise**(4) and then pull it out of the Safety-Multify® Needle / Multify® Needle (**disconnect**). The Safety-Multify® Needle/ Multify® Needle remains in the vein while blood is collected and when the blood collection tubes are changed. For multiple specimens – as described above – the subsequent S-Monovettes are connected to the Safety-Multify® Needle/ Multify® Needle and additional specimens are taken. At the end of blood collection, the last S-Monovette® is disconnected from the Safety-Multify® Needle / Multify® Needle and step A or step B are followed for the Safety-Multify® Needle and step C for the Multify® Needle.


Safety-Multify® Needle

Multify® Needle


**For the Safety-Multify® Needle:**

- A. The Safety-Multify® Needle is withdrawn from the vein (5) and the needle protector is activated outside the vein (6) by pushing the needle protector over the needle with the thumb and index finger, while holding the hose in place with the other fingers. The needle must perceptibly engage in the needle protector with an audible "click" (6).
- B. The needle protector is activated in parallel with withdrawing the Safety-Multify® Needle from the vein (5). The needle must perceptibly engage in the needle protector with an audible "click" (6).

Activating the safety mechanism is designed to be done with one hand.


**Activating the needle protector – needle outside the vein:** Withdraw the Safety-Multify® Needle from the vein and activate the needle protector.

**Activating the needle protector – needle in the vein:** Activate the needle protector in **parallel** with withdrawing the Safety-Multify® Needle from the vein.

**For the Multify® Needle:**

- C. The Multify® Needle is withdrawn from the vein (5).

4. Dispose of the Multify® Needle (6, pictogram of Multify® Needle) and the safely locked Safety-Multify® Needle immediately in a container for sharp objects (7, pictogram of Safety-Multify® Needle).

Handling for blood collection: see our handling video for the Safety-Multify® Needle: [www.sarstedt.com/en/download/videos/ Diagnostics / Venous blood](http://www.sarstedt.com/en/download/videos/Diagnostics/Venous%20blood)




**Short infusion with the Safety-Multify® Needle / Multify® Needle**

1. Prepare the puncture site with a suitable disinfectant. Do not touch the puncture site after disinfection. Position the patient to allow easy access to the vein.
2. Open the blister packaging of the Safety-Multify® Needle / Multify® Needle aseptically using the peel-off technique with thumb and index finger. The transparent film must be facing upwards. Remove the multi-adapter at the end of the hose by twisting and pulling it. After the disinfectant has dried, remove the Luer lock of the infusion system and connect the female Luer connection at the end of the Safety-Multify® / Multify® needle hose with male Luer connection of the infusion system. The air must be vented from the infusion system according to the guidelines of the institution (priming) (1,2).

**Maximum venting volume** of the Safety-Multify® Needles / Multify® Needles with needle diameters 20G, 21G, 23G, 25G:

|              |          |
|--------------|----------|
| 240 mm hose: | 0.450 ml |
| 200 mm hose: | 0.300 ml |
| 80 mm hose:  | 0.120 ml |

3. Remove the protective sleeve from the needle and ensure that the bevelled side of the Safety-Multify® Needle / Multify® Needle is facing upwards. Perform the venous puncture by pressing the coloured wings together upwards. Then a short infusion (max 24h) can be carried out.
4. To end the short infusion, the infusion system of the Safety-Multify® Needle / Multify® Needle is disconnected (3) and the Safety-Multify® Needle / Multify® Needle is then withdrawn from the vein by following step A or step B for the Safety-Multify® Needle and step C for the Multify® Needle.

5. Dispose of the Multify® Needle (5, pictogram of Multify® Needle) and the safely locked Safety-Multify® Needle (6, pictogram of Safety-Multify® Needle) immediately in a container for sharp objects.

Storage

Store the Safety-Multify® Needle / Multify® Needle at room temperature.  
 NOTE: Do not expose to direct sunlight. Exceeding the maximum recommended storage temperature may affect the quality of the Safety-Multify® Needle / Multify® Needle.

International / ISO standards

CLSI\* GP41 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.  
 \*CLSI (Clinical and Laboratory Standards Institute)

Standards for needles:

DIN 13097-4: Hypodermic needles - Part 4: Point geometry, requirements and testing.  
 DIN EN ISO 6009: Hypodermic needles for single use - Colour coding for identification.  
 DIN EN 10088-1: Stainless steels - Part 1: List of stainless steels.  
 DIN EN ISO 9626: Stainless steel needle tubing for the manufacture of medical devices - Requirements and test methods.

Standards for labels:

DIN EN 1041: Information supplied by the manufacturer of medical devices  
 DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

Standards for sterilisation:

DIN EN ISO 11135: Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices  
 DIN EN 556-1: Sterilization of medical devices - Requirements for medical devices to be designated „STERILE“ - Part 1: Requirements for terminally sterilized medical devices

Standards for packaging:

DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems  
 DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

System-related standards:

DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes  
 MDSAP: Medical Device Single Audit Program in accordance with ISO 13485  
 DIN EN ISO 14971: Medical devices - Application of risk management to medical devices  
 MDR 2017/745: EU Medical Device Directive

Key for icons and labels:

|                                                                                       |                                 |                                                                                     |                                  |
|---------------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------|----------------------------------|
|  REF | Article number                  |  | Use before                       |
|  LOT | Batch number                    |  | Follow the instructions for use  |
|      | CE marking                      |  | Store away from sunlight         |
|      | Manufacturer                    |  | Store in a dry place             |
| Store at room temperature                                                             |                                 |  | Do not use if package is damaged |
|      | Risk of contamination if reused |  | Pyrogen free                     |
|      | Do not resterilise              |  | DEHP free                        |

Technical modifications reserved

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

 SARSTEDT AG & Co. KG  
 Sarstedtstr. 1  
 D-51588 Nümbrecht  
 www.sarstedt.com

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Norme internazionali / ISO

CLSI\* GP41 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture.  
\*CLSI (Clinical and Laboratory Standards Institute)

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DIN 13097-4: Hypodermic needles - Part 4: Point geometry, requirements and testing.  
DIN EN ISO 6009: Hypodermic needles for single use - Colour coding for identification.  
DIN EN 10088-1: Stainless steels - Part 1: List of stainless steels.  
DIN EN ISO 9626: Stainless steel needle tubing for the manufacture of medical devices - Requirements and test methods.

Normen zu Etiketten:

DIN EN 1041: Information supplied by the manufacturer of medical devices  
DIN EN ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

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DIN EN ISO 11135: Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices  
DIN EN 556-1: Sterilization of medical devices - Requirements for medical devices to be designated „STERILE“ - Part 1: Requirements for terminally sterilized medical devices

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DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems  
DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

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DIN EN ISO 13485: Medical devices - Quality management systems – Requirements for regulatory purposes  
MDSAP: Medical Device Single Audit Program secondo la norma ISO 13485  
DIN EN ISO 14971: Medical devices - Application of risk management to medical devices  
MDR 2017/745: Regolamento Dispositivi Medici UE

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DIN 13097-4: Hypodermic needles – Part 4: Point geometry, requirements and testing.  
DIN EN ISO 6009: Hypodermic needles for single use – Colour coding for identification.  
DIN EN 10088-1: Stainless steels – Part 1: List of stainless steels.  
DIN EN ISO 9626: Stainless steel needle tubing for the manufacture of medical devices – Requirements and test methods.

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DIN EN 1041: Information supplied by the manufacturer of medical devices  
DIN EN ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements

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DIN EN ISO 11135: Sterilization of health-care products – Ethylene oxide – Requirements for the development, validation and routine control of a sterilization process for medical devices  
DIN EN 556-1: Sterilization of medical devices – Requirements for medical devices to be designated „STERILE“ – Part 1: Requirements for terminally sterilized medical devices

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DIN EN ISO 11607-1: Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems  
DIN EN ISO 11607-2: Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes

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DIN EN ISO 13485: Medical devices – Quality management systems – Requirements for regulatory purposes  
MDSAP: ISO 13485 에 따른 Medical Device Single Audit Program  
DIN EN ISO 14971: Medical devices – Application of risk management to medical devices  
MDR 2017/745: EU 의료기기 규정

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MDR 2017/745: AB tıbbi ürün yönetmeliği