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DE	Gebrauchshinweis – SARSTEDT CryoPure Gefäße	2 - 4
EN	Instructions for use – SARSTEDT CryoPure Collection Tubes	5 - 7
BG	Инструкции за употреба – SARSTEDT CryoPure епруветки	8 - 10
CS	Návod k použití – zkumavky SARSTEDT CryoPure	11 - 13
DA	Bruksanvisning – SARSTEDT CryoPure beholdere	14 - 16
EL	Οδηγίες χρήσεως – Φιαλίδια CryoPure SARSTEDT	17 - 19
ES	Instrucciones de uso – Recipientes CryoPure SARSTEDT	20 - 22
ET	Kasutusjuhend – SARSTEDT'i CryoPure'i katsutid	23 - 25
FR	Mode d'emploi – Cryotubes CryoPure SARSTEDT	26 - 28
HR	Uputa za upotrebu – SARSTEDT CryoPure epruvete	29 - 31
HU	Használati utasítás – SARSTEDT CryoPure fagyasztó csövek	32 - 34
IT	Istruzioni d'uso – Contenitori CryoPure SARSTEDT	35 - 37
KO	사용 설명서 – SARSTEDT CryoPure 튜브	38 - 40
NL	Gebruiksaanwijzing – SARSTEDT CryoPure-buisje	41 - 43
NO	Bruksanvisning – SARSTEDT CryoPure-beholdere	44 - 46
PL	Instrukcja obsługi – próbówki CryoPure SARSTEDT	47 - 49
PT	Instruções de utilização – Recipientes CryoPure SARSTEDT	50 - 52
RU	Инструкция по применению – пробирки SARSTEDT CryoPure	53 - 55
SK	Návod na Použitie – Skúmavky SARSTEDT CryoPure	56 - 58
SL	Navodila za uporabo – epruvete CroPure SARSTEDT	59 - 61
SV	Bruksanvisning – SARSTEDT CryoPure-rör	62 - 64
TH	คำแนะนำในการใช้งาน – หลอด CryoPure ของ SARSTEDT	65 - 67
TR	Kullanım için talimatlar – SARSTEDT CryoPure Tüpleri	68 - 70

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
Intended use

CryoPure collection tubes are single-use collection tubes for the storage and cryopreservation of solid and liquid biological samples (e.g. tissue) at temperatures as low as -196 °C. CryoPure collection tubes must be stored in the gas phase of liquid nitrogen. They may only be used by appropriately trained specialist personnel.

Product description

CryoPure collection tubes are available in four different designs which differ in capacity or screw cap colour (see Figure 1). All collection tubes consist of a clear plastic container with a round base and feature a screw cap, with three versions being available with an external thread and one with an internal thread and silicone seal. All tubes have a labelling area and a clear label specifying the fill quantity. The screw caps can have so-called colour-coded inserts of different colours pressed into them to allow samples to be coded.

Art. no.	Description	Screw cap	Nominal quantity
A 72.377.xxx	CryoPure collection tube 1.2 ml, sorted by colour	external thread	1.0 ml
B 72.379.xxx	CryoPure collection tube 2 ml, sorted by colour	external thread	1.8 ml
C 72.380.xxx	CryoPure collection tube 2 ml, sorted by colour	Internal thread	1.6 ml
D 72.383.xxx	CryoPure collection tube 5 ml, sorted by colour	external thread	4.5 ml

Accessories	
65.386.xxx	Colour-coded insert, sorted by colour
93.856.040	CryoRack 40

Handling

Freezing biological samples

1. Open the CryoPure collection tube by turning the screw cap.
2. Fill the CryoPure collection tube with your sample until you reach the specified fill mark.
Important! Over- or underfilling will increase the risk of contamination and explosions (see safety information / warnings).
3. Close the tube again with the screw cap.
4. Freeze your sample while observing the necessary cooling rate. When storing in a nitrogen tank, please ensure you observe the safety precautions.
5. If necessary, a colour-coded insert can be pressed into the screw cap prior to freezing to allow sample coding.

Thawing biological samples

1. Remove the CryoPure collection tube from the nitrogen tank in accordance with the safety precautions.
2. Allow your sample to thaw in a covered container, e.g., a water bath with a lid. Covering it will provide protection in the event of an explosion, which may occur in the case of improper use. Please observe the thawing speed necessary for your sample.
3. Open the CryoPure collection tube by turning the screw cap.
4. Remove the sample for further work steps.

Sample storage:

While the CryoPure collection tubes are designed to store samples at a temperature as low as -196 °C, they should be stored in the gas phase of liquid nitrogen, not in liquid nitrogen itself, for safety reasons.

IMPORTANT!

1. The storage duration and temperature of a filled CryoPure collection tube depend on the shelf life of the sample being examined and the analyte to be determined
2. The stability of the analyte should be assessed by the laboratory in question or on the basis of specialist literature or the instructions for use of the analytical equipment.
3. CryoPure collection tubes are not autoclavable.
4. CryoPure collection tubes are not intended to be centrifuged.

Transport:

CryoPure collection tubes are not suitable for transport. They are not primary containers for transport according to the ADR (Packaging Instruction P650) and the IATA Regulation.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled CryoPure collection tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and burned.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Safety information/Warnings

In the event of improper use, if the CryoPure collection tubes are stored in the liquid phase of nitrogen, liquid nitrogen may enter the tubes and contaminate the samples as falling temperatures form a vacuum in the tube. Nitrogen that has entered the sample may enter the gas phase during thawing and produce increased internal pressure, which could result in an explosion and the release of infectious material.

The following safety measures should always be observed when handling liquid nitrogen:

- Only store CryoPure collection tubes in the gas phase of liquid nitrogen.
- Fill CryoPure collection tubes up to the fill mark. Over- or underfilling will increase the risk of contamination and explosions.
- Tighten the CryoPure collection tubes by hand, avoiding overtightening them and refrain from using aids.
- Wear personal protective equipment. Liquid nitrogen can result in cold burns, frostbite and severe damage to eyes, and other injuries.
- Thaw CryoPure collection tubes in a covered container (e.g. water bath).

Additional safety measures:

- General precautions: Wear gloves and other general personal protective equipment to protect against possible exposure to potentially infectious specimen material and transmitted pathogens.
- Handle all biological specimens according to the guidelines and procedures in your facility. In case of direct contact with biological specimens, consult a doctor due to the risk of transmission of HIV, HCV, HBV or other infectious diseases.
- Do not use CryoPure collection tubes after their shelf life has expired.
- CryoPure collection tubes are intended for single use. After use they should be disposed of according to the applicable regulations. There is a risk of contamination if reused.

Purity

CryoPure collection tubes have received "Cryo Performance Tested" certification and fulfil the following purity criteria:

Successfully passing this defined testing procedure awards CryoPure tubes the following certifications:

- ✓ Sterile
- ✓ Pyrogen-free/endotoxin-free
- ✓ Non-cytotoxic
- ✓ Non-mutagenic
- ✓ DNA-free
- ✓ free from DNase/RNase
- ✓  

We guarantee compliance with the following limits:

- Sterility validated in accordance with the ISO 11137 series of standards
- Pyrogens/endotoxins <0.06 EU/ml
- Non-cytotoxic in accordance with the ISO 10993 series of standards
- Mutagen-free according to Ames Test II
- Human DNA <0.5 pg/µl, bacterial DNA <0.02 pg/µl
- DNase <1x10⁻⁵ U/µl, RNase <1x10⁻⁹ Kunitz units/µl

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE mark
	<i>In vitro</i> diagnostic device
	Pyrogen-free/non-pyrogenic
	Follow the instructions for use
	If reused: Risk of contamination
	Keep away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

The following is also applicable to sterile products:

	Sterile fluid path
	Single sterile barrier system with protective packaging outside
	Do not resterilise

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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