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DE	Gebrauchshinweis – SARSTEDT Archivierungsstopfen	2 - 3
EN	Instructions for use – SARSTEDT archiving stopper	4 - 5
BG	Инструкции за употреба – SARSTEDT Тапа за архивиране	6 - 7
CS	Návod k obsluze – archivační zátka SARSTEDT	8 - 9
DA	Brugsanvisning – SARSTEDT-lagringspropper	10 - 11
EL	Οδηγίες χρήσεως – Πώμα διατήρησης SARSTEDT	12 - 13
ES	Instrucciones de uso – Tapones para archivo SARSTEDT	14 - 15
ET	Kasutusjuhend – SARSTEDT arhiveerimiskork	16 - 17
FI	Käyttöohje – SARSTEDT näyteastian tulppa	18 - 19
FR	Mode d'emploi – Bouchon d'archivage SARSTEDT	20 - 21
HR	Uputa za upotrebu – SARSTEDT čepovi za arhiviranje	22 - 23
HU	Használati utasítás – SARSTEDT záródugó	24 - 25
IT	Istruzioni d'uso – Tappi per l'archiviazione SARSTEDT	26 - 27
KO	사용 설명서 – SARSTEDT 보관 마개	28 - 29
LT	Naudojimo instrukcijos – SARSTEDT archyvavimo kamščiai	30 - 31
LV	Lietošanas norādes – SARSTEDT arhivēšanas aizbāznis	32 - 33
NL	Gebruiksaanwijzing – SARSTEDT-bewaardoppen	34 - 35
NO	Bruksanvisning – SARSTEDT oppbevaringsplugger	36 - 37
PL	Instrukcja obsługi – zatyczki do archiwizacji SARSTEDT	38 - 39
PT	Instruções de utilização – rolha de arquivamento SARSTEDT	40 - 41
RO	Instrucțiuni de utilizare – dopuri de arhivare SARSTEDT	42 - 43
RU	Инструкция по применению – Архивные пробки для пробирок SARSTEDT	44 - 45
SK	Návod na Použitie – SARSTEDT archivačný zatlačací uzáver	46 - 47
SL	Navodila za uporabo – zamaški za arhiviranje SARSTEDT	48 - 49
SV	Bruksanvisning – SARSTEDT arkiveringspluggar	50 - 51
TH	คำแนะนำในการใช้งาน – จุกปิดหลอดเก็บตัวอย่างของ SARSTEDT	52 - 53
TR	Kullanım için talimatlar – SARSTEDT Arşivleme tapası	54 - 55
ZH	使用说明 – SARSTEDT存档塞	56 - 57

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Intended use

The archiving stoppers are push-in stoppers used to close specimen vessels and secondary vessels (Ø 13, 15, 16 mm) to reduce evaporation during storage. The vessels may contain human specimen material (e.g. serum, plasma). Storage may be at room temperature or refrigerator temperature depending on the stability of the parameters to be tested.

The archiving stoppers can be used for (fully) automatic or manual processing of specimens. The specifications of the device manufacturer must be noted.

The products are intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The archiving stoppers are plastic and are suitable for tubes with a diameter of 13, 15, or 16 mm; they are available with or without a filter. The filter seals on contact with liquid to prevent leakage of the liquid through the ventilation hole.

Safety information and warnings

1. General precautions: Wear gloves and general personal protective equipment to protect yourself from possible exposure to potentially infectious specimen material and transmitted pathogens.
2. Handle all biological specimens and collection utensils according to the guidelines and procedures in your institution. In case of direct contact with biological specimens, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. The archiving stoppers are intended for single use only. Used stoppers must not be put back onto the tubes. Dispose of all products in the hazardous biological waste containers.
4. The archiving stopper is for protection against evaporation and is not liquid tight.
5. Do not press the archiving stopper into overfilled tubes or specimen fluid.
6. Always store the closed tubes in an upright position, do not lie them down or mix them.
7. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year as indicated.

Storage


Store the products at room temperature.

Handling

READ THIS DOCUMENT COMPLETELY BEFORE STARTING HANDLING.

1. Use gloves, gown, goggles or other suitable protective gear to protect against pathogens or potentially infectious materials transmitted by blood.
2. Open the bag and remove the archiving stopper.
 - a. Place the archiving stopper onto the tube to be closed and push it in.
 - b. For (fully) automatic use, follow the device manufacturer's instructions.
3. Store vessels that have been closed with archiving stoppers in an upright position and in accordance with the specified storage conditions for the specimen material.
4. When necessary, pull out the archiving stopper to open the vessels for further steps. This can be done manually or (fully) automatically.

Stopper position according to tube diameter.


The push-in depth of the archiving stopper depends on the diameter of the specimen vessel to be closed (see illustration).

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and burned.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Key for symbols and labels:

	Article number
	Batch number
	CE mark
	For <i>in-vitro</i> diagnostic use
	Follow the instructions for use
	Risk of contamination if reused
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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