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Intended use


The prepared tube with boric acid is used as a sample tube and for collecting, transporting and storing urine samples from a primary tube for microbiological *in vitro* diagnostic determinations. The boric acid stabilises the microbial growth at room temperature for up to 48 hours.

The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The prepared tube with boric acid is available in two variants for the nominal volumes 25 ml or 10 ml and consists of a transparent plastic tube with a conical bottom and coloured screw cap. The version for 25 ml nominal volume also has a skirt.

Product overview

Type	Name	Length/Ø
A	Prepared tube with boric acid 10 ml	100x16 mm
B	Prepared tube with boric acid 25 ml	92x25 mm

The product is prepared with boric acid. For a prepared tube filled to the nominal volume, the average concentration of boric acid is 1.5%. This inhibits the growth of urinary pathogenic microorganisms in the urine at room temperature for up to 48 hours without impairing their viability. Studies were conducted for selected strains of bacteria:

- *Escherichia coli*
- *Klebsiella pneumoniae*
- *Proteus mirabilis*
- *Streptokokus faecalis*
- *Candida albicans*

Information on product-specific studies is available on the SARSTEDT homepage: www.sarstedt.com/en/products/diagnostic/urine/urine-monovetter/

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself from urine and possible exposure to pathogens transmitted by biological samples.
2. Handle all biological specimens according to the guidelines and procedures in your facility. In case of direct contact with biological specimens, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. The product is intended for single use. Dispose of the product and all collection utensils in containers for hazardous biological waste.
4. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year indicated.

Please note the safety data sheet at www.sarstedt.com/en/download/safety-data-sheets/.

Storage

Store the product at room temperature.

Transport

The product is a primary receptacle according to the ADR (Packing instruction P650) and IATA regulations.

Limitations

1. Significant underfilling can lead to false results.
NOTE: Filling up to the fill level mark ensures an optimal urine-preservative ratio. A deviation of +/-10% is allowed.
2. The use of the stabilised urine sample for strip tests and other non-microbiological analysis methods has not been tested and may lead to false results.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING A URINE TRANSFER.

Material required for the urine transfer

1. Gloves, gown, goggles or other suitable protective wear to protect against pathogens or potentially infectious materials transmitted by urine.
2. Required number of prepared tubes with boric acid.
3. Transfer pipettes or other suitable laboratory equipment for sample transfer.
4. Block rack or other suitable device for placing the prepared tubes with boric acid.
5. Disposal container for hazardous biological waste.

Handling for urine transfer into a prepared tube containing boric acid

1. Prepare the urine specimen by gently swirling the primary tube several times.
2. Hold the prepared tube of boric acid vertically with the closure facing upwards and make sure that the additive is on the bottom of the tube or the inner wall of the tube.
3. Open the prepared tube with boric acid and set the cap aside.
4. Place the opened prepared tube with boric acid in a suitable device in a stable position.
5. Carefully transfer the sample material up to the nominal volume into the prepared tube with boric acid. Pay attention to the "FILL TO LINE" marking.
6. Close the prepared tube with boric acid again with the screw cap previously set aside.
7. Invert the prepared tube with boric acid and swirl it about 5 times or until the boric acid is completely dissolved.
8. Place the prepared tube of boric acid in a suitable device.
9. Close the primary tube.
10. Dispose of all contaminated consumables in the hazardous biological waste containers.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled prepared sample tubes must be disposed of in suitable containers for hazardous biological waste.
4. Dispose of potentially contaminated material according to the regulations and guidelines of the facility.
5. Please take note of the information in the safety data sheet.

Product-specific standards and guidelines as amended / updated

DIN EN 14254: *In-vitro* diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans
 CLSI* GP16 "Urinalysis Approved Guideline"
 *CLSI (Clinical and Laboratory Standards Institute)

Further literature:

European Urinalysis Guidelines: Scand J. Clin. Lab. Invest 2000; 60: 1–96.

Information on product-specific studies is available on the SARSTEDT website:

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Manufacturer

Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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