

S-Monovette® Hirudin 1,6 ml

Jetzt mit reduziertem Volumen für die gewohnt gute Probenqualität

NEU

**Neu / sofort verfügbar:
1,6 ml**



Art.-Nr. 04.1959.001

**Bisher / nicht mehr verfügbar:
2,7 ml**



Art.-Nr. 04.1944.001

Nur so viel wie nötig – 1,6 ml Blut

- Ausreichendes Probenvolumen für erforderliches Analysenvolumen von 3 x 300µl
- Reduktion des Blutverlustes um 40 %
- Optimierte Geometrie gewährleistet immer optimale Probenqualität
- Vereinfachte Pipettierung durch hohen Füllstand

SARSTEDT AG & Co.
Postfach 12 20 · D-51582 Nümbrecht
Telefon +49 22 93 305-0
Telefax +49 22 93 305-3450
Service 0800 (Deutschland)
Telefon (0800) 0 83 305-0
info@sarstedt.com
www.sarstedt.com

 **SARSTEDT**

S-Monovette® Hirudin 1.6 ml

Now with reduced volume and consistently good sample quality

NEW

**new / available now:
1.6 ml**



Art. No. 04.1959.001

**old / no longer available:
2.7 ml**



Art. No. 04.1944.001

Only as much as necessary – 1.6 ml blood

- Sufficient sample volume for the required analysis volume of 3 x 300µl
- Blood loss reduced by 40%
- Optimised geometry ensures optimal sample quality at all times
- Easy pipetting thanks to high filling level

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com

 **SARSTEDT**