

Automatic-Sarpette®

New design for even more convenient pipetting



The SARSTEDT Automatic-Sarpette® in a new ergonomic design enables convenient and fatigue-free use even for prolonged pipetting.





Features and benefits:

- Ergonomic design and optimal weight distribution
- One-handed operation for convenient, efficient use
- Users can choose from three different pump speeds
- Precise and responsive control of filling and dispensing speed using two pipetting triggers
- All plastic and glass pipettes fit perfectly in silicone pipette holder with graduated steps
- A 50 ml pipette can be filled in under 10 seconds at maximum motor power
- Low-noise pump motor for pleasant working environment
- NiMH battery with a service life of up to 8 hours
- LED display for monitoring battery charge

Ordering information – Automatic-Sarpette®

Order no	Description	Packaging
90.189.200	Automatic-Sarpette® incl. EU charger, table stand, wall bracket, 2 replacement filters (0.45 µm and 0.20 µm)	1 pcs./case
90.189.202	Automatic-Sarpette® incl. GB charger, table stand, wall bracket, 2 replacement filters (0.45 µm and 0.20 µm)	1 pcs./case
90.189.203	Automatic-Sarpette® incl. US charger, table stand, wall bracket, 2 replacement filters (0.45 µm and 0.20 µm)	1 pcs./case
90.189.204	Automatic-Sarpette® incl. AUS charger, table stand, wall bracket, 2 replacement filters (0.45 µm and 0.20 µm)	1 pcs./case
90.189.205	Automatic-Sarpette® incl. KR charger, table stand, wall bracket, 2 replacement filters (0.45 µm and 0.20 µm)	1 pcs./case

Ordering information – replacement parts for Automatic-Sarpette®

	Order no	Description	Packaging
	90.189.220	Replacement filter for Automatic-Sarpette®, pore size: 0.45 µm	5 pcs./bag
	90.189.221	Replacement filter for Automatic-Sarpette®, pore size: 0.2 µm	5 pcs./bag
	90.189.222	Silicone pipette holder for Automatic Sarpette®	1 pcs./bag
	90.189.223	Replacement batteries for Automatic-Sarpette®	2 pcs./bag












Serological pipettes

- Manufactured from transparent polystyrene
- Larger pipetting volume due to negative graduations
- Variable working methods due to reverse scaling
- Optimized mouth for universal fit in the most common pipetting aids
- Guide ribs on the filter end piece of the 25 ml pipettes provide a secure fit in the retaining adapter of pipetting aids
- Simple volume identification with international color code on each individual blister pack and pipette
- Easy to open, anti-static paper/plastic packaging
- Available in individual sterile* packaging or in a bag of 25 units



* Individually packaged pipettes are sterile and certified pyrogen-free/endotoxin-free and non-cytotoxic.

Ordering information – Serological pipettes 1 ml, 2 ml, 5 ml, 10 ml, 25 ml, 50 ml

Order no	Total volume/ graduations		Design	Color code	Packaging unit/case
86.1251.001	1 ml	1/100 ml	Plugged, ind. wrapped sterile		1/1,000
86.1251.025	1 ml	1/100 ml	Plugged, sterile, in 25 units		25/1,000
86.1252.001	2 ml	1/100 ml	Plugged, ind. wrapped sterile		1/1,000
86.1252.025	2 ml	1/100 ml	Plugged, sterile, in 25 units		25/1,000
86.1253.001	5 ml	1/10 ml	Plugged, ind. wrapped sterile		1/500
86.1253.025	5 ml	1/10 ml	Plugged, sterile, in 25 units		25/500
86.1254.001	10 ml	1/10 ml	Plugged, ind. wrapped sterile		1/500
86.1254.025	10 ml	1/10 ml	Plugged, sterile, in 25 units		25/500
86.1685.001	25 ml	2/10 ml	Plugged, ind. wrapped sterile		1/200
86.1685.020	25 ml	2/10 ml	Plugged, sterile, in 20 units		20/200
86.1256.001	50 ml	1/2 ml	Plugged, ind. wrapped sterile		1/90

SARSTEDT, Inc.
1025 St. James Church Road
P.O. Box 468
Newton, NC 28658 - 0468
Toll free: +1 800 257 5101
Fax: +1 828 465 4003
customerservice.us@sarstedt.com
www.sarstedt.com