

Cell culture products
by SARSTEDT



**READY.
SET.
GROW!**

SARSTEDT: Excellence in cell culture

Try it out for
free & with no
obligation!



cellculture.
sarstedt.com/en

Set new standards in cell cultivation with SARSTEDT!

Our innovative product line, consisting of flasks, dishes, plates and a wide range of complements, provides many different growth surfaces for all *in vitro* conditions and is tailored to the requirements of most cell types. Our precise production and highest quality standards pave the way for excellent science.

The SARSTEDT colour codes make it easy to find the right receptacle at a glance. **Put your trust in the scientific expertise and product portfolio of SARSTEDT.**

Cell culture dishes

- ✓ SARSTEDT SUREGrip – Outer ring enables a **secure grip through haptic feedback**
- ✓ Clear arrows on the base and lid for **uniform positioning**
- ✓ **Gas exchange** through air vents in the lid
- ✓ **Safe stacking** thanks to distinct rings
- ✓ Dishes in \varnothing 35 mm and \varnothing 60 mm with a **grid for cloning experiments**

Cell culture flasks

- ✓ **Accessible** for pipettes and scrapers
- ✓ Large **labelling fields**
- ✓ **Engraved and printed** graduated scales
- ✓ **Good stability** with stacking rim
- ✓ **Reduction of contamination risk** due to angled neck and anti-drip rim

Cell culture plates

- ✓ Alphanumerical labelling on the edge and between the wells for **better orientation** when filling
- ✓ **Reduced risk of contamination** when pipetting thanks to free-standing wells
- ✓ **Secure hold** due to non-slip side grids
- ✓ Transparent side walls allow **visual inspection** of the medium
- ✓ Integrated air vents and condensation rings in the lid for **gas exchange and minimisation of evaporation**

The **SARSTEDT colour code** indicates the three surface properties of the culture vessels that create optimal conditions for cell growth.

RED – The standard surface for adherent cells

YELLOW – The SARSTEDT Cell⁺ surface for difficult adherent cells and for serum-reduced/serum-free conditions of cultivated cells

GREEN – The SARSTEDT suspension culture surface



Find the right surface for your cell culture here:



 **SARSTEDT**

SARSTEDT AG & Co. KG

P.O. Box 12 20

D-51582 Nümbrecht

Phone: +49 2293 305 0

Fax: +49 2293 305 3992

export@sarstedt.com

www.sarstedt.com

We support the entire cell culture workflow with an extensive product portfolio.

For an overview, see cellculture.sarstedt.com/en

Let's grow! SARSTEDT



SARSTEDT