

Число лунок	96 проб. установл.	96	96	96	96	384	96	384
Половинная/полная кромка	нет	нет	половинная	полная	половинная	полная	половинная	полная
Профиль	высок.	высок.	высок.	низкий	низкий	-	"Lightcycler"	"Lightcycler"
Кат. № ПЦР-планшета	72.985	72.1978 72.1978.010	72.1979 72.1979.010 72.1979.003 72.1979.201 72.1979.700 72.1979.102 72.1979.132	72.1980 72.1980.010 72.1980.201 72.1980.600 72.1980.700	72.1981 72.1981.010	72.1984.202	72.1982.202	72.1985.202

**BioRad® / MJ Research®**

CFX96 Touch™ Real-Time PCR				●				
CFX384 Touch™ Real-Time PCR						●		
CFX Automation System II								
T100™ Thermal Cycler	●	●		●	●			
T100™ Thermal Cycler		●	●	●		●		
C1000Touch™ Thermal Cycler		●	●	●	●	●		
iCycler iQ™ Thermal Cycler	●	●	●					
iQ4™ Thermal Cycler	●	●	●					
iQ5™ Thermal Cycler	●	●	●					
MyCycler™ Thermal Cycler	●	●	●					
Chromo4™		●		●				
Opticon™, Opticon2™				●				
BaseStation™				●				

**Corbett Research® / Qiagen®**

Palm Cycler 96 лунки		●						
Palm Cycler 384 лунки				●				

**Eppendorf®**

Mastercycler® nexus	●	●						
Mastercycler® ep realplex		●	●			●		
Mastercycler® gradient	●	●		●				
Mastercycler® ep gradient	●	●	●	●				
Mastercycler® pro	●	●		●				

**Ericom®**

Deltacycler			●	●				
SingleBlock			●	●				
TwinBlock			●	●				

**MWG®**

Primus 96 лунки		●		●				
Primus 384 лунки						●		
The Q-Lifecycler		●	●	●				

**Roche®**

Система Lightcycler® 96							●	
Система Lightcycler® 480							●	●

**Stratagene® / Agilent®**

Система ПЦР в реальном времени AriaMx				●	●			
Mx3000P™	●	●						
Mx3005P™	●	●	●	●				
Mx4000™	●	●	●					
Gradient Cycler		●		●				
Robocycler 384 лунки						●		

**Teche®**

Cyclogene		●		●				
Flexigene		●	●	●				
Genius / Genius Quad		●	●	●				
OMN-E		●	●	●				
PCR Express	●	●	●			●		
Primus 96		●						
Px2 / PxE		●	●			●		
Quantica			●	●				
TC412 / TC512		●		●		●		
Touchgene / Touchgene Gradient		●	●	●	●			