



### Product description

Order number	93.854.125
Product description	S-Monovette® rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, orange, material: PP

### Product characteristics

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

### Size

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**S-Monovette® rack D12, Ø opening: 12 mm, 5 x 10, orange**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product orange

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 20  
Piece(s) / pallet 1560  
Depth of case 428 mm  
Width of case 323 mm  
Height of case 130 mm  
Case volume 0.018 cbm  
Weight of product 0.1214 kg  
Weight of case 2.8 kg  
EAN of inner box 4038917122689  
EAN case 4038917018524

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2023-11-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**