



### **Product description**

Order number	78.684
Product description	Absorbent liner, suitable for mailing container 84 x 30 mm, (LxW): 75 x 65 mm, 100 piece(s)/case

### **Product characteristics**

Application	Mailing of samples, Multi-Safe disposal boxes
Suitable for	mailing container 84 x 30 mm

### **Size**

Width of product	65 mm
Length of product	75 mm

### **Material & colours**

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Absorbent liner, suitable for mailing container 84**  
**x 30 mm, (LxW): 75 x 65 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      2000  
Type of smallest  
subpackaging                              case  
Piece(s) / inner box                      100  
Piece (s) / outer case                      2000  
Piece(s) / pallet                            120000  
Depth of case                              410 mm  
Width of case                                260 mm  
Height of case                               305 mm  
Case volume                                0.0325 cbm  
Weight of product                         0.0007 kg  
Weight of case                              1.9 kg  
EAN of inner box                          4038917083362  
EAN case                                      4038917002295

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

