


Boxen für 15 & 50ml Röhren




- Passend für 15 & 50 ml Röhren
- Numerische Codierung zur einfachen Probenidentifizierung
- Einsetzbar im Temperaturbereich von -80°C bis 121°C
- Abmessung L x B x H (ohne Deckel): 144 x 144 x 127 (102) mm

Bestellinformation

Bestell Nummer	Bezeichnung	Farbe	Kapazität	Verpackung
93.880	Lagerbox für 15 ml Röhren	Grün	36 (6 x 6)	2/Beutel 12/Karton
93.881	Lagerbox für 50 ml Röhren	Blau	16 (4 x 4)	

SARSTEDT AG & Co.
Postfach 12 20 · D-51582 Nümbrecht
Telefon (+49) 0 22 93 30 50
Telefax (+49) 0 22 93 305-282
☎ Service 0800 (Deutschland)
Telefon (0800) 0 83 30 50
info@sarstedt.com
www.sarstedt.com

Boxes for 15 & 50 ml Tubes



Technical modifications reserved

This document may contain information on products that may not be available in particular countries

- Designed to accommodate 15 & 50 ml tubes
- Each space numbered for easy sample identification
- Ideal for use in temperatures ranging from -80°C to 121°C
- L x W x H (without lid): 144 x 144 x 127 (102) mm

Ordering Information

Art. No.	Description	Colour	Capacity	Packaging
93.880	Storage box for 15 ml tubes	green	36 (6 x 6)	2/bag
93.881	Storage box for 50 ml tubes	blue	16 (4 x 4)	12/case

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-122
info@sarstedt.com
www.sarstedt.com