

To make your
workflow excellent.




Conversion to Multi-Safe REC

Safe and sustainable disposal!

Brand NEW

Also available as a
sustainable version
made from recycled
plastic!



 Pre-analytics

Multi-Safe made from virgin PP

Box type	Vol.	Tests			Order no.	Pieces/ box	
		TRBA 250	UN 3291	ISO 23907			
Multi-Safe vario 2000	2.0 l	✓	✓	✓	77.3893.020	50	
Multi-Safe twin plus 500	0.5 l	✓	✓	✓	77.3884.005	126	
Multi-Safe twin plus 500 BHZ	0.5 l	✓	✓	✓	77.3884.105	126	
Multi-Safe twin plus 1500	1.5 l	✓	✓	✓	77.3884.015	94	
Multi-Safe twin plus 1500 BHZ	1.5 l	✓	✓	✓	77.3884.115	94	
Multi-Safe twin plus 2000	2.0 l	✓	✓	✓	77.3884.020	80	
Multi-Safe twin plus 2000 BHZ	2.0 l	✓	✓	✓	77.3884.120	80	
Multi-Safe twin plus 2500	2.5 l	✓	✓	✓	77.3884.025	70	
Multi-Safe twin plus 2500 BHZ	2.5 l	✓	✓	✓	77.3884.125	70	
Multi-Safe twin plus 3000	3.0 l	✓	✓	✓	77.3884.030	62	
Multi-Safe twin plus 3000 BHZ	3.0 l	✓	✓	✓	77.3884.130	62	
Multi-Safe twin plus 5000	5.0 l	✓	✓	✓	77.3884.050	34	
Multi-Safe twin plus 5000 BHZ	5.0 l	✓	✓	✓	77.3884.150	34	

Multi-Safe REC made from recycled PP

Up to **64%** recycled content
in the entire box!
Up to **51%**
CO₂e savings compared
to new material.



Box type	Vol.	Tests			Order no.	Pieces/ box	
		TRBA 250	UN 3291	ISO 23907			
Multi-Safe vario 2000 REC	2.0 l	✓	✓	✓	77.3893.220	50	
Multi-Safe twin plus 500 REC	0.5 l	✓	✓	✓	77.3884.205	126	
Multi-Safe twin plus 500 BHZ REC	0.5 l	✓	✓	✓	77.3884.305	126	
Multi-Safe twin plus 1500 REC	1.5 l	✓	✓	✓	77.3884.215	94	
Multi-Safe twin plus 1500 BHZ REC	1.5 l	✓	✓	✓	77.3884.315	94	
Multi-Safe twin plus 2000 REC	2.0 l	✓	✓	✓	77.3884.220	80	
Multi-Safe twin plus 2000 BHZ REC	2.0 l	✓	✓	✓	77.3884.320	80	
Multi-Safe twin plus 2500 REC	2.5 l	✓	✓	✓	77.3884.225	70	
Multi-Safe twin plus 2500 BHZ REC	2.5 l	✓	✓	✓	77.3884.325	70	
Multi-Safe twin plus 3000 REC	3.0 l	✓	✓	✓	77.3884.230	62	
Multi-Safe twin plus 3000 BHZ REC	3.0 l	✓	✓	✓	77.3884.330	62	
Multi-Safe twin plus 5000 REC	5.0 l	✓	✓	✓	77.3884.250	34	
Multi-Safe twin plus 5000 BHZ REC	5.0 l	✓	✓	✓	77.3884.350	34	

Multi-Safe made from virgin PP

Box type	Vol.	Tests			Order no.	Pieces/ box	
		TRBA 250	UN 3291	ISO 23907			
Multi-Safe sani 200	0.2 l	✓	✓	✓	77.3890.002	110	
Multi-Safe mini 300	0.3 l	✓	✓	✓	77.3890.003	135	
Multi-Safe euroMatic® 3000	3.0 l	✓	✓	✓	77.3897.030	55	
Multi-Safe euroMatic® 3000 BHZ	3.0 l	✓	✓	✓	77.3897.230	55	
Multi-Safe euroMatic® 4000	4.0 l	✓	✓	✓	77.3897.040	50	
Multi-Safe euroMatic® 4000 BHZ	4.0 l	✓	✓	✓	77.3897.240	50	
Multi-Safe euroMatic® 6000	6.0 l	✓	✓	✓	77.3897.060	40	
Multi-Safe euroMatic® 6000 BHZ	6.0 l	✓	✓	✓	77.3897.260	40	
Multi-Safe euroMatic® 7000	7.0 l	✓	✓	✓	77.3897.070	34	
Multi-Safe euroMatic® 7000 BHZ	7.0 l	✓	✓	✓	77.3897.270	34	

Multi-Safe REC made from recycled PP

Up to **64%** recycled content
in the entire box!
Up to **51%**
CO₂e savings compared
to new material.



Box type	Vol.	Tests			Order no.	Pieces/ box	
		TRBA 250	UN 3291	ISO 23907			
Multi-Safe sani 200 REC	0.2 l	✓	✓	✓	77.3890.202	110	
Multi-Safe mini 300 REC	0.3 l	✓	✓	✓	77.3890.203	135	
Multi-Safe euroMatic® 3000 REC	3.0 l	✓	✓	✓	77.3897.530	55	
Multi-Safe euroMatic® 3000 BHZ REC	3.0 l	✓	✓	✓	77.3897.630	55	
Multi-Safe euroMatic® 4000 REC	4.0 l	✓	✓	✓	77.3897.540	50	
Multi-Safe euroMatic® 4000 BHZ REC	4.0 l	✓	✓	✓	77.3897.640	50	
Multi-Safe euroMatic® 6000 REC	6.0 l	✓	✓	✓	77.3897.560	40	
Multi-Safe euroMatic® 6000 BHZ REC	6.0 l	✓	✓	✓	77.3897.660	40	
Multi-Safe euroMatic® 7000 REC	7.0 l	✓	✓	✓	77.3897.570	34	
Multi-Safe euroMatic® 7000 BHZ REC	7.0 l	✓	✓	✓	77.3897.670	34	

From blood collection to “Multi-Safe” disposal – Simple and safe



Safely through the day

As a renowned manufacturer of laboratory products and medical devices, we consider it our duty to provide safe disposal options for our single-use products.

As large and small waste is disposed of differently, our wide range of products offers a multitude of disposal options for our users. A variety of choice that is reflected in the different volumes, openings and many variants, even through to the final closure.

We ensure that all boxes meet the TRBA 250 (German regulations) as a minimum requirement so that we can guarantee the safe handling of the disposal box. Naturally our Multi-Safe disposal boxes also comply with international regulations, such as UN 3291 and ISO 23907.

As an independent family business, we are committed to the environment and future generations. Guided by this responsibility, we continuously optimise our actions in order to promote social growth, minimise our ecological footprint and consistently protect human health.

Our commitment to sustainable practices for our product lines is shown through actions, not merely words. When developing, manufacturing and optimising our products, we attach great importance to resource-efficient processes and products, using recyclable materials and promoting reusable systems.

The new product line of Multi-Safe REC disposal boxes made from recycled plastic (Post-Consumer Recycled = PCR) is one current example of our sustainable approach, along with other sustainable product solutions. With a recycled content of up to 64% on the entire box and a CO₂e reduction of up to 51%, it fulfils the requirement of both ISO 23907 and UN 3291.

Brand NEW

Also available as a sustainable version made from recycled plastic!




If you have any questions,
we'll be happy to help!

Visit our website:
www.sarstedt.com

Pre-analytical workflow made by SARSTEDT

Utilise the synergy of our
coordinated systems.



SARSTEDT AG & Co. KG

Sarstedtstraße 1
D-51588 Nümbrecht

Phone: +49 2293 305 0

export@sarstedt.com
www.sarstedt.com

Discover the 360°
solutions for your
pre-analytics from
SARSTEDT:



workflow.sarstedt.com