

# Pressure Infusion Cuffs

Improved "Comfort" with each application!





## Pressure Infusion Cuff "Comfort"

### Improved "Comfort" with each application!

- For invasive pressure measurement, pressure infusions and pressure transfusions
- Latex-free & suitable for MRI
- Washable cover
- Large tube diameter and large valve outlet for easy and fast pumping up; quick deflating
- Design of the valve-gauge-unit allows a comfortable one-hand operation of bellows and valve
- Optimised fit also for longish infusion bags
- Modular design allows the exchange of single components

Latex-free  
&  
suitable for MRI



### Insert pocket – for recording patient, ward, medication etc.:

- Improves identification and documentation
- Prevents writing on the cover



## Pressure Infusion Cuff "Classic"

### Quality that pays

- For invasive pressure measurement, pressure infusions and pressure transfusions
- Modular design allows the exchange of single components
- Latex-free versions available
- Washable cover



### Pressure Infusion Cuff 500 For single use

- For invasive pressure measurement, pressure infusions and pressure transfusions
- For transfusion/infusion bags up to 500 ml volume
- With pressure gauge
- Cost-efficient and hygienic
- Suitable for single patient use



# Pressure Infusion Cuffs

## Comfort-Line – latex free

Order No.	Description	Pck.
97.8100.600	<b>Pressure Infusion Cuff “S-Comfort”</b> - Up to 500 ml volume - With fixed insertion	1
97.8100.610	<b>Pressure Infusion Cuff “K-Comfort”</b> - Up to 1000 ml volume - Variable fit due to velcro	1
97.8100.620	<b>Pressure Infusion Cuff “FR-Comfort”</b> - Up to 1000 ml volume - Solid back and velcro for variable fit - With additional insert pocket for noting down patient, ward, medication etc.	1

## Spare Parts

Order No.	Description	Pck.
97.8100.630	<b>Rubber ballon latex-free</b> • Inner tube diameter 8 mm	5
97.8100.631	<b>Valve-Gauge-Unit</b> • Valve-Gauge-Unit for inner tube diameter of 8 mm	1
97.8100.632	<b>Bellows for valve-gauge-unit</b> • Bellows for the connection to the valve-gauge-unit, red, latex-free	1
97.8100.633	<b>Cover for pressure infusion cuff “S-Comfort”</b>	1
97.8100.634	<b>Cover for pressure infusion cuff “K-Comfort”</b>	1
97.8100.635	<b>Cover for pressure infusion cuff “FR-Comfort”</b>	1

## Classic-Line

Order No.	Description	Pck.
97.8100.000	<b>Pressure infusion cuff K</b> - Up to 1000 ml volume - Variable fit due to velcro	1
97.8100.011	<b>Pressure infusion cuff K latex-free</b> - Up to 1000 ml volume - Variable fit due to velcro	1
97.8100.100	<b>Pressure infusion cuff Standard</b> - Up to 500 ml volume - With fixed insertion	1
97.8100.111	<b>Pressure infusion cuff Standard latex-free</b> - Up to 500 ml volume - With fixed insertion	1
97.8100.200	<b>Pressure infusion cuff FR</b> - Up to 1000 ml volume - Solid back and velcro for variable fit	1
97.8100.211	<b>Pressure infusion cuff FR latex-free</b> - Up to 1000 ml volume - Solid back and velcro for variable fit	1

## Spare Parts

Order No.	Description	Pck.
97.8100.010	<b>Rubber ballon, latex free</b>	5
97.8100.001	<b>Rubber ballon, containing latex</b>	5
97.8100.210	<b>Cover for pressure infusion cuff FR</b>	1
97.8100.002	<b>Cover for pressure infusion cuff K</b>	1
97.8100.102	<b>Cover for pressure infusion cuff Standard</b>	1
97.8100.003	<b>Pressure gauge</b>	5
97.8100.004	<b>Bellows with valve, latex-free</b>	5

## Disposable Pressure Infusion Cuff

Order No.	Description	Pck.
97.8100.410	<b>Pressure Infusion Cuff 500</b> - Up to 500 ml volume - With fixed insertion	10

Distributed by:

SARSTEDT AG & Co. KG  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone +49 2293 305-0  
Fax +49 2293 305-3992  
export@sarstedt.com  
www.sarstedt.com

