

To make your
workflow excellent.



Multi-Safe disposal boxes

For safe disposal!



Pre-analytics

From blood collection to "Multi-Safe" disposal – Simple and safe



Safely through the day

As a renowned manufacturer of laboratory products and medical devices, we consider it our duty to provide safe disposal options for our single-use products.

As large and small waste is disposed of differently, our wide range of products offers a multitude of disposal options for our users. A variety of choice that is reflected in the different volumes, openings and many variants, even through to the final closure.

We ensure that all boxes meet the TRBA 250 (German regulations) as a minimum requirement so that we can guarantee the safe handling of the disposal box. Naturally our Multi-Safe disposal boxes also comply with international regulations, such as UN 3291 and ISO 23907.

As an independent family business, we are committed to the environment and future generations. Guided by this responsibility, we continuously optimise our actions in order to promote social growth, minimise our ecological footprint and consistently protect human health.

Our commitment to sustainable practices for our product lines is shown through actions, not merely words. When developing, manufacturing and optimising our products, we attach great importance to resource-efficient processes and products, using recyclable materials and promoting reusable systems.

The new product line of Multi-Safe REC disposal boxes made from recycled plastic (Post-Consumer Recycled = PCR) is one current example of our sustainable approach, along with other sustainable product solutions. With a recycled content of up to 64% on the entire box and a CO₂e reduction of up to 51%, it fulfils the requirement of both ISO 23907 and UN 3291.

Brand NEW

Also available as a sustainable version made from recycled plastic!





Table of Contents

Multi-Safe sani & mini / Multi-Safe sani REC , Multi-Safe mini REC	Page 6 - 7
Multi-Safe vario / Multi-Safe vario REC	Page 8 - 9
Multi-Safe twin plus / Multi-Safe twin plus REC	Page 10 - 13
Multi-Safe euroMatic® / Multi-Safe euroMatic® REC	Page 14 - 15
Multi-Safe medi	Page 16 - 17
Multi-Safe eco	Page 18 - 19
Multi-Safe steri	Page 20 - 21
Multi-Safe "Custom products"	Page 22
Safety tray	Page 23
Multi-Safe accessories	Page 24 - 25

Multi-Safe sani & mini

Multi-Safe for on the go – the "small ones" with a "big" job

- Safe needle box for on the go
- Suitable for outpatient use by doctors, nurses and rescue services or for doctors and medical cases
- One-handed opening and closing of the day and final closure
- Lid already assembled
- Compact

Multi-Safe sani 200

Brand NEW

Also available as a sustainable version made from recycled plastic!



Up to **64%** recycled content!
Up to **51%** CO₂e reduction!



Multi-Safe mini 300

Up to **64%**
recycled content!
Up to **51%**
CO₂e reduction!

160 mm
77 mm
41 mm



Luer systems



Multi-Safe sani 200

Volume	Material	Container/lid colour	Design	Piece(s)/carton	Piece(s)/pallet	Case dimensions	Order no
0.2 ml	Polypropylene		Multi-Safe sani 200	110	3,300	39.7 x 39.7 x 32.5 cm	77.3890.002
0.2 ml	Recycled polypropylene		Multi-Safe sani 200 REC	110	3,300	39.7 x 39.7 x 32.5 cm	77.3890.202

Multi-Safe mini 300

Volume	Material	Container/lid colour	Design	Piece(s)/carton	Piece(s)/pallet	Case dimensions	Order no
0.3 ml	Polypropylene		Multi-Safe mini 300	135	2,700	60 x 40 x 31 cm	77.3890.003
0.3 ml	Recycled polypropylene		Multi-Safe mini 300 REC	135	2,700	60 x 40 x 31 cm	77.3890.203

Multi-Safe vario 2000

Large opening 2.0

- Easy assembly of the lid with snap hook
- The size of the disposal opening can be adjusted by breaking out sections as required. Please break out before use.
- Day closure and final closure enables one-handed operation
- QR code on the box enables quick and direct access to the instructions for use
- Various removal and twist-off aids for needles, Luer systems and insulin pens
- Integrated stamp in the lid enables the easy and safe disconnection of the insulin pen when closing
- Compatible with commercially available trays
- Certified in line with DIN ISO 23907, UN 3291 and TRBA 250


Brand NEW

Also available as a sustainable version made from recycled plastic!



Compatible
with commercially available trays





Multi-Safe vario 2000

Volume	Material	Container/lid colour	Design	Piece(s)/carton	Piece(s)/pallet	Case dimensions	Order no
2.0 l	Polypropylene		Multi-Safe vario 2000	50	1,000	60 x 40 x 35.9 cm	77.3893.020
2.0 l	Recycled polypropylene		Multi-Safe vario 2000 REC	50	1,000	60 x 40 x 35.9 cm	77.3893.220



Multi-Safe twin plus

The disposal box with added safety!

- Optimised assembly of the lid with snap hook
- Innovative lid: the day closure and final closure enable one-handed operation
- Disposal opening can be gradually enlarged by breaking out the sections. Please break out before use.
- IML label (in-mould labelling) stamped directly during production with product labelling in 9 languages
- The QR code on the box enables quick and direct access to the instructions for use
- Various removal and turning aids
- In particular, the integrated stamp for insulin pens introduces significant benefits as, by pressing down the lid, the stamp pushes the insulin pens completely into the box
- Various table and wall holders provide the option of fastening the boxes in different areas and prevent toppling (see accessories p. 24)
- Also available as version with biohazard labelling
- Certified in line with DIN ISO 23907, UN 3291 and TRBA 250

Brand NEW

Also available as a sustainable version made from recycled plastic!



Up to **64%**
recycled content!
Up to **51%**
CO₂e reduction!





With innovative day closure and final closure with one-handed operation and "Press for final closure" imprint




Integrated QR code for direct online access to the operating instructions



Also available with biohazard labelling




Directly integrated label in 9 languages – label also available as customised solution



twin plus accessories (see page 24)





**Multi-Safe
twin plus 500 REC**



**Multi-Safe
twin plus 1500 REC**



**Multi-Safe
twin plus 2000 REC**



**Multi-Safe
twin plus 5000 REC**

**Multi-Safe
twin plus 2500 REC**




**Multi-Safe
twin plus 3000 REC**





Multi-Safe twin plus

Volume	Material	Container/ lid colour	Design	Piece(s)/ carton	Piece(s)/ pallet	Case dimensions	Order no
0.5 l	Polypropylene		Multi-Safe twin plus 500	126	3,024	59.9 x 26.6 x 45.1 cm	77.3884.005
0.5 l	Polypropylene		Multi-Safe twin plus 500 with biohazard labelling	126	3,024	59.9 x 26.6 x 45.1 cm	77.3884.105
0.5 l	Recycled polypropylene		Multi-Safe twin plus 500 REC	126	3,024	59.9 x 26.6 x 45.1 cm	77.3884.205
0.5 l	Recycled polypropylene		Multi-Safe twin plus 500 REC with biohazard labelling	126	3,024	59.9 x 26.6 x 45.1 cm	77.3884.305
1.5 l	Polypropylene		Multi-Safe twin plus 1500	94	1,128	80 x 40 x 38 cm	77.3884.015
1.5 l	Polypropylene		Multi-Safe twin plus 1500 with biohazard labelling	94	1,128	80 x 40 x 38 cm	77.3884.115
1.5 l	Recycled polypropylene		Multi-Safe twin plus 1500 REC	94	1,128	80 x 40 x 38 cm	77.3884.215
1.5 l	Recycled polypropylene		Multi-Safe twin plus 1500 REC with biohazard labelling	94	1,128	80 x 40 x 38 cm	77.3884.315
2.0 l	Polypropylene		Multi-Safe twin plus 2000	80	960	80 x 40 x 38 cm	77.3884.020
2.0 l	Polypropylene		Multi-Safe twin plus 2000 with biohazard labelling	80	960	80 x 40 x 38 cm	77.3884.120
2.0 l	Recycled polypropylene		Multi-Safe twin plus 2000 REC	80	960	80 x 40 x 38 cm	77.3884.220
2.0 l	Recycled polypropylene		Multi-Safe twin plus 2000 REC with biohazard labelling	80	960	80 x 40 x 38 cm	77.3884.320
2.5 l	Polypropylene		Multi-Safe twin plus 2500	70	840	80 x 40 x 38 cm	77.3884.025
2.5 l	Polypropylene		Multi-Safe twin plus 2500 with biohazard labelling	70	840	80 x 40 x 38 cm	77.3884.125
2.5 l	Recycled polypropylene		Multi-Safe twin plus 2500 REC	70	840	80 x 40 x 38 cm	77.3884.225
2.5 l	Recycled polypropylene		Multi-Safe twin plus 2500 REC with biohazard labelling	70	840	80 x 40 x 38 cm	77.3884.325
3.0 l	Polypropylene		Multi-Safe twin plus 3000	62	744	80 x 40 x 38 cm	77.3884.030
3.0 l	Polypropylene		Multi-Safe twin plus 3000 with biohazard labelling	62	744	80 x 40 x 38 cm	77.3884.130
3.0 l	Recycled polypropylene		Multi-Safe twin plus 3000 REC	62	744	80 x 40 x 38 cm	77.3884.230
3.0 l	Recycled polypropylene		Multi-Safe twin plus 3000 REC with biohazard labelling	62	744	80 x 40 x 38 cm	77.3884.330
5.0 l	Polypropylene		Multi-Safe twin plus 5000	34	408	80 x 40 x 38 cm	77.3884.050
5.0 l	Polypropylene		Multi-Safe twin plus 5000 with biohazard labelling	34	408	80 x 40 x 38 cm	77.3884.150
5.0 l	Recycled polypropylene		Multi-Safe twin plus 5000 REC	34	408	80 x 40 x 38 cm	77.3884.250
5.0 l	Recycled polypropylene		Multi-Safe twin plus 5000 REC with biohazard labelling	34	408	80 x 40 x 38 cm	77.3884.350

Multi-Safe euroMatic®

Innovative design and engineering – the premium model

Up to **64%**
recycled content!
Up to **51%**
CO₂e reduction!

Automatic lid opening by **pressing the red button**



- Square container shape
- Tough container and lid
- Anti-tip mounting with adhesive adapter and/or wall holder (see accessories on page 24)
- Innovative day and final closures allow one-handed operation of the lid (pressing the red “button” automatically opens the lid)
- Hinges on the day lid enable multiple opening/closing operations, while the secure final closure is made visible by folding out the yellow ring
- One-handed operation and large disposal opening enable fast and safe disposal of bulky objects
- Side grooves at the top edge visually confirm the closure of the container
- Side handle
- Various pull-off and twist-off aids

Brand NEW

Also available as a sustainable version made from recycled plastic!




Integrated QR code for direct online access to the operating instructions




Also available with **biohazard labelling**



euroMatic accessories® (see page 24)



euroMatic® 7000 REC**euroMatic® 6000 REC****euroMatic® 4000 REC****euroMatic® 3000 REC****Multi-Safe euroMatic®**

Volume	Material	Container/lid colour	Design	Piece(s)/carton	Piece(s)/pallet	Case dimensions	Order no
7 l	Polypropylene		Multi-Safe euroMatic® 7000	34	306	80.2 x 39.2 x 51.4 cm	77.3897.070
7 l	Polypropylene		Multi-Safe euroMatic® 7000 with biohazard labelling	34	306	80.2 x 39.2 x 51.4 cm	77.3897.270
7 l	Recycled polypropylene		Multi-Safe euroMatic® 7000	34	306	80.2 x 39.2 x 51.4 cm	77.3897.570
7 l	Recycled polypropylene		Multi-Safe euroMatic® 7000 with biohazard labelling	34	306	80.2 x 39.2 x 51.4 cm	77.3897.670
6 l	Polypropylene		Multi-Safe euroMatic® 6000	40	360	80.2 x 39.2 x 51.4 cm	77.3897.060
6 l	Polypropylene		Multi-Safe euroMatic® 6000 with biohazard labelling	40	360	80.2 x 39.2 x 51.4 cm	77.3897.260
6 l	Recycled polypropylene		Multi-Safe euroMatic® 6000	40	360	80.2 x 39.2 x 51.4 cm	77.3897.560
6 l	Recycled polypropylene		Multi-Safe euroMatic® 6000 with biohazard labelling	40	360	80.2 x 39.2 x 51.4 cm	77.3897.660
4 l	Polypropylene		Multi-Safe euroMatic® 4000	50	450	80.2 x 39.2 x 51.4 cm	77.3897.040
4 l	Polypropylene		Multi-Safe euroMatic® 4000 with biohazard labelling	50	450	80.2 x 39.2 x 51.4 cm	77.3897.240
4 l	Recycled polypropylene		Multi-Safe euroMatic® 4000	50	450	80.2 x 39.2 x 51.4 cm	77.3897.540
4 l	Recycled polypropylene		Multi-Safe euroMatic® 4000 with biohazard labelling	50	450	80.2 x 39.2 x 51.4 cm	77.3897.640
3 l	Polypropylene		Multi-Safe euroMatic® 3000	55	495	80.2 x 39.2 x 51.4 cm	77.3897.030
3 l	Polypropylene		Multi-Safe euroMatic® 3000 with biohazard labelling	55	495	80.2 x 39.2 x 51.4 cm	77.3897.230
3 l	Recycled polypropylene		Multi-Safe euroMatic® 3000	55	495	80.2 x 39.2 x 51.4 cm	77.3897.530
3 l	Recycled polypropylene		Multi-Safe euroMatic® 3000 with biohazard labelling	55	495	80.2 x 39.2 x 51.4 cm	77.3897.630

Multi-Safe medi

Large volume box with handle and label


- Round container shape
- The large, round footprint eliminates the need for additional holders
- Various designs (with or without handle, additional label)
- Closures have been developed so that they automatically fall backwards on opening and uncover the disposal opening
- Side grooves at the top edge visually confirm the closure of the container
- Various pull-off and twist-off aids

Multi-Safe **medi 6**

In the future it will also be available as a sustainable version made from recycled plastic!



Multi-Safe **medi 11** with handle and label



Multi-Safe medi

Volume	Material	Container/ lid colour	Design	Piece(s)/ carton	Piece(s)/ pallet	Case dimensions	Order no
11 l	Polypropylene	■ ■	Multi-Safe medi 11	-	250	-	77.3898.011
11 l	Polypropylene	■ ■	Multi-Safe medi 11	10	160	39.5 x 43.0 x 58.5 cm	77.3898.511
11 l	Polypropylene	■ ■	Multi-Safe medi 11 with handle	-	250	-	77.3898.111
6 l	Polypropylene	■ ■	Multi-Safe medi 6	-	480	-	77.3898.006
6 l	Polypropylene	■ ■	Multi-Safe medi 6	21	336	39.5 x 43.0 x 58.5 cm	77.3898.506
6 l	Polypropylene	■ ■	Multi-Safe medi 6 with handle	-	480	-	77.3898.106
6 l	Polypropylene	■ ■	Multi-Safe medi 6 with handle and label	21	336	39.5 x 43.0 x 58.5 cm	77.3898.206

*UN testing only for medi 6, medi 6 with handle, medi 6 with handle and label

Multi-Safe eco

A safe and stable container for the disposal of clinical waste

- Rectangular container shape
- No additional holders needed thanks to the large footprint
- Various designs
- Multi-Safe 25 & 50 - lid-in-lid system provides variable disposal opening
- Liquid-tight final closure thanks to the adhesive strip fastened on the side of the lid
- Side grooves at the edge visually confirm the closure of the container
- Can be stacked empty and full
- Suitable for oncological waste, cytostatic drugs
- Collection and transport container for small disposal boxes

Multi-Safe 25





Multi-Safe eco

Volume	Material	Container/lid colour	Design	Piece(s)/pallet	Order no
50 l	Polypropylene	■ ■	Multi-Safe eco 50	125	77.3899.050
25 l	Polypropylene	■ ■	Multi-Safe eco 25	160	77.3899.025
50 l	Polypropylene	■ ■	Multi-Safe 50	125	77.3899.150
25 l	Polypropylene	■ ■	Multi-Safe 25	160	77.3899.125

Multi-Safe steri

Huge, practical, (autoclave-)safe

- Rectangular container
- No additional holders needed thanks to the large footprint
- Safe and stable for the disposal of infectious waste
- Autoclave-safe at 121°C to 134 °C
- Small lid melts at around 90°C and releases the autoclave process through a grid
- The small lids can be purchased to label boxes that have already been autoclaved (77.3888.800 - yellow; other colours are possible)
- Easy to transport with handles in the bottom and top part
- Side grooves on the lid visually confirm the closure
- Can be stacked empty and full
- Collection and transport container for small disposal boxes

Multi-Safe **steri 35**



Multi-Safe **steri 60**



Multi-Safe **steri**

Volume	Material	Container/lid colour	Design	Piece(s)/pallet	Order no
60 l	Polypropylene	■ ■	Multi-Safe steri 60	120	77.3899.060
35 l	Polypropylene	■ ■	Multi-Safe steri 35	150	77.3899.035

Custom products

Customised Multi-Safe disposal box

- Above a certain quantity, Multi-Safe disposal boxes can be ordered as custom products (e.g. different colour, personalised labels, etc.)
- Additional certification with your company address exclusively for your customised Multi-Safe disposal box can also be requested.
- Multi-Safe and Multi-Safe eco are not available in special colours.
- If you have any questions about your customised Multi-Safe disposal box, your personal SARSTEDT contact or our helpline will be happy to help you.



Safety tray



Version I

1 x S-Monovette® block rack D17 with handle
2x dispenser box with lid
1x Multi-Safe twin plus 500
3x multi-purpose cup
1x disinfection spray



Version II

2x dispenser box with lid
1x Multi-Safe vario 2000
6x multi-purpose cup
1x disinfection spray



Version III

2x dispenser box with lid
1x Multi-Safe twin plus 500
8x multi-purpose cup
1x disinfection spray



Version IV

1 x S-Monovette® block rack D17 with handle
2x dispenser box with lid
1x Multi-Safe vario 2000
1x multi-purpose cup
1x disinfection spray



Version V

1x S-Monovette® rack D17 with handle
1x dispenser box with lid
1x Multi-Safe Box twin plus 2000
3x multi-purpose cup
1x disinfection spray

Safety tray

Design	Dimensions	Piece(s)/carton	Order no
Safety tray without insertion plate	429 x 283 x 80 mm	1/box	93.1650
Version I	-	1/box	93.1650.010
Version II	-	1/box	93.1650.020
Version III	-	1/box	93.1650.030
Version IV	-	1/box	93.1650.040
Version V	-	1/box	93.1650.050

Accessories

Material	Design	Piece(s)/carton	Order no	Order no. recycled
Multi-Safe twin plus 500	0.5 l suitable for versions I and III	126	77.3884.005 77.3884.105	77.3884.205 77.3884.305
Multi-Safe vario 2000	2.0 l suitable for versions II and IV	50	77.3893.020 77.3893.120	77.3893.220 77.3893.420
Multi-Safe twin plus 2000	2.0 l suitable for version V	80	77.3884.020 77.3884.120	77.3884.220 77.3884.320
Dispenser box	With lid for Safety Tray	35	95.813.001	
Multi-purpose cup	125 ml for Safety Tray	1,000	75.1355	
Block D17	With handle, grey	12	93.851	
Block D17	With handle, red	12	93.851.171	
Block D17	With handle, yellow	12	93.851.172	
Block D17	With handle, green	12	93.851.173	
Block D17	With handle, blue	12	93.851.174	
Block D17	With handle, orange	12	93.851.175	

Accessories

For simple and secure mounting on walls, tables or smooth surfaces

Wall holder



Adhesive adapter



Wall holder twin plus*



Wall holder vario 2000*



Wall holder euroMatic®



Fastening plate



Round tube holder

Wall holder

Order no	Article description	Material	Piece(s)/carton	Suitable for series
95.965.001	Adhesive adapter	Polypropylene	150	twin plus, euroMatic®
95.963	Wall holder with adhesive adapter	Stainless steel	30	twin plus
95.963.001	Wall holder euroMatic®	stainless steel	20	euroMatic®
95.963.006	Round tube holder	POM	10	twin plus, euroMatic®, vario
95.963.005	Fastening plate	Stainless steel	20	twin plus, euroMatic®, vario
95.963.007	Wall holder twin plus	Stainless steel	20	twin plus
95.963.020	Wall holder vario 2000	Stainless steel	5	vario

*Round tube holder and fastening plate also suitable for the twin plus and vario wall holders

Multi-Safe label



This label can be ordered in addition to the disposal boxes.

Absorbent materials



ISO 23907 requires the following absorbent materials to be placed in the disposal box to prevent the leakage of liquid residue (does not apply for twin plus series). It is important to ensure that adequate absorbent material is used to cover the entire fill volume of a box.

Label

Article description	Piece(s)/inner box	Units/box	Order no
Label 70 x 80 mm	2,000	8,000	84.1390.022

Absorbent materials (ISO 23907)

Absorption	Minimum order quantity (units)	for Multi-Safe	Order no
16 ml per unit	1,000	euroMatic®	78.669
9 ml per unit	2,000	sani, mini	78.684
1,200 ml per unit	100	euroMatic® (6.0 and 7.0 litres), medi	95.1712



Accessories information
www.sarstedt.com/en/products/diagnostic/multi-safe/accessories/

OVERVIEW

Volume	Box type	Tests			Height	Ø of lid	Ø of base	
		TRBA 250	UN 3291	ISO 23907				
0.2 l	Multi-Safe sani 200	✓	✓	✓	131 mm	75 x 47 mm	57 x 32 mm	
0.3 l	Multi-Safe mini 300	✓	✓	✓	160 mm	98 mm	41 x 77 mm	
2.0 l	Multi-Safe vario 2000	✓	✓	✓	140 mm	243 x 102 mm	95 x 240 mm	
0.5 l	Multi-Safe twin plus 500	✓	✓	✓	162 mm	89 mm	57 mm	
1.5 l	Multi-Safe twin plus 1500	✓	✓	✓	125 mm	179 x 176 mm	112 mm	
2.0 l	Multi-Safe twin plus 2000	✓	✓	✓	162 mm	179 x 176 mm	112 mm	
2.5 l	Multi-Safe twin plus 2500	✓	✓	✓	197 mm	179 x 176 mm	112 mm	
3.0 l	Multi-Safe twin plus 3000	✓	✓	✓	233 mm	179 x 176 mm	112 mm	
5.0 l	Multi-Safe twin plus 5000	✓	✓	✓	379 mm	179 x 176 mm	112 mm	
3.0 l	Multi-Safe euroMatic® 3000	✓	✓	✓	185 mm	166 x 43 mm	132 x 132 mm	
4.0 l	Multi-Safe euroMatic® 4000	✓	✓	✓	245 mm	166 x 43 mm	132 x 132 mm	
6.0 l	Multi-Safe euroMatic® 6000	✓	✓	✓	330 mm	166 x 43 mm	132 x 132 mm	
7.0 l	Multi-Safe euroMatic® 7000	✓	✓	✓	390 mm	166 x 43 mm	132 x 132 mm	
6.0 l	Multi-Safe medi 6	✓	✓	✓	195 mm	230 mm	190 mm	
11.0 l	Multi-Safe medi 11	✓	-	✓	220 mm	290 mm	257 mm	
6.0 l	Multi-Safe medi 6	✓	✓	✓	195 mm	230 mm	190 mm	
6.0 l	Multi-Safe medi 6 with handle	✓	✓	✓	195 mm	230 mm	190 mm	
11.0 l	Multi-Safe medi 11 with handle	✓	-	✓	220 mm	290 mm	257 mm	
6.0 l	Multi-Safe medi 6, handle, label	✓	✓	✓	195 mm	230 mm	190 mm	
11.0 l	Multi-Safe medi 11	✓	-	✓	220 mm	290 mm	257 mm	
25.0 l	Multi-Safe eco 25	-	✓	-	260 mm	295 x 395 mm	230 x 330 mm	
50.0 l	Multi-Safe eco 50	-	✓	-	500 mm	295 x 395 mm	230 x 330 mm	
35.0 l	Multi-Safe steri 35	-	✓	-	330 mm	393 x 393 mm	340 x 340 mm	
60.0 l	Multi-Safe steri 60	-	✓	-	520 mm	393 x 393 mm	340 x 340 mm	
25.0 l	Multi-Safe 25	-	✓	-	260 mm	295 x 395 mm	230 x 330 mm	
50.0 l	Multi-Safe 50	-	✓	-	500 mm	295 x 395 mm	230 x 330 mm	

	Piece(s)/carton	Piece(s)/pallet	Order no	Order no. with biohazard	Order no. recycled	Order no. recycled, with biohazard
110	3,300	77.3890.002	Not available, as no label	77.3890.202	Not available, as no label	
135	2,700	77.3890.003	Not available, as no label	77.3890.203	Not available, as no label	
50	1,000	77.3893.020		77.3893.220		
126	3,024	77.3884.005	77.3884.105	77.3884.205	77.3884.305	
94	1,128	77.3884.015	77.3884.115	77.3884.215	77.3884.315	
80	960	77.3884.020	77.3884.120	77.3884.220	77.3884.320	
70	840	77.3884.025	77.3884.125	77.3884.225	77.3884.325	
62	744	77.3884.030	77.3884.130	77.3884.230	77.3884.330	
34	408	77.3884.050	77.3884.150	77.3884.250	77.3884.350	
55	495	77.3897.030	77.3897.230	77.3897.530	77.3897.630	
50	450	77.3897.040	77.3897.240	77.3897.540	77.3897.640	
40	360	77.3897.060	77.3897.260	77.3897.560	77.3897.660	
34	306	77.3897.070	77.3897.270	77.3897.570	77.3897.670	
-	480	77.3898.006				
-	250	77.3898.011				
21	336	77.3898.506				
-	480	77.3898.106				
-	250	77.3898.111				
21	336	77.3898.206				
10	160	77.3898.511				
-	160	77.3899.025				
-	125	77.3899.050				
-	150	77.3899.035				
-	120	77.3899.060				
-	160	77.3899.125				
-	125	77.3899.150				

Handling

Safe handling of euroMatic®, steri, medi and twin plus disposal boxes



Centre the lid on the lower section of the box. The red day lid and the handle face forwards.



Simultaneously press on the opposite corners of the lid to join to the lower section of the box.



The edge of the lid must be firmly inserted around the entire circumference, until it audibly and visibly locks into place.

Safe handling of the medi disposal boxes



Centre the lid on the lower section of the box. The yellow day lid faces forwards.



Simultaneously press on the opposite corners of the lid to join to the lower section of the box.



The edge of the lid must be firmly inserted around the entire circumference, until it audibly and visibly locks into place.

Safe handling of eco/Multi-Safe disposal boxes (lid-in-lid)



Centre the lid on the lower section of the box.



Simultaneously press on the opposite corners of the lid to join to the lower section of the box.



The edge of the lid must be firmly inserted around the entire circumference, until it audibly and visibly locks into place.

Day closure

1 x lock into place
(for reuse)



Final closure

2 x lock into place
(completely sealed)



Multi-Safe Instructions for Use

Intended use: The disposal container is intended exclusively for the disposal of pointed and sharp objects, such as needles and scalpel blades from the medical field.

Safety instructions

- Only use containers of a size that is suitable for holding the items to be disposed of.
- Store the box out of the reach of children.
- The lid must be on the container and locked into place before it is filled.
- Connect containers with the recommended adhesive adapter by rotating or hanging in the wall holder to prevent accidents.
- Do not use the day lid to press down the items to be disposed of.
- Scalpels must be disposed of in the container with particular care. If too much force is used when throwing scalpels in or if other items are placed on top, there is a risk of the angle changing and damaging the container walls or the base of the container.
- Only place items to be disposed of upright in the container.
- Do not press items into the container with force.
- Do not place any liquids in the container.
- Do not fill the box above the fill line (risk of injury!).
- Do not place hands or any other objects in the container (risk of injury).
- Do not throw the box. Do not shake the box. Do not drop the box.
- Before sealing the container, ensure that no items are projecting through the opening.
- Close the box by hand, without applying too much force.
- Do not push the box against your body (e.g. chest) for the final closure.
- Do not compress a full box.
- Use the handle.
- Before disposing of the box, check carefully that the lid is tightly closed.

ISO 23907

- One of the absorbent liners from page 25 must be used to ensure the compliant use of the box.

Disposal and the environment

- Plastics used in our products are environmentally neutral in terms of disposal.
- The plastic can be sent to landfill, as they do not release any harmful substances into the air, earth or water.
- When disposed of in incinerators, the plastics PS, PP, PC and PE serve as valuable sources of energy to maintain the temperature in the combustion chamber – the more plastic in the waste, the less heating oil is used. Burning these plastics does not release any hazardous emissions.

Recommendation:

Fill the Multi-Safe container only to about 2/3 of its volume.
Do not overfill the Multi-Safe: Risk of injury!

Fill line
observed



Overfilled



Material info

	Polystyrene	Polypropylene	Polycarbonate	Polyethylene
Abbreviation	PS	PP	PC	HD-PE High Density LD-PE Low Density
Optical properties	transparent, shiny surface with 90% light transmission (at 400-800 nm)	translucent, shiny surface	transparent, 88% light transmission (at 400-800 nm)	translucent to opaque, waxy surface
General mechanical properties	Low elongation at break and thermal resistance, very good electrical insulating properties, not suitable for high centrifugal accelerations	high breaking strength, resistant to stress cracking, stable form and high rigidity	very good mechanical, optical, electrical and thermal properties, can be autoclaved and sterilised	Relatively low tensile strength and surface hardness, high toughness, soft to firm, sensitive to stress cracking, water-repellent
Autoclavability	Not suitable	Products made of PP can be autoclaved at up to 121 °C without any appreciable loss of mechanical properties. The user must check whether other product characteristics are affected in terms of the desired use.	Products made of PC can be autoclaved up to 121 °C without any appreciable loss of mechanical properties. The user must check whether other product characteristics are affected in terms of the desired use.	Not suitable
Max. Continuous operating temperature* ¹	60 - 70 °C	100 - 110 °C	115 - 125 °C	HD-PE 70 - 80 °C LD-PE 60 - 75 °C
Temporary maximum operating temperature* ¹	75 - 80 °C	120 - 140 °C	125 - 140 °C	HD-PE 90 - 120 °C LD-PE 80 - 90 °C
Use at sub-zero temperatures* ²	Not well-suited	Limited suitability* ¹	down to -80 °C	Limited suitability* ¹
Density g/cm ³	1.05	0.90	1.19	HD-PE 0.95 LD-PE 0.92
Flammability	Flammable	Flammable	Flammable	Flammable
	300 - 400 °C	300 - 360 °C	380 - 450 °C	350 - 360 °C
Moisture absorption	< 0.1%	< 0.1%	0.1 - 0.3%	< 0.1%
General resistance to chemicals	PS is resistant to saline solutions, bases, non-oxidising acids, alkalis and alcohols. Benzene, volatile oils, strongly oxidising agents and flavourings attack PS, forming stress cracks.	PP is resistant to aqueous solutions of inorganic salts, acids, alkalis and to organic solvents up to 60 °C. Alcohols, esters and ketones also do not attack PP. Aromatic and halogenated hydrocarbons, oxidising substances such as concentrated nitric acid and fats, oils and waxes at high temperatures cause PP to swell.	PC is resistant to higher concentrations of mineral acids, many organic acids (e.g. carbonic acid, oleic acid and citric acid), oxidants, reducing agents, neutral and acidic saline solutions, a range of fats and oils, saturated, aliphatic and cycloaliphatic hydrocarbons and alcohols (except methanol). Bases, ammonia gas and its soluble form, and amines attack PC. PC is soluble in a range of industrial solvents. Other organic compounds like benzene, acetone and carbon tetrachloride eventually lead to expansion or swelling.	PE has a high resistance to chemicals. The chemical resistance of HD-PE is generally higher than that of LD-PE. Diluted acids, alkalis, alcohol, oil and water and saline solutions do not affect PE. Concentrated, oxidising acids, such as nitric acid and halogens have a corrosive effect.
Disposal	PS is a pure hydrocarbon and thus environmentally neutral in terms of disposal. The controlled burning of this plastic does not release any hazardous emissions.	PP is a pure hydrocarbon and thus environmentally neutral in terms of disposal. The controlled burning of this plastic does not release any hazardous emissions.	PC is a pure hydrocarbon and thus environmentally neutral in terms of disposal. The controlled burning of this plastic does not release any hazardous emissions.	PE is a pure hydrocarbon and thus environmentally neutral in terms of disposal. The controlled burning of this plastic does not release any hazardous emissions.

*¹ Suitability depends on the type of plastic used and the type of stress


*² Warning: The plastics become brittle at sub-zero temperatures. The use of products at sub-zero temperatures should be tested in advance for the relevant application. This information is intended as a guide and does not guarantee product properties.

If you have any questions,
we'll be happy to help!

Visit our website:
www.sarstedt.com

Pre-analytical workflow made by SARSTEDT

Utilise the synergy of our
coordinated systems.



SARSTEDT AG & Co. KG

Sarstedtstraße 1
D-51588 Nümbrecht

Phone: +49 2293 305 0

export@sarstedt.com
www.sarstedt.com

Discover the 360°
solutions for your
pre-analytics from
SARSTEDT:



workflow.sarstedt.com