

Ordering information

Blood collection systems



S-Monovette®

Serum Gel CAT

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1602.001	01.1602	01.1602.088
	4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1935.001	04.1935	04.1935.088
	4.7 ml	75/91 x 15 mm	■	50 / 50 / 500	03.1524.001	03.1524	
	4.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1925.001	04.1925	04.1925.088
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1923.001		04.1923.088
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1905.001		
	1.1 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1667.001		

Serum CAT

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1601.100	01.1601.200	01.1601.288
	5.5 ml	75/91 x 15 mm	■	50 / 50 / 500	03.1397.100	03.1397.200	
	4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1934.100	04.1934.200	04.1934.288
	4.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1924.100	04.1924.200	
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1943.100	04.1943.200	
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1904.100		
	1.6 ml	75/91 x 13 mm	■	50 / 50 / 500		04.2911.200	
	1.2 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1663.100		

Lithium Heparin Gel LH / Lithium Heparin Gel+ LH

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1634.100		01.1634.288
	4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1940.100	04.1940.200	04.1940.288
	4.7 ml	75/91 x 15 mm	■	50 / 50 / 500	03.1631.100	03.1631.200	
	4.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1927.100		
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1928.100	04.1928.200	04.1928.288
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1907.100		
	1.1 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1669.100		
	Lithium Heparin Gel+ 4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1954.100	04.1954.200	04.1954.288
	4.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1953.100	04.1953.200	04.1953.288
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1952.100	04.1952.200	04.1952.288

Lithium Heparin LH / Lithium Heparin LH liquid

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1604.100	01.1604.200	01.1604.288
	5.5 ml	75/91 x 15 mm	■	50 / 50 / 500	03.1628.100	03.1628.200	
	4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1936.100	04.1936.200	04.1936.288
	4.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1920.100	04.1920.200	
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1929.100		04.1929.288
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1906.100	04.1906.200	
	1.6 ml	75/91 x 13 mm	■	50 / 50 / 500		04.2913.200	
	1.2 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1666.100		
liquid	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1608.100		

Sodium Heparin NH

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1613.100		
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1913.100		

Packaging: All S-Monovettes listed above are wrapped in a sterile packaging.
 SP: Smallest sub-packaging of a product
 IB: Inner box, the SP is packaged in the IB
 OC: Outer case, the outer case is usually also the minimum order quantity
 MOQ: Minimum order quantity

Ammonium Heparin AH

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500		01.1603.200	

Citrate 9NC, 3.2%

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	5.4 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1930.100		
	4.3 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1922.100	04.1922.200	
	3.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1919.100	04.1919.200	04.1919.288
	3.0 ml	66/82 x 11 mm	■	50 / 50 / 500	05.1165.100	05.1165.200	
	2.9 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1902.100	04.1902.200	04.1902.288
	1.8 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1955.100	04.1955.200	04.1955.288
	1.4 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1668.100		

EDTA K3E

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1605.100		
	4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1931.100	04.1931.200	
	4.0 ml	75/91 x 15 mm	■	50 / 50 / 500	03.1068.100		
	4.0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1921.100		04.1921.288
	3.4 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1914.100	04.1914.200	
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1917.100	04.1917.200	
	2.7 ml	66/82 x 11 mm	■	50 / 50 / 500	05.1167.100	05.1167.200	
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1901.100	04.1901.200	04.1901.288
	1.8 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1951.100		04.1951.288
	1.6 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1956.100		
	1.6 ml	66/82 x 11 mm	■	50 / 50 / 500	05.1081.100	05.1081.200	
	1.2 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1664.100		

EDTA K2E

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1915.100	04.1915.200	
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1900.100	04.1900.200	
	1.2 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1660.100		

EDTA Gel K2E

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1621.200		
	4.9 ml	90/106 x 13 mm	■	50 / 50 / 500	04.1932.101		

Fluoride EDTA FE

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	5.5 ml	75/91 x 15 mm	■	50 / 50 / 500		03.1626.200	
	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1918.100	04.1918.200	04.1918.288
	2.7 ml	66/82 x 11 mm	■	50 / 50 / 500	05.1073.100	05.1073.200	
	2.6 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1903.100	04.1903.200	
	1.2 ml	66/82 x 8 mm	■	50 / 50 / 500	06.1665.100		

Fluoride Heparin FH

Preparation	Nominal volume	Length without/with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label	Order no. colored not writable label	Order no. pre-barcoded label
	2.7 ml	66/82 x 11 mm	■	50 / 50 / 500		05.1076.200	

S-Monovette®

Neutral Z

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored not writable label
	7.5 ml	92/108 x 15 mm	□	50 / 50 / 500	01.1728.001
	4.9 ml	90/106 x 13 mm	□	50 / 50 / 500	04.1926.001
	2.7 ml	66/82 x 11 mm	□	50 / 50 / 500	05.1729.001

Special S-Monovette® & S-Sedivette®

Special S-Monovette®

Description	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no.
S-Monovette® Serum Gel LightPROTECT	4,0 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1925.003
S-Monovette® GlucoEXACT	3.1 ml	66/82 x 11 mm	■	50 / 50 / 500	05.1074.201
S-Monovette® GlucoEXACT	3.1 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1945.201
S-Monovette® GlucoEXACT	3.1 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1945.288
S-Monovette® Metal Analysis*	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1604.400
S-Monovette® Metal Analysis*	7.5 ml	92/108 x 15 mm	■	50 / 50 / 500	01.1609.201
S-Monovette® ThromboExact	2.7 ml	66/82 x 11 mm	■	50 / 50 / 500	05.1168.001
S-Monovette® Hirudin	1.6 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1959.001
S-Monovette® PFA, Citrate 9NC, 3.8%	3.8 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1910.001
S-Monovette® CTAD	2.9 ml	65/81 x 13 mm	■	50 / 50 / 500	04.1909.200
S-Monovette® CTAD	4.3 ml	75/91 x 13 mm	■	50 / 50 / 500	04.2909.200
S-Monovette® ACD-A	5.7 ml	90/106 x 13 mm	■	50 / 50 / 500	01.1624.001
S-Monovette® Homocystein HCY-Z-Gel	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1908.001
S-Monovette® DNA Exact	2.7 ml	75/91 x 13 mm	■	50 / 50 / 500	04.1948.001
S-Monovette® RNA Exact	≤ 2.4 ml	100/116 x 15 mm	■	50 / 50 / 500	01.2048.001
S-Monovette® cfDNA Exact	9.2 ml	100/116 x 15 mm	■	50 / 50 / 500	01.2040.001

* To use with S-Monovette® Safety needle for Metal Analysis 85.1162.600

S-Sedivette® Citrate 4NC

Preparation	Nominal volume	Length without/ with cap x Ø	Cap color	Packaging (SP/IB/OC)	Order no. colored writable label
	3.5 ml	130/146 x 8 mm	■	50 / 50 / 500	06.1690.100

Sedimentation rack for S-Sedivette®

Illustration	Description	Packaging (SP/IB/OC)	Order no.
	Sedimentation rack with graduated back wall	1 / 1 / 1	90.1090

Needles and adapters


Safety-Multifly® needle for S-Monovette®

Illustration	Ø/Length	Color code	Description	Tube length	Packaging (SP/IB/OC)	MOQ	Order no.
	0.9 mm / 19 mm	■	20G x 3/4" Multi-adapter	80 mm	1 / 120 / 480	480	85.1637.205
	0.8 mm / 19 mm	■	21G x 3/4" Multi-adapter	80 mm	1 / 120 / 480	480	85.1638.205
	0.6 mm / 19 mm	■	23G x 3/4" Multi-adapter	80 mm	1 / 120 / 480	480	85.1640.205
	0.5 mm / 19 mm	■	25G x 3/4" Multi-adapter	80 mm	1 / 120 / 480	480	85.1642.205
	0.9 mm / 19 mm	■	20G x 3/4" Multi-adapter	200 mm	1 / 120 / 480	480	85.1637.235
	0.8 mm / 19 mm	■	21G x 3/4" Multi-adapter, screw-on guide sleeve	200 mm	1 / 120 / 480	480	85.1638.200
	0.8 mm / 19 mm	■	21G x 3/4" Multi-adapter	200 mm	1 / 120 / 480	480	85.1638.235
	0.6 mm / 19 mm	■	23G x 3/4" Multi-adapter	200 mm	1 / 120 / 480	480	85.1640.235
	0.5 mm / 19 mm	■	25G x 3/4" Multi-adapter	200 mm	1 / 120 / 480	480	85.1642.235

S-Monovette® Safety Needle

Illustration	Ø/Length	Color code	Description	Packaging (SP/IB/OC)	MOQ	Order no.
	0.9 mm/25 mm	■	20G x 1"	1 / 50 / 500	500	85.1372.200
	0.8 mm/25 mm	■	21G x 1"	1 / 50 / 500	500	85.1373.200
	0.7 mm/25 mm	■	22G x 1"	1 / 50 / 500	500	85.1441.200
	0.6 mm/25 mm	■	23G x 1"	1 / 50 / 500	500	85.1442.200
	0.9 mm/38 mm	■	20G x 1½"	1 / 50 / 500	500	85.1160.200
	0.8 mm/38 mm	■	21G x 1½"	1 / 50 / 500	500	85.1162.200
	0.8 mm/38 mm	■	21G x 1½"	1 / 50 / 500	500	85.1162.600 for S-Monovette® Metal Analysis 01.1604.400
	0.7 mm/38 mm	■	22G x 1½"	1 / 50 / 500	500	85.1440.200

Adapters


Illustration	No.	Description	Packaging (SP/IB/OC)	MOQ	Order no.
	1	Membrane adapter	1 / 100 / 1,000	1,000	14.1112
	2	Multi-Adapter Luer	1 / 100 / 1,000	1,000	14.1205
	3	Multi-Adapter Luer-Lock	1 / 100 / 1,000	1,000	14.1205.100
	4	Blood culture adapter LongNeck	1 / 50 / 250	250	14.1207
	5	Blood culture adapter Universal	1 / 20 / 100	100	14.1209

Caps

Screw caps for recapping S-Monovettes

Illustration	Screw cap outer Ø	Length without/ with cap x Ø	Packaging (SP/IB/OC)	MOQ	Order no.
	18 mm	92/108 x 15 mm	1,000 / 1,000 / 5,000	5,000	65.729
	18 mm	75/91 x 15 mm			
	16 mm	90/106 x 13 mm			
	16 mm	91/75 x 13 mm	1,000 / 1,000 / 5,000	5,000	65.728
	16 mm	81/65 x 13 mm			
	13 mm	108/92 x 11 mm			
	13 mm	82/66 x 11 mm	1,000 / 1,000 / 5,000	5,000	65.1121
	13 mm	146/130 x 8 mm			
	13 mm	82/66 x 8 mm			

Accessories

Illustration	Description	Packaging (SP/IB/OC)	MOQ	Order no.
	Haemo-Diff® with smear edge	500 / 500 / 2,000	2,000	14.1218
	Haemo-Diff® without smear edge	500 / 500 / 2,000	2,000	14.1217
	Seraplas® V11, suitable for S-Monovette® with Ø 11 mm	100 / 500 / 5,000	5,000	53.677
	Seraplas® V13, suitable for S-Monovette® with Ø 13 mm	100 / 500 / 1,000	1,000	53.419
	Seraplas® V15, suitable for S-Monovette® with Ø 15 mm	100 / 100 / 1,000	1,000	53.428
	Seraplas® V16, suitable for S-Monovette® with Ø 16 mm	100 / 100 / 1,000	1,000	53.420
	Color-coded cap, neon red	100 / 100 / 1,000	1,000	14.1513
	Color-coded cap, blue	100 / 100 / 1,000	1,000	14.1513.001
	Color-coded cap, yellow	100 / 100 / 1,000	1,000	14.1513.002
	Color-coded cap, neon green	100 / 100 / 1,000	1,000	14.1513.005
	Color-coded cap, grey	100 / 100 / 1,000	1,000	14.1513.008
	tournistrip® disposable paper tourniquet	200 / 200 / 4,800	2,000	95.1006
	Single-patient tourniquet TourniSTRETCH®	25 / 25 / 500	25	95.1008
	Reusable tourniquet	1 / 1 / 200	1	95.1007
	Reusable tourniquet, latex-free	1 / 1 / 200	1	95.1007.001

Micro needles & micro sample tubes

Micro Needle

Illustration	Description	Ø/Length	Color code	Packaging (SP/IB/OC)	MOQ	Order no.
	Micro Needle, 21 G x 3/4", No 2	0.8 mm/36 mm		1 / 50 / 500	50	85.921
	Micro Needle, 23 G x 3/4", No 16	0.6 mm/36 mm		1 / 50 / 500	50	85.923
	Micro Needle, 25 G x 3/4", No 18	0.5 mm/36 mm		1 / 50 / 500	50	85.925

Prepared micro sample tube with attached push cap and transparent label

Illustration	Preparation	Nominal volume	Ø/Length	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	1.3 ml	10.8 mm/40 mm		100 / 100 / 1,000	1,000	41.1501.105
	Serum Gel CAT	1.1 ml	10.8 mm/40 mm		100 / 100 / 1,000	1,000	41.1500.005
	Lithium Heparin LH	1.3 ml	10.8 mm/40 mm		100 / 100 / 1,000	1,000	41.1503.105
	EDTA K3E	1.3 ml	10.8 mm/40 mm		100 / 100 / 1,000	1,000	41.1504.105
	Citrate 9NC 3.2%	1.3 ml 1.0 ml 0.5 ml	10.8 mm/40 mm		100 / 100 / 1,000	1,000	41.1506.105 41.1506.100 41.1506.102

Prepared micro sample tube with screw cap and paper label

Illustration	Preparation	Nominal volume	Ø/Length	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	1.3 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.1392.105
	Serum Gel CAT	1.1 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.1378.005
	Lithium Heparin LH	1.3 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.1393.105
	EDTA K3E	1.3 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.1395.105
	EDTA K3E (Membrane in cap for direct adaptation)	1.3 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.3395.105
	Fluoride Heparin FH	1.3 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.1394.105
	Citrate 9NC 3.2%	1.3 ml	10.8 mm/45 mm		100 / 100 / 1,000	1,000	41.1350.105

Carrier tubes for Pediatric S-Monovette®, micro sample tubes and Microvette® and Multivette®

Illustration	Suitable for	Ø/Length	Color code	Packaging (SP/IB/OC)	MOQ	Order no.
	S-Monovette® with Ø 8 mm	54 x 11 mm	□	250 / 250 / 1,000	1,000	55.1570
	S-Monovette® with Ø 8 mm	54 x 11 mm	■	250 / 250 / 1,000	1,000	55.1570.002
	S-Monovette® with Ø 8 mm, Microvette® and micro tubes	65 x 13 mm	□	250 / 250 / 1,000	1,000	55.1571
	S-Monovette® with Ø 8 mm, Microvette® and micro tubes	65 x 13 mm	■	250 / 250 / 1,000	1,000	55.1571.051
	S-Monovette® with Ø 8 mm	55 x 13 mm	□	250 / 250 / 1,000	1,000	55.1572
	Microvette® and micro tube with screw cap	75 x 13 mm	□	250 / 250 / 2,000	2,000	55.525.003
	Microvette® and micro tube with screw cap	75 x 13 mm	■	250 / 250 / 2,000	2,000	55.525.002
	Multivette® 600	60 x 11.5 mm	□	250 / 250 / 1,000	1,000	55.1577
	Multivette® 600	60 x 11.5 mm	■	250 / 250 / 1,000	1,000	55.1577.002

Microvette® APT & Microvette® 100/200

Microvette® APT with round bottom inner tube

Illustration	Preparation	Nominal volume	Length/∅	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	EDTA K2E	250 µl	75 mm x 13 mm	■	50 / 50 / 500	500	20.1331.100
	EDTA K2E	500 µl	75 mm x 13 mm	■	100 / 100 / 500	500	20.1330.100

Microvette® 100 with conical inner tube

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	100 µl	■	100 / 100 / 1,000	1,000	20.1280.100
	Lithium Heparin LH	100 µl	■	100 / 100 / 500	500	20.1282.100

Microvette® 100 with round bottom inner tube

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	EDTA K3E	100 µl	■	100 / 100 / 1,000	1,000	20.1278.100

Microvette® 200 with conical inner tube

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	200 µl	■	100 / 100 / 1,000	1,000	20.1290.100
	Lithium Heparin LH	200 µl	■	100 / 100 / 1,000	1,000	20.1292.100
	Fluoride Heparin FH	200 µl	■	100 / 100 / 1,000	1,000	20.1289.100

Microvette® 200 with round bottom inner tube


Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum Gel CAT	200 µl	■	100 / 100 / 1,000	1,000	20.1291
	EDTA K3E	200 µl	■	100 / 100 / 1,000	1,000	20.1288.100

Microvette® 300/500 & Microvette® CB 200 BSG

Microvette® 300 with conical inner tube

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	300 µl	■	100 / 100 / 1,000	1,000	20.1308.100
	Lithium Heparin LH	300 µl	■	100 / 100 / 1,000	1,000	20.1309.100
	Fluoride Heparin FH	300 µl	■	100 / 100 / 1,000	1,000	20.1307.100

Microvette® 500 with round bottom inner tube

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	500 µl	■	100 / 100 / 1,000	1,000	20.1343.100
	Serum Gel CAT	500 µl	■	100 / 100 / 1,000	1,000	20.1344
	Lithium Heparin LH	500 µl	■	100 / 100 / 1,000	1,000	20.1345.100
	Lithium Heparin LH	500 µl	■	100 / 100 / 1,000	1,000	20.1346.100
	EDTA K3E	500 µl	■	100 / 100 / 1,000	1,000	20.1341.100
	EDTA K3E	500 µl	■	100 / 100 / 1,000	1,000	20.1341.102*
	EDTA K2E	500 µl	■	100 / 100 / 1,000	1,000	20.1339.100

* for lead determination

Microvette® CB 200 ESR

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Citrate	200 µl	■	50 / 50 / 500	500	18.1325.100

Sedimentation rack for Microvette® CB 200 ESR

Illustration	Description	Packaging (SP/IB/OC)	Order no.
	ESR rack with graduated back wall	1 / 1 / 1	90.1091

Multivette® 600 & Safety Lancet

Multivette® 600

Illustration	Preparation	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Serum CAT	600 µl	■	100 / 100 / 1,000	1,000	15.1670.100
	Serum Gel CAT	600 µl	■	100 / 100 / 1,000	1,000	15.1674
	Lithium Heparin LH	600 µl	■	100 / 100 / 1,000	1,000	15.1673.100
	Lithium Heparin Gel LH	600 µl	■	100 / 100 / 1,000	1,000	15.1675.100
	EDTA K3E	600 µl	■	100 / 100 / 1,000	1,000	15.1671.100

Safety Lancet

Illustration	Description	Needle size/ Incision length	Penetration depth	Color code	Packaging (SP/IB/OC)	MOQ	Order no.
	Safety Lancet Mini	28 G	1.6 mm	■	200 / 200 / 1,200	1,200	85.1015
	Safety Lancet Normal	21 G	1.8 mm	■	200 / 200 / 1,200	1,200	85.1016
	Safety Lancet Extra	18 G	1.8 mm	■	200 / 200 / 1,200	1,200	85.1017
	Safety Lancet Super	Blade 1.5 mm	1.6 mm	■	200 / 200 / 1,200	1,200	85.1018
	Safety Lancet Neonatal	Blade 1.5 mm	1.2 mm	■	200 / 200 / 1,200	1,200	85.1019

Incision Lancet - Safety-Heel®

Illustration	Description	Incision length	Penetration depth	Color code	Packaging (SP/IB/OC)	MOQ	Order no.
	Premie	1.75 mm	0.85 mm	■	50 / 50 / 200	200	85.1025
	Newborn	2.5 mm	1.0 mm	■	50 / 50 / 200	200	85.1026

Blood gas analysis

Blood Gas Monovette®

Illustration	Description	Length x Ø	Nominal volume	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Blood Gas Monovette®	66 x 11 mm	1 ml	□/■	50 / 50 / 500	500	05.1146
	Blood Gas Monovette®	66 x 11 mm	1 ml	□/■	1 / 100 / 500	500	05.1146.020*
	Blood Gas Monovette®	66 x 11 mm	2 ml	□/■	50 / 50 / 500	500	05.1147
	Blood Gas Monovette®	66 x 11 mm	2 ml	□/■	1 / 100 / 500	500	05.1147.020*
	Blood Gas Monovette®, membrane adapter assembled	66 x 11 mm	2 ml	□/■	1 / 100 / 500	500	05.1147.028
	Blood Gas, ventilator assembled	66 x 11 mm	2 ml	□/■	1 / 100 / 500	500	05.1147.048
	Ventilator	66 x 11 mm		☒/■	500 / 500 / 1,000	1.000	14.1148

* individually wrapped

Plastic blood gas capillaries

Illustration	Description	Nominal volume	Packaging (SP/IB/OC)	MOQ	Order no.
	Plastic blood gas capillaries	100µl	250 / 250 / 16,000	250	19.930.100
	Plastic blood gas capillaries	125µl	250 / 250 / 16,000	250	19.930.125
	Plastic blood gas capillaries	140µl	200 / 200 / 12,800	200	19.930.140
	Plastic blood gas capillaries	175µl	200 / 200 / 12,800	200	19.930.175

Plastic blood gas capillary sets

Illustration	Description	Nominal volume	Packaging (SP/IB/OC)	MOQ	Order no.
	Blood gas capillary set, caps, and mixing wires	100µl	1 / 1 / 16	1	51.931.100
	Blood gas capillary set, caps, and mixing wires	125µl	1 / 1 / 16	1	51.931.125
	Blood gas capillary set, caps, and mixing wires	140µl	1 / 1 / 16	1	51.931.140
	Blood gas capillary set, caps, and mixing wires	175µl	1 / 1 / 16	1	51.931.175

Accessories

Illustration	Description	Cap color	Packaging (SP/IB/OC)	MOQ	Order no.
	Closure cap, OD 2.05 mm, for 100µl / 125µl capillaries	□	500 / 500 / 32,000	500	65.935.205
	Closure cap, OD 2.3 mm, for 140µl / 175µl capillaries	■	400 / 400 / 32,000	400	65.935.230
	Mixing wires, Ø 1 mm, length 9 mm	-	250 / 250 / 41,250	250	95.936
	Bell-shaped magnet, red	■	10 / 10 / 160	10	95.937

Notes

[illegible]

Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Legend


Packaging: SP: Smallest sub-packaging of a product
IB: Inner box, the SP is packaged in the IB
OC: Outer case, the outer case is usually also the minimum order quantity
MOQ: Minimum order quantity

If you have any questions,
we'll be happy to help!

Visit our website: www.sarstedt.com

Pre-analytical workflow powered by SARSTEDT

Take advantage of the synergy
of our synchronized systems.



SARSTEDT, Inc.

1025 St. James Church Road
P.O. Box 468
Newton, NC 28658 - 0468

Toll free: +1 800 257 5101

customerservice.us@sarstedt.com
www.sarstedt.com

Discover the 360°
solutions for your
pre-analytics from
SARSTEDT



workflow.sarstedt.com