

# 5 ml Micro Tubes

The intermediate size for laboratories



The micro tube for medium volumes:

- Easy handling of sample volumes between 2 and 5 ml
- Reliable SafeSeal locking closure
- PCR Performance Tested quality for purity requirements
- Biosphere® Plus for maximum protection against contamination



# 5 ml micro tubes

## Reliable handling of sample volumes between 2 and 5 ml

Routine applications in the laboratory such as dilution series, cell lysis or sample storage often require a filling volume greater than 2 ml. Larger volume tubes that are often narrow and have screw caps are not always suitable for these applications.

The SARSTEDT 5 ml micro tube combines the advantages of a SafeSeal reaction tube and a 5 ml screw cap tube in one. Easy to open with one hand thanks to the ergonomic design, the snap-on lid with a sealing lip and holding tab closes securely even under thermal stress. The wide hinge allows for easy and precise closure, and cap mix-ups associated with screw caps are eliminated.



## Advantages of the 5 ml micro tube

- Robust polypropylene – can be centrifuged up to 21,000 x g, autoclavable
- Molded graduation – easy to check volume
- Writing space and frosted cap – easy marking for quick identification
- Low height – minimizes contamination risk from pipette cone contact
- Compact size – space-saving storage of samples, media or cultures

## Ordering information

Order no.	Name	Packaging
72.701	5 ml micro tube	100 / bag 800 / case
72.701.400	5 ml micro tube, PCR Performance Tested	100 / bag 800 / case
95.1298	5 ml micro tube centrifuge adapter	8 / bag 96 / case



## Accessories

Order no.	Name	Format	External dimensions (mm)	Packaging
93.850.001	Polycarbonate rack for 18 x 5 ml micro tubes	6 x 3	137 x 70 x 40 (L x W x H)	1 / box
93.1431.001	Polycarbonate rack for 48 x 5 ml micro tubes	12 x 4	257 x 90 x 40 (L x W x H)	1 / box

Technical modifications reserved

This document may contain information on products that may not be available in particular countries

SARSTEDT, Inc.  
1025 St. James Church Road  
P.O. Box 468  
Newton, NC 28658 - 0468  
Toll free: +1 800 257 5101  
Fax: +1 828 465 4003  
customerservice.us@sarstedt.com  
www.sarstedt.com

