

Ordering information

10 μ l PIPETTE TIPS



SMALL
AND FINE

Precise **DE-M**

10 µl pipette tips are suitable for precise pipetting of even the smallest volumes of liquid thanks to their compact size and their resulting small dead volume.

Innovative

The tip is suitable for any application thanks to its sophisticated new packaging system. Our system adapts to your needs. Whether it's ultra-clean working with the SingleRefill, high throughput with the StackPack or simply a filled box – there is a suitable solution for any requirement.

TOP TIP Identification at a glance! In addition to the familiar colour code of the trays, you will now also find the most important information about the tip directly on each tray.



Ordering information

Vol.	Order no	Description	Filter	Purity	Design	Packaging	Units per box
10 µl	70.3010.275	Low retention filter tip, 10 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Low retention		1,920
	70.3010.375	Low retention filter tip, 10 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Low retention		1,920
	70.3010.285	Low retention filter tip, 2.5 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Low retention		1,920
	70.3010.385	Low retention filter tip, 2.5 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Low retention		1,920
	70.3010.255	Filter tip, 10 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Standard		1,920
	70.3010.355	Filter tip, 10 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Standard		1,920
	70.3010.265	Filter tip, 2.5 µl, Biosphere® plus, transparent, 96 pcs/Box	+		Standard		1,920
	70.3010.365	Filter tip, 2.5 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	+		Standard		1,920
	70.3010.205	Pipette tip, 10 µl, Biosphere® plus, transparent, 96 pcs/Box	-		Standard		1,920
	70.3010.305	Pipette tip, 10 µl, Biosphere® plus, transparent, 96 pcs/SingleRefill	-		Standard		1,920
	70.3010.210	Low retention pipette tip, 10 µl, PCR performance tested, transparent, 96 pcs/Box	-		Low retention		1,920
	70.3010.110	Low retention pipette tip, 10 µl, PCR performance tested, transparent, 5 x 96 pcs/StackPack	-		Low retention		1,920
	70.3010.200	Pipette tip, 10 µl, PCR performance tested, transparent, 96 pcs/box	-		Standard		1,920
	70.3010.100	Pipette tip, 10 µl, PCR performance tested, transparent, 5 x 96 pcs/StackPack	-		Standard		1,920
	70.3010	Pipette tip, 10 µl, PCR Performance Tested, transparent, 1000 pcs/bag	-		Standard		10,000

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com



Further information
can be found at
www.refillrevolution.tips