

# S-Monovette® Serum Gel LightPROTECT

Light protection –  
from blood collection to archiving



- Light protection from the outset
- Light protection across the UV range
- More sustainability at the workplace by avoiding aluminium foil
- Usual convenient specimen identification as labels are always legible
- Process-compliant pre-analytics

The S-Monovette® Serum Gel LightPROTECT eliminates the need to wrap the blood sample in aluminium foil. Direct legibility and thus usual convenient specimen identification are maintained. This also facilitates an effortless contribution to more sustainability at the workplace by conserving resources (aluminium foil).

Light protection is guaranteed in the relevant UV range (190 nm – 380 nm). Average light permeability in this range is less than 1%.




Fig. 1: Transmission of the S-Monovette® Serum Gel LightPROTECT tube wall



Art. no.	Preparation	Volume	Length / Ø	Case
04.1925.003	S-Monovette® Serum Gel LightPROTECT	4.0 ml	75x13 mm	500 / box 50 / inner box

SARSTEDT AG & Co. KG  
 Sarstedtstraße 1  
 D-51588 Nümbrecht  
 Phone: +49 2293 305 0  
 Fax: +49 2293 305 3992  
 export@sarstedt.com  
 www.sarstedt.com