



# Pre-barcode order overview






















## EU colour code

Preparation	Order no. Pre-Barcode	Prefix	Volume (ml)	Length/Ø	Colour code Cap/label
 Serum	04.1934.088	001	4.9	90mm x 13mm	
	04.1943.088	003	2.7	75mm x 13mm	
	01.1601.088	005	7.5	92mm x 15mm	
 Serum Gel	04.1923.088	007	2.7	75mm x 13mm	
	04.1925.088	008	4.0	75mm x 13mm	
	04.1935.088	009	4.9	90 mm x 13mm	
	01.1602.088	010	7.5	92mm x 15mm	
 K3 EDTA	04.1901.088	011	2.6	65mm x 13mm	
	04.1951.088	013	1.8	65mm x 13mm	
	04.1917.088	049	2.7	75mm x 13mm	
	04.1921.088	051	4.0	75mm x 13mm	
 Citrate 1:10	04.1902.088	015	2.9	65mm x 13mm	
	04.1919.088	017	3.0	75mm x 13mm	
	04.1955.088	019	1.8	75mm x 13mm	
 Fluoride	04.1918.088	021	2.7	75mm x 13mm	
 Fluoride/citrate	04.1945.288	023	3.1	75mm x 13mm	
 Lithium heparin	04.1929.088	024	2.7	75mm x 13mm	
	04.1936.088	026	4.9	90mm x 13mm	
	01.1604.088	032	7.5	92mm x 15mm	
 Lithium Heparin Gel	01.1634.088	036	7.5	92mm x 15mm	
	04.1940.088	038	4.9	90mm x 13mm	
	04.1928.088	040	2.7	75mm x 13mm	
 Lithium Heparin Gel*	04.1953.088	028	4.0	75mm x 13mm	
	04.1954.088	030	4.9	90mm x 13mm	
	04.1952.088	034	2.7	75mm x 13mm	
 Urine Monovette® without preparation	10.252.088	042	10.0	102 mm x 15mm	
	10.252.388*	043	10.0	102 mm x 15mm	
	10.258.088	046	8.5	92mm x 15mm	
	10.258.388*	047	8.5	92mm x 15mm	
	10.250.088	048	3.2	75mm x 13mm	
 Urine Monovette® with boric acid	10.253.088	044	10.0	102 mm x 15mm	
	10.253.388*	045	10.0	102 mm x 15mm	

\* individually sterile-packed

## ISO colour code

Preparation	Order no. Pre-Barcode	Prefix	Volume (ml)	Length/Ø	Colour code Cap/label
 Serum	04.1934.288	002	4.9	90mm x 13mm	
	04.1943.288	004	2.7	75mm x 13mm	
	01.1601.288	006	7.5	92mm x 15mm	
 Serum Gel	04.1923.088	007	2.7	75mm x 13mm	
	04.1925.088	008	4.0	75mm x 13mm	
	04.1935.088	009	4.9	90 mm x 13mm	
	01.1602.088	010	7.5	92mm x 15mm	
 K3 EDTA	04.1901.288	012	2.6	65mm x 13mm	
	04.1951.288	014	1.8	65mm x 13mm	
	04.1917.288	050	2.7	75mm x 13mm	
	04.1921.288	052	4.0	75mm x 13mm	
 Citrate 1:10	04.1902.288	016	2.9	65mm x 13mm	
	04.1919.288	018	3.0	75mm x 13mm	
	04.1955.288	020	1.8	75mm x 13mm	
 Fluoride	04.1918.288	022	2.7	75mm x 13mm	
 Fluoride/citrate	04.1945.288	023	3.1	75mm x 13mm	
 Lithium heparin	04.1929.288	025	2.7	75mm x 13mm	
	04.1936.288	027	4.9	90mm x 13 mm	
	01.1604.288	033	7.5	92mm x 15 mm	
 Lithium Heparin Gel	01.1634.288	037	7.5	92mm x 15mm	
	04.1940.288	039	4.9	90mm x 13mm	
	04.1928.288	041	2.7	75mm x 13mm	
 Lithium Heparin Gel*	04.1954.288	031	4.9	90mm x 13mm	
	04.1953.288	029	4.0	75mm x 13mm	
	04.1952.288	035	2.7	75mm x 13mm	
 Urine Monovette® without preparation	10.252.088	042	10.0	102 mm x 15mm	
	10.252.388*	043	10.0	102 mm x 15mm	
	10.258.088	046	8.5	92mm x 15mm	
	10.258.388*	047	8.5	92mm x 15mm	
	10.250.088	048	3.2	75mm x 13mm	
 Urine Monovette® with boric acid	10.253.088	044	10.0	102 mm x 15mm	
	10.253.388*	045	10.0	102 mm x 15mm	

\* individually sterile-packed

## SARSTEDT B.V.

Rithmeesterpark 23 A  
4838 GZ Breda  
Tel +31 76 501 75 50  
Fax +31 76 501 76 26  
info.nl@sarstedt.com

## SARSTEDT BV

Uitbreidingsstraat 84/3  
2600 Berchem  
Tel +32 3 541 76 92  
Fax +32 3 541 81 03  
info.be@sarstedt.com

www.sarstedt.com

# Structure of the unique barcode

The diagram shows a label for 'S-Monovette® Serum Gel CAT/4.9ml'. The label features a barcode with the number '00901000001' printed below it. Below the barcode, there are several fields: 'REF 04.1935.088', 'STERILE R', 'IVD CE', and 'SARSTEDT' (vertical text). Below these, there are 'LOT 1234567' and 'XXXX-XX' (expiry date). The label also has a small 'F404030' code at the bottom right.

**Prefix** for unambiguous differentiation of the sample tubes

**Unique barcode numbering** for unambiguous assignment to the patient order

Imprinted **article number, batch** and **expiry date** for visual inspection



Supplementary product information such as article number, name, preparation, batch or expiry date are linked to the barcode. This product information is available digitally through the SARSTEDT WebService.