


For the use of SARSTEDT consumables in cryopreservation

Cryo Performance Tested items fulfil the following certified purity criteria:

- ✓ Sterile
- ✓ Pyrogen-free/endotoxin-free
- ✓ DNA-free
- ✓ DNase/RNase-free
- ✓ Non-cytotoxic
- ✓ Non-mutagenic
- ✓ CE 

We guarantee that the following limit values are met:

Sterility validated according to ISO 11137

Pyrogens/endotoxins <0.06 EU/ml

Human DNA <0.5 pg/μl, bacterial DNA <0.02 pg/μl

DNase <1x10⁻⁵ U/μl, RNase <1x10⁻⁹ Kunitz-units/μl

Non-cytotoxic according to ISO 10993

Non-mutagenic according to the Ames Test II