

Nombre de puits	96 prémontés	96	96	96	96	384	96	384
Demi-Juue/Juue intégrale	sans	sans	Demi	Intégrale	Demi	Intégrale	Demi	Intégrale
Profile	High	High	High	Low	Low	-	« Lightcycler »	« Lightcycler »
Réf. plaques de PCR	72.985	72.1978 72.1978.010	72.1979 72.1979.010 72.1979.003 72.1979.201 72.1979.700 72.1979.102 72.1979.132	72.1980 72.1980.010 72.1980.201 72.1980.600 72.1980.700	72.1981 72.1981.010	72.1984.202	72.1982.202	72.1985.202

BioRad® / MJ Research®

CFX96 Touch™ PCR en temps réel				●				
CFX384 Touch™ PCR en temps réel						●		
Système d'automatisation CFX II								
Cycleur thermique T100™	●	●		●	●			
Cycleur thermique S1000™		●	●	●		●		
Cycleur thermique C1000Touch™		●	●	●	●	●		
Cycleur thermique iCycler IQ™	●	●	●					
Cycleur thermique IQ4™	●	●	●					
Cycleur thermique IQ5™	●	●	●					
Cycleur thermique MyCycler™	●	●	●					
Chromo4™		●		●				
Opticon™, Opticon2™				●				
BaseStation™				●				

Corbett Research® / Qiagen®

Palm Cyclyer 96 puits		●						
Palm Cyclyer 384 puits				●				

Eppendorf®

Mastercycler® nexus	●	●						
Mastercycler® ep realplex		●	●			●		
Mastercycler® gradient	●	●		●				
Mastercycler® ep gradient	●	●	●	●				
Mastercycler® pro	●	●		●				

Ericom®

Deltacycler			●	●				
SingleBlock			●	●				
TwinBlock			●	●				

MWG®

Primus 96 puits		●		●				
Primus 384 puits						●		
Q-Lifecycler		●	●	●				

Roche®

Système Lightcycler® 96							●	
Système Lightcycler® 480							●	●

Stratagene® / Agilent®

Système AriaMx PCR en temps réel				●	●			
Mx3000P™	●	●						
Mx3005P™	●	●	●	●				
Mx4000™	●	●	●					
Cycleur à gradient		●		●				
Robocycler 384 puits						●		

Teche®

Cyclogene		●		●				
Flexigene		●	●	●				
Genius / Genius Quad		●	●	●				
OMN-E		●	●	●				
PCR Express	●	●	●			●		
Primus 96		●						
Px2 / PxE		●	●			●		
Quantica			●	●				
TC412 / TC512		●		●		●		
Touchgene / Touchgene Gradient		●	●	●	●			