

Número de pocillos	96 preinsertada	96	96	96	96	384	96	384
Medio faldón / faldón	sin	sin	medio	con	medio	con	medio	con
Perfil	alto	alto	alto	bajo	bajo	-	«Lightcycler»	«Lightcycler»
Referencia placas de PCR	72.985	72.1978 72.1978.010	72.1979 72.1979.010 72.1979.003 72.1979.201 72.1979.700 72.1979.102 72.1979.132	72.1980 72.1980.010 72.1980.201 72.1980.600 72.1980.700	72.1981 72.1981.010	72.1984.202	72.1982.202	72.1985.202

#### BioRad® / MJ Research®

CFX96 Touch™ Real-Time PCR				●				
CFX384 Touch™ Real-Time PCR						●		
CFX Automation System II								
T100™ Thermal Cycler	●	●		●	●			
S1000™ Thermal Cycler		●	●	●		●		
C1000Touch™ Thermal Cycler		●	●	●	●	●		
iCycler iQ™ Thermal Cycler	●	●	●					
iQ4™ Thermal Cycler	●	●	●					
iQ5™ Thermal Cycler	●	●	●					
MyCycler™ Thermal Cycler	●	●	●					
Chromo4™		●		●				
Opticon™, Opticon2™				●				
BaseStation™				●				

#### Corbett Research® / Qiagen®

Palm Cycler 96 pocillos		●						
Palm Cycler 384 pocillos				●				

#### Eppendorf®

Mastercycler® nexus	●	●						
Mastercycler® ep realplex		●	●			●		
Mastercycler® gradient	●	●		●				
Mastercycler® ep gradient	●	●	●	●				
Mastercycler® pro	●	●		●				

#### Ericom®

Deltacycler			●	●				
SingleBlock			●	●				
TwinBlock			●	●				

#### MWG®

Primus 96 pocillos		●		●				
Primus 384 pocillos						●		
The Q-Lifecycler		●	●	●				

#### Roche®

Lightcycler® 96 System							●	
Lightcycler® 480 System							●	●

#### Stratagene® / Agilent®

AriaMx Real-Time PCR System				●	●			
Mx3000P™	●	●						
Mx3005P™	●	●	●	●				
Mx4000™	●	●	●					
Gradient Cycler		●		●				
Robocycler 384 pocillos							●	

#### Techné®

Cyclogene		●		●				
Flexigene		●	●	●				
Genius / Genius Quad		●	●	●				
OMN-E		●	●	●				
PCR Express	●	●	●				●	
Primus 96		●						
Px2 / PxE		●	●				●	
Quantica			●	●				
TC412 / TC512		●		●			●	
Touchgene / Touchgene Gradient		●	●	●	●			