

# Product Specification Multivette<sup>®</sup> 600 Lithium heparin gel LH, 600 μl, cap orange, screw cap



### **Product description**

Order number	
Product description	

#### 15.1675

Multivette<sup>®</sup> 600 Lithium heparin gel LH, venous and capillary blood samples, preparation: Lithium heparin gel, nominal volume: 600  $\mu$ l, (LxØ) with cap: 53.3 x 8 mm, with print, cap orange, screw cap, flat base, 100 piece(s)/case

## **Product characteristics**

Type of collection	Venous and capillary
Type of preparation	Lithium heparin gel
Colour code	EU
Label/ Print	with print
Colour of print/label	orange
Сар	standard screw cap
Closure type	screw cap
Application	venous and capillary blood samples
Base shape	flat base
Recommended centrifugation	10000 x g
Centrifugation time	5 Min
Application Base shape Recommended centrifugation	venous and capillary blood samples flat base 10000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This document was prepared by EDP support and is valid without signature.





Page 2	2
--------	---

## Size

Sample volume	600 µl
Diameter	8 mm
Length including cap	53.3 mm
Length excluding cap	49.8 mm
Material & colours	
Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Purity & certification	
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes
Packaging	
Minimum order qty. Type of smallest subpackaging Piece(s) / inner box Piece (s) / outer case Piece(s) / pallet Depth of box Width of box Height of box Depth of case Width of case	1000 case 100 1000 63000 201 mm 81 mm 81 mm 423 mm 213 mm
Height of case	185 mm
Case volume	0.0167 cbm
Weight of product	0.0024 kg
Weight of case	3.03 kg
EAN of inner box	4038917087438
EAN case	4038917005869

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

