

## Product Specification Micro sample tube Serum Gel CAT, 1.1 ml, screw cap, EU/ISO

Page 1



## **Product description**

Order number 41.1378.005

Product description Micro sample tube Serum Gel CAT, venous blood

collection, preparation: Clotting activator/gel, nominal volume: 1.1 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap,

cap: brown, colour code EU/ISO, with paper label,

label/print: brown, skirted conical base, 100 piece(s)/bag

## **Product characteristics**

Type of collection venous

Type of preparation Clotting activator/gel

Colour code EU/ISO

Label/ Print with paper label

Colour of print/label brown Graduation yes

Closure type screw cap

Application venous blood collection
Base shape skirted conical base

Volume of work 1,1 ml Recommended centrifugation 10000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification Micro sample tube Serum Gel CAT, 1.1 ml, screw cap, EU/ISO

Page 2

Centrifugation time 5 Min

Size

Sample volume 1.1 ml
Diameter 10.8 mm
Length of product 47 mm
Length including cap 47 mm
Length excluding cap 45 mm

**Material & colours** 

Colour of product transparent

Closure material Polypropylene (PP)

Colour of cap brown

Tube material Polypropylene (PP)

**Purity & certification** 

Satisfies the requirement IATA, ADR

Product category In-vitro diagnostic (IVD)

CE certified CE - manufacturer's self-declaration

Batched yes

**Packaging** 

Minimum order qty. 1000 Type of smallest bag

subpackaging

Piece(s) / inner box 100 Piece (s) / outer case 1000 96000 Piece(s) / pallet Depth of case 348 mm Width of case 278 mm Height of case 135 mm Case volume 0.0131 cbm Weight of product  $0.0019 \, kg$ Weight of case 2.17 kg

EAN of inner box 4038917127967 EAN case 4038917021500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

