



### Product description

Order number	72.730.106
Descrição do artigo	Microtubo com tampa de rosca, volume de trabalho: 0,5 ml, fundo cônico com base plana, com recartilhado, transparente, tampa: natural, tampa montada acoplada, com campo de escrita impresso, com escala, estéril, 100 unid./pacote

### Product characteristics

Rótulo/impressão	com impressão
Cor impressão/rótulo	branca
Graduação	sim
Tampa	tampa com anel
Ranhura	sim
Forma do fundo	fundo cônico com base plana
Volume de trabalho	0,5 ml
Centrifugação máxima (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Diâmetro	10.8 mm
Comprimento do produto	44 mm

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	Polypropylene (PP)
Cor da tampa	natural
Material do O-ring	Borracha de etileno-propileno-dieno (EPDM)

### Purity & certification

Corresponde às exigências	IATA
Categoria do produto	Diagnóstico in-vitro
Certificação CE	CE
Padrão de pureza	estéril
Esterilização	Irradiação de elétrons
Em lotes	sim

### Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	96000
Profundidade da embalagem externa	348 mm
Largura da embalagem externa	278 mm
Altura da embalagem externa	140 mm
Volume da embalagem externa	0.0135 cbm
Peso do produto	0.0016 kg
Peso da embalagem externa	1.9 kg
EAN da caixa interna	4038917119481
EAN da embalagem externa	4038917016674

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

