



### Product description

Order number	62.554.502
Descrição do artigo	Tubo de rosca, volume de trabalho: 15 ml, (CxØ): 120 x 17 mm, material: PP, fundo cônico, transparente, HD-PE, tampa: vermelha, tampa montada, com impressão, rótulo/impressão: branca/azul, com escala, livre de DNA/DNase/RNase, livre de pirogênio/endotoxinas, não citotóxico, estéril, 50 unid./pacote

### Product characteristics

Rótulo/impressão	com impressão
Cor impressão/rótulo	branca/azul
Graduação	sim
Tampa	tampa de rosca padrão
Tipo de fechamento	tampa de rosca
Forma do fundo	fundo cônico
Volume de trabalho	15 ml
Centrifugação máxima (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Diâmetro	17 mm
Comprimento do produto	120 mm

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	vermelha

### Purity & certification

Corresponde às exigências	IATA
Categoria do produto	Diagnóstico in-vitro
Certificação CE	CE
Padrão de pureza	livre de DNA/DNase/RNase, livre de pirogênio/endotoxinas, não citotóxico, estéril
Esterilização	Irradiação de elétrons
Isenção de pirogênio	sim
Em lotes	sim

### Packaging

Quantidade mínima do pedido	500
Tipo de embalagem interior	pacote
Unidades por caixa interna	50
Unidades por caixa	500
Unidades por palete	20000
Profundidade da embalagem externa	598 mm
Largura da embalagem externa	398 mm
Altura da embalagem externa	165 mm
Volume da embalagem externa	0.0393 cbm
Peso do produto	0.0063 kg
Peso da embalagem externa	3.74 kg
EAN da caixa interna	4038917113182
EAN da embalagem externa	4038917013932

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

