



## Product Specification

**Blood Gas Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)**

**Page 1**



### Product description

Order number	05.1146.020
Product description	Blood Gas Monovette® , preparation: calcium-balanced lithium heparin, nominal volume: 1 ml, cap white/orange, label/print: transparent/orange, connection: Luer (m), (LxØ) without cap: 66 x 11 mm, 1 piece(s)/blister, sterile, 100 piece(s)/case

### Product characteristics

Type of collection	Arterial and venous
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	50 IU/ml blood
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Application	blood gas analysis
Connection	Luer (m)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Blood Gas Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)**

**Page 2**

### Size

Sample volume	1 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	100 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white/orange
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of box	351 mm
Width of box	111 mm
Height of box	171 mm
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**

**Blood Gas Monovette® , calcium-balanced lithium  
heparin, 1 ml, cap white/orange, connection: Luer  
(m)**

Case volume	0.0415 cbm
Weight of product	0.0043 kg
Weight of case	3.42 kg
EAN of inner box	4038917247252
EAN case	4038917247245

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

