



**Product Specification**  
**Blood gas capillaries, calcium-balanced lithium  
heparin, 125 µl**

**Page 1**



**Product description**

Order number	19.930.125
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container

**Product characteristics**

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

**Size**

Sample volume	125 µl
Diameter	2.05 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	container
Piece(s) / inner box	250
Piece (s) / outer case	16000
Piece(s) / pallet	1280000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0.0196 cbm
Weight of product	0.0003 kg
Weight of case	6.64 kg
EAN of inner box	4038917234283
EAN case	4038917234276

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

