



**Product Specification**  
**Liquor/CSF-Monovette® , 10 ml, cap blue, (LxØ):**  
**102 x 15 mm, 1 piece(s)/blister**

**Page 1**

**Product description**

Order number	10.252.027
Product description	Liquor/CSF-Monovette® , nominal volume: 10 ml, without preparation, collection and storage of diagnostic specimens, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap blue, with paper label, "Liquor/Cerebrospinal Fluid", sterile, 1 piece(s)/blister

**Product characteristics**

Type of preparation	without preparation
Label/ Print	with paper label
Colour of print/label	blue
Graduation	yes
Cap	Luer screw cap
Closure type	screw cap
Light protection	no
Application	collection and storage of diagnostic specimens
Example of use	CSF, collection, transport and storage of CSF/liquor
Base shape	round base
Volume of work	10 ml
Connection	Luer (m)

**Size**

Sample volume	10 ml
Diameter	15 mm
Length including cap	117 mm
Length excluding cap	102 mm
Length including cap and plunger	127 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Liquor/CSF-Monovette® , 10 ml, cap blue, (LxØ):**  
**102 x 15 mm, 1 piece(s)/blister**

**Page 2**

Piston plunger material                      Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0077 kg
Weight of case	5.46 kg
EAN of inner box	4038917104487
EAN case	4038917040334

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**