




**Product Specification**  
**Absorbent liner, suitable for Mailing container 114  
x 44 mm, (LxW): 135 x 95 mm**

**Page 1**



**Product description**

Order number	78.669
Product description	Absorbent liner, suitable for Mailing container 114 x 44 mm, (LxW): 135 x 95 mm, 100 piece(s)/case

**Product characteristics**

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	Mailing container 114 x 44 mm

**Size**

Width of product	95 mm
Length of product	135 mm

**Material & colours**

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for Mailing container 114**  
**x 44 mm, (LxW): 135 x 95 mm**

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	385 mm
Width of case	385 mm
Height of case	160 mm
Case volume	0.0237 cbm
Weight of product	0.0026 kg
Weight of case	3.04 kg
EAN of inner box	4038917083249
EAN case	4038917002196

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

