



**Product Specification**  
**Absorbent liner, suitable for Mailing container 126  
x 30 mm, (LxW): 75 x 90 mm**

**Page 1**



**Product description**

Order number 78.670  
Product description Absorbent liner, suitable for Mailing container 126 x 30 mm,  
(LxW): 75 x 90 mm, 1,000 piece(s)/case

**Product characteristics**

Application Mailing of samples, Multi-Safe disposel boxes  
Suitable for Mailing container 126 x 30 mm

**Size**

Width of product 90 mm  
Length of product 75 mm

**Material & colours**

Product material cellulose

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for Mailing container 126**  
**x 30 mm, (LxW): 75 x 90 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      4000  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      4000  
Piece(s) / pallet                              72000  
Depth of case                              515 mm  
Width of case                                360 mm  
Height of case                               360 mm  
Case volume                                0.0667 cbm  
Weight of product                        0.001 kg  
Weight of case                              4.36 kg  
EAN of inner box                        4038917377492  
EAN case                                    4038917159937

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

