



**Product Specification**  
**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, transparent**

**Page 1**



**Product description**

Order number	75.562.105
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: white, with safety label, cap assembled, sterile, 5 piece(s)/bag

**Product characteristics**

Label/ Print	with safety label
Colour of print/label	white/red
Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Volume of work	100 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, transparent**

**Page 2**

**Size**

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	200
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	200
Piece(s) / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0155 kg
Weight of case	3.78 kg
EAN of inner box	4038917113632
EAN case	4038917014052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

