



Product Specification
Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent

Page 1



Product description

Order number	75.9922.744
Product description	Multi-purpose container, max. working volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, 250 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Base shape	flat base
Volume of work	70 ml

Size

Sample volume	70 ml
Max. volume	70 ml
Diameter	44 mm
Diameter of opening	44 mm
Height of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent

Page 2

Length of product 55 mm
Length excluding cap 55 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap yellow

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 500
Piece(s) / pallet 8000
Depth of case 548 mm
Width of case 408 mm
Height of case 373 mm
Case volume 0.0834 cbm
Weight of product 0.0123 kg
Weight of case 7 kg
EAN of inner box 4038917152754
EAN case 4038917027793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT